



Apr 28, 2025

Human Pluripotent Stem Cell Culture

 [Cell Genomics](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.kqdg334pqq25/v1>

Anita Adami¹

¹Laboratory of Molecular Neurogenetics, Department of Experimental Medical Science, Wallenberg Neuroscience Center and Lund Stem Cell Center, BMC A11, Lund University, 221 84 Lund, Sweden.



Laura Castilla-Vallmanya

Lund University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kqdg334pqq25/v1>

External link: <https://doi.org/10.1016/j.xgen.2025.100979>

Protocol Citation: Anita Adami 2025. Human Pluripotent Stem Cell Culture. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.kqdg334pqq25/v1>

**Manuscript citation:**

Adami A, Garza R, Gerdes P, Johansson PA, Dorazehi F, Koutounidou S, Castilla-Vallmanya L, Atacho DA, Sharma Y, Johansson JG, Tam O, Kirkeby A, Barker RA, Hammell MG, Douse CH, Jakobsson J (2025) LINE-1 retrotransposons mediate cis-acting transcriptional control in human pluripotent stem cells and regulate early brain development. Cell Genomics 5(10). doi: [10.1016/j.xgen.2025.100979](https://doi.org/10.1016/j.xgen.2025.100979)

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 24, 2025

Last Modified: April 28, 2025

Protocol Integer ID: 152744

Keywords: ASAPCRN, human pluripotent stem cell culture this protocol, human pluripotent stem cell culture, pluripotent stem cell, induced pluripotent stem cell, stem cell, thawing

Funders Acknowledgements:

Aligning Science Across Parkinson's through the Michael J. Fox Foundation for Parkinson's Research

Grant ID: ASAP-000520

Aligning Science Across Parkinson's through the Michael J. Fox Foundation for Parkinson's Research

Grant ID: ASAP-025170

Abstract

This protocol described how to perform routine maintenance of human induced pluripotent stem cells (iPSCs), including freezing and thawing.



Coating Plates with Laminin-521

- 1 Coat your plates with Laminin-521 ($0.7 \mu\text{g}/\text{cm}^2$; Biolamina) diluted in DPBS +/- at 37°C for at least 2 hours before use or at 4°C overnight.

Thaw Frozen Cells

- 2 Take your vials of frozen cells and keep them on dry ice until ready.
- 3 Prepare one 15 ml tube with wash buffer (9.5 ml DMEM/F-12 + 0.5 ml KSR) per vial.
- 4 Put the vials in 37°C water bath. Check the vial periodically. Take vial out when there is only a small chunk of ice left.
- 5 Spray the vial with 70% EtOH. Wipe dry. In the hood, open the vial and take the cells out of the vial with a P1000 into your tubes with wash media.
- 6 Pellet cells in the conical tube at 400g for 5 min. Take off supernatant. Resuspend the cell pellet in appropriate amount of iPS media supplemented with Rock inhibitor (StemMACS iPS-Brew XF, 0.5% penicillin/streptomycin and $10 \mu\text{M}$ of Y27632).
- 7 Plate the cells onto a Laminin-521 coated plate and put the plate back into the incubator.

Regular Culture of Human Pluripotent Stem Cells

- 8 Feed cells with pre-warmed iPS media every day (StemMACS iPS-Brew XF with 0.5% penicillin/streptomycin).
- 9 Split the cells when they are 70-90% confluent.
- 9.1 Take off media and rinse with 1 ml PBS -/-.



- 9.2 Add 500 μ l of Accutase and incubate at 37°C for 7-10 minutes.
- 9.3 Monitor cell detaching under the microscope.
- 9.4 Harvest your cell suspension into a 15 ml tube with pre-warmed wash buffer.
- 9.5 Pellet cells in the conical tube at 400g for 5 min. Take off supernatant. Resuspend the cell pellet in appropriate amount of iPS media supplemented with Rock inhibitor (StemMACS iPS-Brew XF, 0.5% penicillin/streptomycin and 10 μ M of Y27632).
- 9.6 Put the new plates in the incubator.

Freezing Cells

- 10 Split cells following the same steps as in the section above.
- 11 Count your cells with a cell counter and transfer the desired amount to a new 15 ml tube (we normally freeze 300,000 cells per vial) and spin at 400xg for 5 min.
- 12 Resuspend the pellet in the appropriate volume of Cryostor (1 ml per vial) and immediately transfer to the labeled cryotubes.
- 13 Put vials into Mr. Frosty and place the Mr. Frosty in a -80°C freezer as quickly as possible.
- 14 The next day, transfer all the vials to -150°C freezer or liquid nitrogen for long-term storage.