

Apr 15, 2025

Human Primary Astrocyte (hPA) CRISPRa reprogramming protocol

DOI

<https://dx.doi.org/10.17504/protocols.io.ewov12477gr2/v1>

Samuel Reisman¹, Charles Gersbach¹

¹Duke University

Gersbach Lab



Andrea R Daniel

Duke University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ewov12477gr2/v1>

Protocol Citation: Samuel Reisman, Charles Gersbach 2025. Human Primary Astrocyte (hPA) CRISPRa reprogramming protocol. protocols.io <https://dx.doi.org/10.17504/protocols.io.ewov12477gr2/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 09, 2025

Last Modified: April 15, 2025

Protocol Integer ID: 126469

Keywords: human primary astrocyte, methods for crispra, hpa, crispra, reprogramming protocol

Funders Acknowledgements:

NIH

Grant ID: HG012053

Abstract

This protocol describes methods for CRISPRa reprogramming of human primary astrocytes.



Materials

	A	B	C
	AGM	unit	500 mL
	AGM (Sciencell 1801)	mL	remove 60ml from bottle
	FBS	mL	50 mL
	AGM supplement (Sciencell 1801 supp)	mL	5 mL
	P/S (Sciencell 1801 supp)	mL	5 mL

	A	B	C	D	E
	NDM (Guo 2014)	unit	100 mL	250 mL	500 mL
	DMEM/F12 (GIBCO)	mL	97.5	243.75	487.5
	0.5% FBS (GIBCO)	mL	0.5	1.25	2.5
	3.5 mM glucose (Sigma)*	mg	63.056	157.64	315.28
	penicillin/str eptomycin (GIBCO)	mL	1	2.5	5
	N2 supplement (GIBCO)	mL	1	2.5	5
	Filter sterilize after combining				

Supplements



Aliquot prep

BDNF (Peprotech 450-02-50UG)

- Dilute 50ug in 250ul H₂O (final conc 200 ug/ml)
- Add 2.25ml 0.1% BSA (final: 2.5ml @ 20 ug/ml).
 - This should be 1000x of what most protocols call for (20 ng/ml)
- Make 50ul aliquots, store at -80C

NT-3 (Peprotech 450-03)

- Dilute 10ug in 100ul H₂O (final conc. 100 ug/ml)
- Add 900ul sterile PBS + 0.1% BSA (final: 1ml @ 10ug/ml).
 - This should be 1000x of what most protocols call for (10 ng/ml)
 - Make 50ul aliquots, store at -80C
- To make BSA PBS
 - 0.1g BSA in 10ml PBS, filter sterilize (0.22um filter)



1

	A	B	C
	Day	Action	Notes
	-1	Coat plates with PDL	PDL is stored at -20 at 100x in 400ul aliquots (dilute to 40ml with sterile H2O)
	0	Prepare plates and transduce cells	
	1	Change media (AGM)	
	2	Change media (AGM + 1ug/ml Puro)	
	3	Change media (NDM + 1ug/ml Puro)	
	4	Check cells, refresh media if there is lots of debris	I usually change media here
	5	Check cells, refresh media if there is lots of debris	I usually change media here
	6	Change media to NDM without puro	
	7	Check cells, refresh media if there is lots of debris	I usually do not change media here



	A	B	C
	8	Change media to NDM + BDNF (20ng/ml) + NT-3 (10ng/ml)	BDNF and NT3 stocks are both 1000x and are stored at -80C in 50ul aliquots
	9	Check cells, refresh media if there is lots of debris	I usually do not change media here
	10	Harvest. Proceed to sorting or RNA extraction	
	Be very gentle with all media changes and rinses		