

Sep 09, 2019

Human skin single cell dissociation

DOI

<https://dx.doi.org/10.17504/protocols.io.ripd4dn>

James Fletcher¹, Rachel Botting¹, Emily Stephenson¹, Peter Vegh¹, Muzlifah Haniffa¹

¹ICM, Newcastle University



James Fletcher

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ripd4dn>

Protocol Citation: James Fletcher, Rachel Botting, Emily Stephenson, Peter Vegh, Muzlifah Haniffa 2019. Human skin single cell dissociation. protocols.io <https://dx.doi.org/10.17504/protocols.io.ripd4dn>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: July 06, 2018



Last Modified: September 09, 2019

Protocol Integer ID: 13615

Keywords: single cell dissociation splitting dermi, single cells from human skin, dissociating single cell, human skin, epidermis, cell

Abstract

Splitting dermis from epidermis, then dissociating single cells from human skin.



- 1 Cut skin into thin strips, ~1.5cm wide
- 2 Stretch the tissue using forceps, and using an 80µm guard on a dermatome cut the top layer of skin. Place this in ice cold PBS
- 3 Cut a mesh of slits into the 200µm skin layer to aid enzymatic access
- 4 Place the skin strips in a 50ml Falcon tube in RPMI + 2U/ml Dispase II
- 5 Invert the tubes to ensure the dispase has sufficient access to all surfaces
- 6 Seal the tubes with parafilm and place in a 37°C water bath for 1 hour. Invert the tubes every 20 minutes to aid digestion
- 7 Dunk the skin sheets in PBS to dilute the dispase
- 8 Use fine forceps to peel the epidermis from the dermis, placing each in PBS in separate containers on ice while peeling
- 9 Place the epidermis in a cell culture dish with 20ml RF10 (RPMI + 10% FCS (v/v) + 100U/ml Penicillin + 100µg/ml Streptomycin + 1% L-Glutamine (v/v). Add 1.6mg/ml type IV Collagenase. Place this in a 37°C 5% CO₂ incubator overnight (~12 hours)
- 10 Repeat step 9 with the dermis
- 11 Collect the supernatant with a serological pipette and pass it through a 100µm cell strainer, keeping the dermis and epidermis separate
- 12 Wash the base of the culture dish with RF10 to collect any settled cells and pass this solution through a 100µm cell strainer



- 13 Centrifuge the collected cells at 500g for 5 minutes (high acceleration/deceleration) and resuspend in PBS to count