

Jul 28, 2023

## Human urine procurement and processing

DOI

<https://dx.doi.org/10.17504/protocols.io.dm6gp3bmjvzp/v1>

hannah.anvari<sup>1</sup>, asma.giornazi<sup>1</sup>, shriya.shah<sup>1</sup>, maryellen.pavone<sup>1</sup>, francesca.e.duncan francesca.duncan<sup>1</sup>

<sup>1</sup>Northwestern Medicine OBGYN, Northwestern University



hannah.anvari

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.dm6gp3bmjvzp/v1>

**Protocol Citation:** hannah.anvari , asma.giornazi , shriya.shah , maryellen.pavone , francesca.e.duncan francesca.duncan 2023. Human urine procurement and processing. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.dm6gp3bmjvzp/v1>

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** July 28, 2023

**Last Modified:** July 28, 2023

**Protocol Integer ID:** 85637

**Keywords:** human urine procurement, storage of human urine, human urine, processing purpose, protocol, follicular fluid, processing, purpose, tissue, storage, sample, use in the collection

**Funders Acknowledgements:**

SenNet

Grant ID: U54AG075932

## Abstract

**Purpose:** This protocol is intended for use in the collection and storage of human urine in a research setting. The protocol details the collection, processing, and long-term storage for downstream analysis alongside other tissues or samples of interest (i.e. follicular fluid, ovary).

## Guidelines

Researchers will adhere to all safety and training protocols required by the institution (Northwestern Medicine/Northwestern University) including but not limited to:

1. Biosafety Certification
2. Bloodborne Pathogens Certification
3. Collaborative Institutional Training Initiative (CITI program) Certification

## Materials

1. Parter Medical Products™ Sterile Specimen Cups (Fisher Scientific, 22-150-265 or equivalent)
2. Corning 5 mL Internal Threaded Polypropylene Cryogenic Vial, Self-Standing with Round Bottom (Fisher Scientific, 09-761-69 or equivalent)


## Safety warnings

- ! Researchers will wear personal protective equipment when working with human specimens, including gloves, masks, and lab coats.

## Ethics statement

Human urine procurement and processing will adhere to the approved IRB protocol approved through NU (NU12G09) for collection of biofluids associated with study participants through Northwestern Medicine.

- 1 A urine sample is collected in a standard specimen cup during pre-op from a nurse. Once the urine sample has been collected it is placed immediately on ice following collection. A research coordinator will bring the urine sample to the lab on ice (2-4 °C). Urine samples should be transported in a blue, biohazard labeled cooler bag (Fig. 1A, B).



**Fig. 1.** Urine collection and processing overview. Urine samples are collected and transported to the research coordinator on ice in a blue biohazard labeled cooler bag (A, B). Upon arrival a photo is taken of the urine biospecimen cup (C) the urine sample is then immediately aliquoted into 5ml cryovials (D) and then stored at -80°C until analyzed.

- 3 On arrival the urine sample biospecimen cup is photographed and then aliquoted from the specimen cup (Fig. 1C) into 5 mL cryovials (Fig. 1D) and then stored at -80°C until analysis.