



Sep 17, 2024

Hybrid selection protocol using 10x Single-Cell RNA-Seq assay library

DOI

<https://dx.doi.org/10.17504/protocols.io.x54v924jzl3e/v1>

Xian Adiconis¹, Joshua Levin¹

¹Broad Institute of MIT and Harvard



xian Adiconis

Broad Institute

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.x54v924jzl3e/v1>

Protocol Citation: Xian Adiconis, Joshua Levin 2024. Hybrid selection protocol using 10x Single-Cell RNA-Seq assay library. [protocols.io https://dx.doi.org/10.17504/protocols.io.x54v924jzl3e/v1](https://dx.doi.org/10.17504/protocols.io.x54v924jzl3e/v1)

Manuscript citation:

Simmons SK, Adiconis X, Haywood N, Parker J, Lin Z, Liao Z, Tuncali I, Al'Khafaji AM, Shin A, Jagadeesh K, Gosik K, Gatzert M, Smith JT, El Kods DN, Kuras Y, Baecher-Allan C, Serrano GE, Beach TG, Garimella K, Rozenblatt-Rosen O, Regev A, Dong X, Scherzer CR, Levin JZ. Experimental and Computational Methods for Allelic Imbalance Analysis from Single-Nucleus RNA-seq Data. bioRxiv [Preprint]. 2024 Aug 16:2024.08.13.607784. doi: 10.1101/2024.08.13.607784. PMID: 39185246; PMCID: PMC11343128.

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 10, 2024

Last Modified: September 17, 2024

Protocol Integer ID: 107222

Keywords: single-cell RNA-seq, 10x Chromium, hybrid selection, targeted sequencing, end multiplexed sequencing library protocol version c3, seq assay library this protocol, target enrichment of cdna library, 10x genomic, hybrid selection protocol, seq assay library, sureselectxt target enrichment system, seq assay, sureselectxt target enrichment system for illumina, cell rna, rna, cdna library, cg000059-demonstratedprotocollexome-revc, target enrichment, hybridization, using 10x single, seq, hybrid

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-000301

Abstract

This protocol is for target enrichment of cDNA libraries generated with 10x Genomics single-cell RNA-seq assays. The protocol consists of parts of vendor provided protocols with minor modification. The hybridization and PCR part is based on "CG000059_DemonstratedProtocolExome_RevC" (10x Genomics) and the capture part is based on "SureSelectXT Target Enrichment System for Illumina Paired-End Multiplexed Sequencing Library Protocol Version C3, September 2019" (Agilent).

Materials

Oligos

	A	B	C
	Oligo Name	Vendor	Sequence
	P5 primer	IDT	5'-AATGATACGGCGACCACCGA-3'
	P7 primer	IDT	5'-CAAGCAGAAGACGGCATACGA-3'

Reagents

	A	B	C
	Reagents	Vendor	Part Number
	xGen® Universal Blocking Oligo – TS-p5, 25 rxn	IDT	1016184
	xGen® Universal Blocking Oligo – TS-p7(8nt), 25 rxn	IDT	1016188
	SureSelectXT Reagent kit for 16 samples, includes library prep and enrichment reagents for post-capture processing on the HiSeq platform. Includes indexes 1 - 16.	Agilent	G9611A
	SureSelect Custom Tier3 3Mb-5.9Mb Probe (up to 180k oligos) sufficient for post-capture processing of 16 samples. With the following configuration: Design/ ELID # : S3333112	Agilent	5191-6910
	Dynabeads MyOne Streptavidin T1, 2ml	Thermo Fisher Scientific	65601
	Amp Mix	10x Genomics	220129
	Agencourt AMPure® XP SPRI beads, 60 ml	Beckman Coulter Genomics	A63881

Safety warnings

⚠ For hazard information and safety warnings, please refer to the MSDSs (Material Safety Data Sheets).

Sample preparation

30m

- 1 Use  750 ng of cDNA library generated with the 10x Chromium single-cell RNA-seq assay.
- 2 Add  1 μL each of TS-p5 and TS-p7 blocking oligos.
- 3 If the total volume exceeds  3.4 μL , use a speed-vac and set the temperature at  30 $^{\circ}\text{C}$ to reduce the volume to the desired volume of  3.4 μL .

20m

Buffer and reaction mix preparation

20m

- 4 Prepare the Hybridization Buffer with reagents from the SureSelectXT kit at  Room temperature
Prepare the volume for at least 5 reactions to ensure accurate pipetting.

10m

	A	B	C
		1x (μL)	5x (μL)
	Hyb1 (orange)	6.63	33.15
	Hyb2 (red)	0.27	1.35
	Hyb3 (yellow)	2.65	13.25
	Hyb4 (black)	3.45	17.25
	Total	13	65

- 5 Prepare the block mix at  4 $^{\circ}\text{C}$

10m

	A	B
		1x (μL)
	Indexing Block1 (green)	2.5



	A	B
	Block 2 (blue)	2.5
	H2O	0.6
	Total	5.6

Sample and block denaturing

5m

- 6 Add  5.6 μL Block mix to the concentrated  3.4 μL sample, mix well and transfer to the 0.2ml PCR tube strip (Axygen), and place into a thermocycler following Thermocycler Conditions: (Lid@  105 $^{\circ}\text{C}$, 100 μl)

 95 $^{\circ}\text{C}$,  00:05:00

 65 $^{\circ}\text{C}$,  00:00:00 hold

Hybridization

1d

- 7 While step 6 mix is incubated at  65 $^{\circ}\text{C}$ for  00:05:00 , prepare the following at  Room temperature

	A	B
		1x (μl)
	Hybridization Buffer (Step 4)	13
	25% RNase Block solution (for $\geq 3\text{Mb}$)	0.5 μl RNase Block (purple)+1.5 μl H ₂ O=2 μl
	Capture library($\geq 3\text{ Mb}$)	5
	Total	20

Mix well by high speed vortexing for  00:00:05 , spin down briefly, add into the sample and block mix at  65 $^{\circ}\text{C}$, mix by pipetting 8-10x, use a new cap strip to seal the tubes

Incubate at  65 $^{\circ}\text{C}$ for  16:00:00 to  24:00:00



Library capture

2h

- 8 In a strip tube, use  50 μL MyOne Streptavidin T1 beads per sample, wash with  200 μL SureSelect Binding buffer 3 times, resuspend in  200 μL SureSelect Binding buffer. 15m
- 9 Bring the washed beads from step 8 near the thermocycler, add the  65 °C reaction mix (~  27 μL left) right into the  200 μL beads, gentle pipetting mix, incubate at  Room temperature for  00:30:00 , pipetting mix 6 times every  00:05:00 35m
- 10 Place on a magnet stand to pellet the beads and remove the supernatant once the solution appears clear. Resuspend the beads in  200 μL SureSelect Wash buffer 1, incubate at  Room temperature for  00:15:00 15m
- 11 Meanwhile, pre-warm SureSelect Wash buffer 2 at  65 °C in strip tubes with  200 μL per well and 3 wells per reaction on a 96-well heating block (or a thermocycler with reaction volume capacity of  200 μL) 5m
- 12 Place step 10 reaction mix on a magnet stand to pellet the beads and remove the supernatant once the solution appears clear.
- 13 Resuspend the beads in  200 μL pre-warmed SureSelect Wash buffer 2, incubate at  65 °C for  00:10:00 Place the reaction mix on a magnet stand to pellet the beads and remove the supernatant once the solution appears clear. 15m
- 14 Repeat step 13 two more times and total three washes with pre-warmed SureSelect Wash buffer 2. And make sure all the wash buffer has been removed in the final wash. 30m
- 15 Resuspend the beads in  30 μL H₂O and keep on ice. 5m

Enrichment PCR

2h



16 On ice, assemble the following mix

30m

A	B
Amp Mix	50
P5 primer (10 μ M)	2
P7 primer (10 μ M)	2
*cDNA containing beads	30
H ₂ O	16
Total	100

Thermocycler Conditions: (Lid@ 🔥 105 °C , 100 μ l)

🔥 98 °C , ⌚ 00:00:45

then 9 cycles of

🔥 98 °C , ⌚ 00:00:15

🔥 60 °C , ⌚ 00:00:30

🔥 72 °C , ⌚ 00:00:30

then

🔥 72 °C , ⌚ 00:01:00

🔥 4 °C , ⌚ 00:00:00 hold

17 Place the PCR reaction tube on a magnet stand, wait until it clears and recover

5m

🧪 100 μ L supernatant to a new PCR tube strip.18 Add 🧪 180 μ L SPRI beads (1.8x) to the recovered PCR product, mix well by pipetting, pellet the beads on a magnet stand, after supernatant removal, wash with 🧪 300 μ L 80% ethanol twice, elute with 🧪 20 μ L buffer EB (10 mM Tris-HCl, pH 8.5).

15m

19 QC with Quant-it (Thermo Fisher Scientific) and BioAnalyzer DNA HS assay (Agilent).

1h



Protocol references

1. "CG000059_DemonstratedProtocolExome_RevC"

https://assets.ctfassets.net/an68im79xiti/Zm2u8VIFa8qGYW4SGKG6e/4bddcc3cd60201388f7b82d241547086/CG000059_DemonstratedProtocolExome_RevC.pdf

2. "SureSelectXT Target Enrichment System for Illumina Paired-End Multiplexed Sequencing Library Protocol Version C3, September 2019"

newer version at <https://www.agilent.com/cs/library/usermanuals/public/G7530-90000.pdf>