

Nov 01, 2019

Hybridization chain reaction (HCR) protocol for tails of mouse embryos

DOI

<https://dx.doi.org/10.17504/protocols.io.7pyhmpw>


Paul Gerald Layague Sanchez¹, Hidenobu Miyazawa¹, Molecular Instruments²

¹European Molecular Biology Laboratory (EMBL); ²Molecular Instruments



Paul Gerald Layague Sanchez

European Molecular Biology Laboratory (EMBL)



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.7pyhmpw>

Protocol Citation: Paul Gerald Layague Sanchez, Hidenobu Miyazawa, Molecular Instruments 2019. Hybridization chain reaction (HCR) protocol for tails of mouse embryos. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.7pyhmpw>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol in our group and it is working.

Created: September 26, 2019

Last Modified: November 01, 2019

Protocol Integer ID: 28120

Keywords: HCR, in situ hybridization, mouse embryos this protocol, situ hybridization chain reaction, multiplexed zoo of mrna expression, mouse embryo, mount mouse embryo, mount mouse embryos available in the website, hybridization chain reaction, generationin situhybridization chain reaction, mrna expression, situ hcr, multiplexed zoo, hcr

Abstract

This protocol is for the in situ hybridization chain reaction (HCR) of intact tails (somites and presomitic mesoderm, PSM) of mouse embryos. This is adapted from the in situ HCR v 3.0 protocol for whole-mount mouse embryos available in the website of Molecular Instruments (for details, please follow this [link](#)).

References

- Choi, H., Schwarzkopf, M., Fornace, M., Acharya, A., Artavanis, G., & Stegmaier, J. et al. (2018). Third-generationin situhybridization chain reaction: multiplexed, quantitative, sensitive, versatile, robust.*Development*,145 (12), dev165753. doi: [10.1242/dev.165753](https://doi.org/10.1242/dev.165753)
- Choi, H., Calvert, C., Husain, N., Huss, D., Barsi, J., & Deverman, B. et al. (2016). Mapping a multiplexed zoo of mRNA expression.*Development*,143 (19), 3632-3637. doi: [10.1242/dev.140137](https://doi.org/10.1242/dev.140137)

Day 1: Recovery and fixation of samples

- 1 Recover samples, wash quickly with cold PBS, and fix in 4% formaldehyde ( 1.08 mL 37% formaldehyde diluted to  10 mL with PBS). Keep in formaldehyde overnight (around  16:00:00) at  4 °C in the cold room, with gentle shaking/rolling.

Safety information

Formaldehyde is toxic. Operate under a fume hood and wear appropriate hand, body, and eye protection.

Note

To reduce autofluorescence, it is recommended to use fresh fixative and PBS with no divalent cations (like Ca^{2+} and Mg^{2+}).

Day 2: Dehydration of samples

- 2 Wash samples with PBST (0.1% Tween20 in 1x PBS).
 - Wash 1/2  00:05:00
 - Wash 2/2  00:05:00
- 3 Dehydrate samples (one sample per tube) with a series of graded methanol : PBST (PBS + 0.1% Tween) washes.
 - 25% methanol : 75% PBST  00:05:00
 - 50% methanol : 50% PBST  00:05:00
 - 75% methanol : 25% PBST  00:05:00
 - 100% methanol  00:05:00
 - 100% methanol  00:05:00
- 4 Incubate samples at  -20 °C overnight (around  16:00:00) or until use.



Day 3: Rehydration of samples and HCR detection stage

5 Rehydrate samples (one sample per PCR tube, capacity =  200 μ L max) with a series of graded methanol : PBST (PBS + 0.1% Tween) washes.

- 100% methanol  00:05:00
- 75% methanol : 25% PBST  00:05:00
- 50% methanol : 50% PBST  00:05:00
- 25% methanol : 75% PBST  00:05:00
- 100% PBST  00:05:00
- 100% PBST  00:05:00

6 Treat samples with  10 μ g/mL proteinase K (Merck, CAS # 38450-01-6) for  00:05:00 at  Room temperature .

7 Wash samples with PBST.

- Wash 1/2  00:05:00
- Wash 2/2  00:05:00

8 Fix samples in 4% formaldehyde ( 2.16 mL 37% formaldehyde diluted to  20 mL with PBST) for  00:20:00 at  Room temperature .

Pre-warm probe hybridization (PH) buffer at  37 $^{\circ}$ C for later use.

Safety information

Formaldehyde is toxic. Operate under a fume hood and wear appropriate hand, body, and eye protection.

Safety information

Probe hybridization (PH) buffer contains formamide, which is toxic.

**Note**

To reduce autofluorescence, it is recommended to use fresh fixative and PBS with no divalent cations (like Ca^{2+} and Mg^{2+}).

9 Wash samples with PBST.

- Wash 1/3 00:05:00
- Wash 2/3 00:05:00
- Wash 3/3 00:05:00

10 Wash samples with PH buffer for 00:05:00 .

Safety information

Probe hybridization (PH) buffer contains formamide, which is toxic.

11 Put fresh PH buffer to the tubes and incubate for 00:30:00 at 37 °C .

Meanwhile, prepare probe solution by adding 2 pmol (2 μL of
 [M] 1 micromolar (μM) stock) of odd probe mixture and 2 pmol of even probe mixture to every 500 μL of pre-warmed PHP buffer (from Step 8).

Safety information

Probe hybridization (PH) buffer contains formamide, which is toxic.

12 Remove PH buffer and add the probe solution. Incubate overnight (around 16:00:00) at 37 °C .

Day 4: HCR amplification stage

13 Pre-warm probe wash buffer at 37 °C , and equilibrate amplification buffer to Room temperature .

**Safety information**

Probe wash (PH) buffer contains formamide, which is toxic.

- 14 In separate tubes, for every  500 μL hairpin mixture, prepare 30 pmol ( 15 μL of  2 micromolar (μM) stock or  10 μL of  3 micromolar (μM) stock) of hairpin h1 (tube 1) and 30 pmol of hairpin h2 (tube 2). Snap-cool hairpins by heating the tubes at  95 °C for  00:01:30 and cool to  Room temperature in the dark for at least  00:30:00 .

- 15 Wash samples with probe wash buffer at  37 °C .

- Wash 1/4  00:15:00
- Wash 2/4  00:15:00
- Wash 3/4  00:15:00
- Wash 4/4  00:15:00

Safety information

Probe wash (PH) buffer contains formamide, which is toxic.

- 16 Wash samples with 5x SSCT (0.1% Tween20) at  Room temperature .

- Wash 1/2  00:05:00
- Wash 2/2  00:05:00

Note

To prepare  500 mL of 5x SSCT, dilute  125 mL of 20x SSC and  0.5 mL Tween20 to  500 mL with distilled H_2O .



- 17 Wash samples with pre-equilibrated amplification buffer (Step 13) for 00:05:00 at Room temperature .

Meanwhile, prepare the hairpin mixture (30 pmol of each hairpin in 500 μ L mixture) by mixing the snap-cooled hairpins h1 and h2 (Step 14) to amplification buffer at Room temperature .

- 18 Remove the amplification buffer from the sample tubes, and add the hairpin mixture. Incubate overnight (around 16:00:00) in the dark at Room temperature .

Note

For dHCR imaging, amplify for a shorter period of time to ensure single-molecule dots are diffraction-limited.

Day 5: Washing, nuclear staining, and imaging

- 19 Wash samples with 5x SSCT for 00:30:00 in the dark at Room temperature .

- 20 Stain nuclei with 1:1000 DAPI (5 μ g/mL : 20 μ L of 5 mg/mL DAPI diluted to 20 mL with 5x SSCT) in the dark at Room temperature .

- Wash + DAPI 1/3 00:30:00
- Wash + DAPI 2/3 00:30:00
- Wash + DAPI 3/3 00:30:00

- 21 Wash samples with 5x SSCT in the dark at Room temperature

- Wash 1/2 00:30:00
- Wash 2/2 00:30:00

- 22 Samples are now ready for imaging.

Prior to imaging, the samples could be cleared using the fructose-glycerol clearing solution described in [Dekkers et al., 2019](#). Incubate samples in clearing solution in the




dark at 🌡️ 4 °C for at least ⌚ 02:00:00 .

Note

To prepare the fructose-glycerol clearing solution, dissolve 🧪 29.72 g fructose in 🧪 33 mL glycerol and 🧪 7 mL water on a magnetic stirrer (start the dissolution at least a day before intended day of clearing because it will take a while for fructose to dissolve). The clearing solution can be stored in the dark at 🌡️ 4 °C for at most 1 month. This solution can be used to mount the sample on microscope slide/cover glass. Take note that the clearing is solution has high viscosity.

Expected result



Tbx18 and Uncx marking the two halves (anterior and posterior, respectively) of somites in tail of a mouse embryo. Transcripts were visualized using HCR version 3.0. Nuclei were stained using DAPI.