

Oct 19, 2019

Hybridization of DNA oligos

DOI

dx.doi.org/10.17504/protocols.io.8guhtww

Laura Sánchez¹

¹AEGIS - Madrid iGEM Team 2019

AEGIS - Madrid iGEM 2019



Manuela De Las Casas

Universidad Complutense de Madrid, AEGIS - Madrid iGEM 2019

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.8guhtww

Protocol Citation: Laura Sánchez 2019. Hybridization of DNA oligos. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.8guhtww>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 19, 2019

Last Modified: October 19, 2019

Protocol Integer ID: 28916

Abstract

This protocol is used to hibridate two complementary DNA chains.



Materials

Equipment

- Thermocycler
- Vortex
- Centrifuge

Reagents

- 1.0 M NaCl
- Sterile H₂O

DNA

- Promoters 5'
- Promoters 3'

Safety warnings

- ❗ Lab coat and gloves should be weared throughout the whole experiment. All working surfaces must be clean and all the reactives should be treated following manufacturer's instructions.

Annealing of oligonucleotides

- 1 Add the required volume of H₂O to the lyophilized oligonucleotides to obtain a concentration of 100 µM.
Vortex both tubes for 30 s and incubate them at RT for 5 min to dissolve them.
- 2 Prepare the annealing mix by adding into a PCR tube:
 - > 45.5 µl of H₂O
 - > 2.5 µl of 1.0 M NaCl
 - > 1 µl of oligo 5' (100 µM)
 - > 1 µl of the oligo 3' (100 µM)
- 3 Place the PCR tube with the mix in a thermocycler with the following annealing programme:
 - 5 min at 95°C
 - 1 min at 95°C
 - ramp down 1°C per cycle for 72 cycles
 - end by keeping the temperature at 10°C
- 4 Take 10 µl of the annealed oligonucleotides and dilute it with 90 µl of H₂O to obtain a 0.2 µM concentration.
- 5 The annealed oligonucleotides stocks can be stored at -20°C for future use.