

May 05, 2025

Version 2

HyDrop v2 Bead Generation & Ligation Barcoding V.2

 [Nature Communications](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.8epv5n11dv1b/v2>

Suresh Poovathingal¹, Marta Wojno¹, Koen Theunis², Florian De Rop^{2,3,4,5}, Stein Aerts^{2,3,4,5}

¹VIB-KU Leuven Center for Brain & Disease Research, CBD Technologies, Single Cell & Microfluidics Expertise Unit, 3000 Leuven, Belgium;

²Laboratory of Computational Biology, VIB Center for AI & Computational Biology, Leuven, Belgium;

³VIB-KU Leuven Center for Brain & Disease Research, Leuven, Belgium;

⁴Department of Human Genetics, KU Leuven, Leuven, Belgium;

⁵Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD 20815, United States



Koen Theunis

Laboratory of Computational Biology, VIB Center for AI & Com...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.8epv5n11dv1b/v2>

External link: <https://doi.org/10.1038/s41467-026-68742-4>

Protocol Citation: Suresh Poovathingal, Marta Wojno, Koen Theunis, Florian De Rop, Stein Aerts 2025. HyDrop v2 Bead Generation & Ligation Barcoding. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.8epv5n11dv1b/v2> Version created by **Koen Theunis**

Manuscript citation:

Dickmanken H, Wojno M, Mahieu L, Theunis K, Ekşi EC, Christiaens V, Kempynck N, Rop FVD, Roels N, Spanier KI, Vandepoel R, Hulselmans G, Poovathingal S, Aerts S (2026) Evaluating single-cell ATAC-seq atlasing technologies using sequence-to-function modeling. Nature Communications 17(). doi: [10.1038/s41467-026-68742-4](https://doi.org/10.1038/s41467-026-68742-4)

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 05, 2025

Last Modified: May 05, 2025

Protocol Integer ID: 217720

Keywords: ASAPCRN, hydrogel, microfluidics, hydrop, HyDrop, scATAC, scRNA, scOmics, hydrogel bead, hydrop v2 bead generation, protocol details hydrogel, single bead mixture, ml of bead, emulsion of acrylamide monomer, acrylamide monomer, hydrop v2, primer concentration within the bead, bead, barcode, bead stock, monomer droplet, carrier oil, emulsion, containing temed, ligation strategy, barcoding, ligation barcoding this protocol details hydrogel, ligation barcoding

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-025179

Aligning Science Across Parkinson's

Grant ID: ASAP-000430

Abstract

This protocol details hydrogel bead generation and barcoding for HyDrop v2. The barcodes are build by split-pool and ligation strategy (3×96). Here, we will create an emulsion of acrylamide monomers in a carrier oil containing TEMED. The monomer droplets will polymerise and form hydrogel beads. Ideally, you have a bead stock of around 3 mL of beads before you barcode, but 2 mL can work as well. It is best to produce the beads in one single run, from a single bead mixture to prevent disparities in sizes and/or primer concentration within the bead from occurring.

Attachments



[HyDrop-buffer-listin...](#)

12KB



[plate-1-Fwd-96.xls](#)

35KB



[plate-2-Fwd-96.xls](#)

36KB



[plate-2-Rev-96.xls](#)

36KB



[plate-3-ATAC-96.xls](#)

37KB



[acrydite primer.xlsx](#)

8KB



[plate-1-Rev-96.xlsx](#)

11KB



Materials

- 96-well plate
- Cassette HyDrop-v2-ATAC-plate
- Tris
- Tris-HCl
- NaCl
- NaOH
- MgCl_2
- Tween-20
- dH_2O
- Oligo-Mix BC1
- ATP
- T4 ligase
- Oligo-Mix B
- Brij-35
- TX-100

Before start

On the days before the barcoding prepare the primer cassette for the stage 1 and stage 2 barcoding.

Hydrogel bead generation

1

Prepare monomer mix. The volume used here will approximately equate the final volume of dense bead stock that you will produce. The beads will be slightly larger than the droplets generated, but you will also generate small losses during the wash steps.

	A	B	C	D	E
		Vol (uL)	Stock	Final	Unit
	TBSET	110	1	0.1	X
	Acrylamide	269.5	40	10	%w/v
	N,N'-Bis(acryloyl)cystamine	45.83	4.8	0.2	%w/v
	Ammonium persulfate	66	10	0.6	%
	Acrydite primer bead handle	44	100	4	uM
	Nuclease free water	564.67			
	TOTAL	1100			

The volumes used here equals to 10% acrylamide (%T) and 2% bisacryloylcystoylamine (%C). Note that N,N'-Bis(acryloyl)cystamine is dissolved w/v in methanol.

- 2 Add 1000 μ L of filtered HFE-7500 Novac oil to the syringe (3 mL) and add the monomer mix and layer it on top of the HFE and prepare the syringe for injection.
- 3 Mix 2900 μ L of the EA008 surfactant with 14.5 μ L of TEMED (=0.5% TEMED). Mix the content well.
- 4 Prepare two 1.5 mL Eppendorf collection tubes and put 200 μ L of mineral oil in the tubes. This will also serve as protection against evaporation during the overnight incubation afterwards.
- 5 Run the droplet generation for the preparation of the hydrogel beads.



Flow rates:

	A	B	C
		start	stable
monomer mix		400	400
oil		400	750

Collect the emulsion for 55 minutes in two 1.5 mL Eppendorf tubes.

- 6 Remove excess oil from the bottom of the collection tubes and double check if the layer of mineral oil is on top of the emulsions.
- 7 Put the emulsions on a 65°C heat block for 14 hours. Leave beads at 12°C until they are washed.
- 8 Remove top layer of mineral oil and bottom layer of EA008.
- 9 Wash the beads:
 - 3 times with 20% PFO in HFE-7500.
 - 3 times with 1% SPAN80-Hexane.
 - 3 times with TBSET buffer.
 - 3 times with TET buffer and combine all aliquots in the last wash.
- 10 Filter the beads through a 70 um filter to exclude contamination of large beads.
- 11 Store the beads at 4°C.

Oligo preparation (done at least a day prior to barcoding)

- 12 Take 3 fresh 96 well plates. And label as:
 - Cassette 1-HyDrop-v2-ATAC-plate (200 µM)
 - Cassette 2-HyDrop-v2-ATAC-plate (200 µM)
 - Plate 3-HyDrop-v2-ATAC-plate (500 µM)
- 13 For **Cassette 1-HyDrop-v2-ATAC-plate**:
Mix  20 µL of primers from plate-1-Fwd-96 of the  400 micromolar (µM) stock primer plate to  20 µL of primers from plate-1-Rev-96 of the

[M] 400 micromolar (μM) stock primer.

14 For **Cassette 2-HyDrop-v2-ATAC-plate**:

Mix  20 μL of primers from plate-2-Fwd-96 of the [M] 400 micromolar (μM) stock primer plate to  20 μL of primers from plate-2-Rev-96 of the [M] 400 micromolar (μM) stock primer.

15 For **Plate 3-HyDrop-v2-ATAC-plate**:

Aliquot  18 μL of primers from plate-3-ATAC-96 of the [M] 500 micromolar (μM) stock primer plate to a new plate.

Barcoding Round 1

3h 38m

- 16
- Spin down the beads stored in TET buffer.
 - Remove supernatant and wash the beads one more time with TET buffer.
 - Centrifuge at  1000 x g, 00:03:00 , brake 8, and discard the supernatant.

3m



- 17 Prepare the pre-Ligation buffer (PL buffer) below and filter with  0.2 μm filter (if not available).

	A	B	C	D	E	F
	Component	Stock	Unit	Final	Unit	Volume
	Tris 8.0	1000	mM	10	mM	5000
	NaCl	5000	mM	30	mM	3000
	MgCl ₂	1000	mM	1	mM	500
	Tween-20	10	%	0.1	%	5000
	dH ₂ O					486500
	TOTAL					500000

- 18
- Wash the beads with PL buffer 2X times.
 - Centrifuge @  1000 x g, 00:01:00 , brake 8, discard the supernatant.

1m





- 19 Using Hamilton, aliquot  22 μL of the compacted beads to each of the wells of 96 well plate.

Total compacted bead volume: $22 \times 1.2 \times 96 =$  2534.4 μL .

- 20 Prepare the 2X Ligation-primer buffer (T4) below:

	A	B	C	D	E	F	G
	Component	Stock	Unit	Final	Unit	Volume	96×1.2
	Tris 7.5	1000	mM	100	mM	2	230.4
	MgCl ₂	1000	mM	20	mM	0.4	46.08
	Cassette 1-HyDrop-v2-ATAC	200	μM	40	μM	4	NO ADD
	dH ₂ O					13.6	1566.7
						20	16 μL per reaction
							1843.2 μL

μL per reaction

- 21 Using Hamilton, aliquot  16 μL of the 2X ligation-primer buffer to the PCR plate.

- 22 ■ Add  4 μL (Hamilton) of **Cassette 1-HyDrop-v2-ATAC-plate**  
 [M] 200 micromolar (μM) , spin and mix well by vortexing and place on PCR block pre-heated to  75 °C (with heated lid to  105 °C).

- 23 Heat the thermoblock to  75 °C and keep it ready.

- 24 ■ Heat the mix at  75 °C for  00:03:00 on the PCR and transfer and switch off the heat block to cool to  Room temperature , to anneal the oligo to bead (~2.5 hours).
- Every 30 mins vortex the plate at max speed (2000x RPM) for  00:00:30 to resuspend the beads. Keep on looking after switching off for  02:30:00 (reaches  33 °C).

2h 33m
30s

- After  02:00:00 , leave the plate at  Room temperature .

Note

Take it from block and keep it at  Room temperature for  00:10:00 .

- 25 In the meantime, prepare the 1X ligase mix below (T4):

A	B	C	D	E	F	G
Component	Stock	Unit	Final	Unit	Volume	96×1.05
Tris 7.5	1000	mM	50	mM	0.5	52.8
MgCl ₂	1000	mM	10	mM	0.1	10.6
ATP	10	mM	1	mM	5.1	538.6
T4 ligase	400	U/μl	20	U/μl	2.5	252
dH ₂ O					1.8	181.4
					10	10 μl per reaction
						1035.4 μl

- 26
- Using Mantis, aliquot  10 μL of the ligase mix to the 96 well plate, after sealing the plate vortex the content in a vortex shaker, centrifuge at  1000 x g, 00:00:30 .
 - On thermoblock shake  2000 rpm, 00:00:30 , to mix the beads.

30s



- 27 Perform ligation at  25 °C for  01:00:00 with a shaking of  1000 rpm (every minute shake at  1600 rpm, 00:00:30).

1h

- 28 Perform  Overnight incubation at  16 °C on a thermoblock with shaking of  1000 rpm, 00:00:30 . <Pause point - overnight>



Barcoding Round 2

3h 59m 30s



29 Remove the plate from 16 °C incubation and leave at 4 °C till the next processing step (at least 30 min).

30 Prepare 2X ligation-primer buffer.

	A	B	C	D	E	F	G
	Component	Stock	Unit	Final	Unit	Volume	96×1.2
	Tris 7.5	1000	mM	100	mM	2	230.4
	MgCl ₂	1000	mM	20	mM	0.4	46.1
	dH ₂ O					13.6	1566.7
						20	16 µl per reaction
							1843.2 µl

31 Aliquot 16 µL of the above 2X ligation-primer buffer mix to a fresh 96 well plate.

- 32
- Add 4 µL of **Cassette 2-HyDrop-v2-ATAC-plate** [M] 200 micromolar (µM) .
 - Mix well (vortex) and place on a PCR block pre-heated to 75 °C (lid 105 °C).



- 33
- Heat the plates at 75 °C for 00:03:00 , transfer the plate to thermoblock at 75 °C and switch off the heat block to cool to Room temperature , to anneal the oligos.
 - Keep on looking after switching off for 02:30:00 .
 - After 2.5 hours, leave the plate at Room temperature .

2h 33m

Note

Take it from block and keep it at Room temperature for 00:10:00 .

- 34
- In the meantime, perform the STOP-25 cleanup on the Hamilton.
 - Transfer to a 50 ml falcon.
 - Incubate for 00:20:00 after collection.

20m





35 Spin the falcon at  3220 x g, 00:03:00 s and braking of 7. Transfer the beads to 15 ml falcon.

3m



36 Wash the beads with Stop-10 buffer. Centrifuge @  1000 x g, 00:01:00 and discard the supernatant.

1m



37 Wash 3X times with TET buffer. Centrifuge @  1000 x g, 00:01:00 and discard the supernatant.

1m



38 Wash the beads with PL buffer 2X times. Centrifuge @  1000 x g, 00:01:00 and discard the supernatant. Set the level of beads to  2534.4 μL .

1m



39 Using Hamilton, aliquot  22 μL of the compacted beads to each of the wells of 96 well plate of a new deep well 96 well plate.

40 After the completion of the oligo annealing, spin the plate to collect condensate and using Hamilton transfer the contents of the annealed oligo cassette plate to the bead plate.

41 Prepare the ligase mix below:

	A	B	C	D	E	F	G
	Component	Stock	Unit	Final	Unit	Volume	96×1.1
	Tris 7.5	1000	mM	50	mM	0.5	52.8
	MgCl ₂	1000	mM	10	mM	0.1	10.6
	ATP	10	mM	1	mM	5.1	538.6
	T4 ligase	400	U/ μl	20	U/ μl	2.5	252
	dH ₂ O					1.8	181.4
						10	10 μl
							1035.4

42 Using Mantis, aliquot  10 μL of the ligase mix to the 96 well plate, after sealing the plate vortex the content in a vortex shaker, centrifuge  1000 x g, 00:00:30 . On thermoblock shake  1600 rpm, 00:00:30 , to mix the beads.

30s



43 Perform ligation at  25 °C for  01:00:00 with  1000 rpm rotation shaking (every minute shake at  1000 rpm, 00:00:30).

1h



44 Perform  Overnight incubation at  16 °C on a thermoblock with  1000 rpm rotation shaking. <Pause point - overnight>



Barcoding Round 3

1h 32m 30s

45 Remove the plate from  16 °C incubation and leave at  4 °C till ready for subsequent processing step.

46 Perform the STOP-25 cleanup on the Hamilton. Transfer to a 50 ml falcon. Incubate for  00:20:00 after collection.

20m



47 Spin the falcon at  3220 x g, 00:03:00 and braking of 7. Transfer the beads to 15 ml falcon.

3m



48 Wash the beads with Stop-10 buffer. Centrifuge @  1000 x g, 00:03:00 and discard the supernatant.

3m



49 Wash 3X times with TET buffer. Centrifuge @  1000 x g, 00:03:00 and discard the supernatant.

3m



50 Wash the beads with PL buffer 2X times. Centrifuge @  1000 x g, 00:03:00 and discard the supernatant. Set the level of beads to  2534.4 µL .

3m



51 Using Hamilton, aliquot  23 µL of the compacted beads to a fresh 96 deep well plate.

52 Prepare the 2X Ligation buffer below:

	A	B	C	D	E	F	G
	Component	Stock	Unit	Final	Unit	Volum e	96×1.3

	A	B	C	D	E	F	G
	Tris 7.5	1000	mM	100	mM	2	249.6
	MgCl ₂	1000	mM	20	mM	0.4	49.92
	plate 3-HyDrop-v2-ATAC-plate	500	μM	30	μM	3	NO ADD
	dH ₂ O					14.6	1822.1
						20	17

53 Using Hamilton, aliquot  17 μL of the 2X ligation-primer buffer to the PCR plate.

54 Add  3 μL of **plate 3-HyDrop-v2-ATAC-plate** [M] 500 micromolar (μM) . 

55 Prepare the ligase mix below:

	A	B	C	D	E	F	G
	Component	Stock	Unit	Final	Unit	Volume	96×1.1
	Tris 7.5	1000	mM	50	mM	0.5	52.8
	MgCl ₂	1000	mM	10	mM	0.1	10.6
	ATP	10	mM	1	mM	5.1	538.6
	T4 ligase	400	U/μl	20	U/μl	2.5	252
	dH ₂ O					1.8	181.4
						10	10 μl
							1035.4

56 Using Mantis, aliquot  10 μL of the ligase mix to the 96 well plate, after sealing the plate vortex the content in a vortex shaker, centrifuge  1000 x g, 00:00:30 . On thermoblock shake  2000 rpm, 00:00:30 , to mix the beads.

30s



57 Perform ligation at  25 °C for  01:00:00 with rotation mixing of  1000 rpm .

1h

Note

Every minute shake at  1600 rpm, 00:00:30 .



58 Perform  Overnight incubation at  16 °C on thermoblock with rotation mixing of  1000 rpm . <Pause point - overnight>



Bead Cleanup

29m

59 Remove the plate from  16 °C incubation and leave at  4 °C till ready for subsequent processing step.

60 Perform the STOP-25 cleanup on the Hamilton. Transfer to a 50 ml falcon. Incubate for  00:20:00 after collection.

20m



61 Spin the falcon at  3220 x g, 00:03:00 and braking of 7. Transfer the beads to 15 ml falcon.

3m



62 Wash the beads with Stop-10 buffer. Centrifuge @  1000 x g, 00:03:00 and discard the supernatant.

3m



63 Wash 3X times with TET buffer. Centrifuge @  1000 x g, 00:03:00 and discard the supernatant.

3m



Denaturation and Final Steps

24m

64 Prepare the following denaturation buffer and filter with  0.2 µm filter.

	A	B	C	D
	Component	Stock	Final	Volume
	Water	-	-	60.5 ml



	A	B	C	D
	NaOH	10 N		938 µl
	Brij-35			1050 µl
				62.5

- 65 Add  15 mL of denaturation solution to beads and incubate for  00:10:00 on a rotator at  Room temperature .  
- 66 Wash the denaturation buffer 3X times with denaturation solution with  15 mL buffer. 
- 67 Add  15 mL of Neutralization solution to beads and incubate for  00:10:00 on a rotator at  Room temperature .  
- 68 Wash the beads with Neutralization buffer 1X more time. 
- 69 Wash 3X times with TET buffer and proceed to QC. 
- 70 Filter using  70 µm filter.
- 71 Prepare the following.
- 71.1 Wash 3X  15 mL falcon 3 times with distilled water to clear of the dust in the tube and keep (3X 15 ml Falcon) 2X Falcon is for TET buffer; 1X for Bead collection. 
- 71.2 Wash 4X  50 mL falcon 3X times with distilled water to clear of the dust in the tube and keep (4X 50 ml Falcon) 2X Falcon is for bead filtering; 1X for Bead refiltering and 1X for filtered lysis buffer collection. 
- 72 Sequential transfer of beads to  70 µm filter for 2 stage filtration.

73 Collect the final filtrate and transfer to the dust clean  15 mL falcon. Keep  On ice .

74 Prepare the 1X lysis mix in fresh 50 ml falcon.

Prepare the 1X Lysis mix

A	B	C	D	E	F
	Stock	Unit	Final drop	Final	Volume
Tris-HCl (pH 7)	1000	mM	25	125	5000
NaCl	5000	mM	30	150	1200
MgCl ₂	1000	mM	2.5	12.5	500
TX-100	100	%	0.25	1.25	500
BSA	30	%	0.08	0.4	533
dH ₂ O	-				32800
					40000

75 Filter the 1X lysis mix to a 2nd pre-clean 50 ml falcon using  0.2 µm disc filter.

76 Spin down the falcon with beads in TET at  800 x g, 00:01:00 , remove as much as supernatant as possible.

1m



77 Add  10 mL of filter 1X lysis mix. Vortex well. Spin down the falcon with beads in TET at  1000 x g, 00:01:00 . Remove and discard the supernatant.

1m



78 Add  10 mL of filter 1X lysis mix. Vortex well. Incubate the beads at  4 °C  Overnight for equilibration.



79 Next day wash the 8-strip Eppendorf tubes 2X times with the DI water. Remove as much liquid as possible with the p200 multichannel.





80

Spin down the bead at  1000 x g, 00:02:00 and aliquot  40 μ L beads to 8-strip tube.

Put the strips to  -80 °C storage.

2m

