

Dec 16, 2024

Hypothermia Induction in Mice

 In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.6qpvr937bvmk/v1>

Mukesh Kumar¹, Catherine L. Pontius², Rose E. Witte³, Kartik N. Rajagopalan², Timothy A. Ryan^{1,4}

¹Department of Biochemistry, Weill Cornell Medicine, New York, NY 10065;

²Department of Internal Medicine, UT Southwestern Medical Center, Dallas, TX 75390, USA;

³Laboratory of Molecular Genetics, Howard Hughes Medical Institute, Rockefeller University, New York, NY 10065, USA;

⁴Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, Maryland 20815, USA



Mukesh Kumar

Weill Cornell Medical College

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.6qpvr937bvmk/v1>

Protocol Citation: Mukesh Kumar, Catherine L. Pontius, Rose E. Witte, Kartik N. Rajagopalan, Timothy A. Ryan 2024. Hypothermia Induction in Mice. protocols.io <https://dx.doi.org/10.17504/protocols.io.6qpvr937bvmk/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: November 26, 2024

Last Modified: December 16, 2024

Protocol Integer ID: 114382

Keywords: ASAPCRN, hypothermia induction in mice, hypothermia in mice, hypothermia induction, measuring hypothermia, thermoregulation, pharmacological effect, using etomoxir, evaluation of pharmacological effect, mice

Funders Acknowledgements:

NIH

Grant ID: NS036942

NIH

Grant ID: NS11739

ASAP

Grant ID: ASAP-024404

KL2 Award by Rockefeller University

Abstract

This protocol describes a detailed method for inducing and measuring hypothermia in mice using Etomoxir and KLH45, enabling the evaluation of pharmacological effects on thermoregulation.

Guidelines

The protocol needs prior approval by the users' Institutional Animal Care and Use Committee (IACUC) or equivalent ethics committee.



Materials

Materials:

-  Etomoxir **Merck MilliporeSigma (Sigma-Aldrich) Catalog #236020** : 30 mg/kg, dissolved in phosphate-buffered saline (PBS).
-  KLH45 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #SML1998** : 40 mg/kg, dissolved in emulphor (18:1:1 saline: ethanol:  Polyethylene Glycol Monooleyl Ether **Fisher Scientific Catalog #P071325G**
- Control vehicles: PBS or emulphor.

Equipments:

- Heating block or water bath for solution preparation
- Syringes and needles for intraperitoneal (IP) injection
- Rectal thermometer for mice (Kent Scientific, RET-3)
- Cage with controlled environment for post-injection monitoring

Preparation of Solutions

- 1 Dissolve Etomoxir in PBS to achieve the desired concentration for 30 mg/kg dosing.
- 2 Prepare KLH45 in emulphor composed of 18:1:1 saline: ethanol: polyethylene glycol monooleyl ether to achieve a 40 mg/kg dosing concentration.
- 3 Prepare control vehicles (PBS for Etomoxir and emulphor for KLH45).

Pre-Injection Procedure

- 4 Place mice on food restriction 3 hours prior to the injection to minimize metabolic variability.
- 5 Weigh each mouse to calculate the appropriate injection volume based on body weight.

Hypothermia Induction

- 6 Administer the following intraperitoneally (IP):
 - 6.1 Etomoxir group: 30 mg/kg in PBS.
 - 6.2 KLH45 group: 40 mg/kg in emulphor.
 - 6.3 Control groups: Corresponding vehicle (PBS or emulphor).
- 7 Ensure proper injection technique to avoid discomfort or improper dosing.

Body Temperature Measurement



- 8 Measure the baseline body temperature (T_b) of each mouse using a rectal probe.
- 9 Record body temperature at the following time points post-injection: 0 hour (immediately after injection), 0.5 hour, 1 hour, and 3 hours.

Post-Injection Monitoring

- 10 Return the mice to their cages immediately after injection.
- 11 Ensure the mice are housed without access to feed during the post-injection monitoring period.
- 12 Monitor for any signs of distress or abnormal behavior during the experiment.

Data Recording

- 13 Record all temperature measurements and note any significant changes or observations during the experiment.
- 14 Analyze temperature data to assess the hypothermic response induced by the treatments.

Note

- Ensure accurate dosing by calibrating the syringe and needle prior to injection.
- Use sterile techniques when preparing solutions and administering injections.
- Mice should be handled gently to minimize stress, which can influence body temperature measurements.
- Include adequate sample size and control groups for statistical reliability.