



Jan 29, 2025

i3N neuronal differentiation

DOI

dx.doi.org/10.17504/protocols.io.n92ldrj98g5b/v1

Alexis Penverne¹

¹Gandhi Lab



Alexis Penverne

Gandhi Lab, University College London, University of London

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n92ldrj98g5b/v1>

Protocol Citation: Alexis Penverne 2025. i3N neuronal differentiation. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.n92ldrj98g5b/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 20, 2025



Last Modified: January 29, 2025

Protocol Integer ID: 118735


Keywords: neuronal differentiation this protocol, neuronal differentiation, integrated doxycycline, protocol from the kampmann lab, differentiation



Abstract


This protocol is based on the protocol from the Kampmann lab (Chen et al., 2020). It describes the how to differentiate i3Ns, iPSCs with stably integrated doxycycline-inducible Ngn2

Troubleshooting

Induction - day 0








1 iPSCs are maintained on Geltrex in mTeSR (STEMCELL) , and passaged using  0.5 micromolar (μM) EDTA. They need to be 70-80% confluent, preferably at a low passage number when starting an induction



2 Coat 6-well cell culture plates or T75 flasks with Geltrex and incubate at  37 °C for  01:00:00


3 Prepare **pre-differentiation media** supplemented with  10 micromolar (μM) **Rock inhibitors**

Note

How to prepare 50mL **Pre-differentiation media**:

-  49.8 mL Knockout DMEM/F12 (Thermo Fisher)
-  500 μL 100X N2 Supplement (Thermo Fisher)
-  500 μL Non-essential amino acids (Thermo Fisher)
-  50 μL BDNF (From 10ug/mL stock)
-  50 μL NT-3 (From 10ug/mL stock)
-  50 μL Laminin (From 1mg/mL stock)
-  50 μL Doxycyclin (From 2mg/mL stock)

4 Wash iPSCs with PBS, add 500ul of Accutase per well and incubate at  37 °C for  00:05:00 until cells have dissociated. Gently tap the plate if they are not detaching

5 After cells have dissociated, lift cells by adding 1mL of Knockout DMEM/F12 in well and pipetting up and down 2-3 times. Add cell suspension to a 15mL Falcon tube containing  4 mL of Knockout DMEM/F12



- 6 Spin cell suspensions at 300xG for 00:05:00
- 7 Aspirate the supernatant and resuspend the cells in 2 mL of **Pre-differentiation media**. Pipette gently up and down to ensure a single-cell suspension.
- 8 Count cells and prepare cell solution in **pre-differentiation media**. Plate 500k cells per well of a 6-well plate or 3-5 million cells in a single T75 flask.

Days 1 - 2

- 9 **Day 1:** remove media and add fresh **pre-differentiation media** (without Rock inhibitors from this point)
- 10 **Day 2:** remove media and add fresh **pre-differentiation media**
- 10.1 **Day 2:** Prepare plates for final plating. Coat plates with PLO (100ug/mL) in Borate Buffer. Incubate plates overnight at 37 °C

Day 3

10m







- 11 Finish coating the plates: Wash 2x with ddH₂O. Leave to dry under the hood for 00:10:00 and coat with 100ug/mL Laminin in DMEM. Incubate at 37 °C for at least 02:00:00
- 12 At this point the cells should appear more neuronal, with neurites. If cells are healthy, prepare for plating or freeze down.
- 13 Prepare **neuronal media** supplemented with **doxycyclin** (100ug/mL) and 10 micromolar (μM) **Rock inhibitors**

10m



**Note**

How to prepare 50mL **neuronal media**:

-  23.4 mL DMEM/F12 (Thermo Fisher)
-  23.4 mL Neurobasal-A (Thermo Fisher)
-  500 μ L Non-essential amino acids (Thermo Fisher)
-  250 μ L 100X Glutamax (Thermo Fisher)
-  250 μ L 100X N2 supplement (Thermo Fisher)
-  500 μ L 50X B27 supplement (Thermo Fisher)
-  50 μ L BDNF (From 10ug/mL stock)
-  50 μ L NT-3 (From 10ug/mL stock)
-  50 μ L Laminin (From 1mg/mL stock)

- 14 Wash cells with PBS and add Accutase ( 500 μ L in 6-well plate or  4 mL in T75 flask), incubate at  37 °C for  00:05:00 until cells have dissociated
- 15 After cells have dissociated, lift cells by adding 1mL of DMEM+Neurobasal solution in well and pipetting up and down 2-3 times. Add cell suspension to a 15mL Falcon tube containing  4 mL of DMEM+Neurobasal
- 16 Spin cell suspensions at 300xG for  00:05:00
- 17 Plate cells into PLO-coated plates in **neuronal media** supplemented with **Doxycycline and Rock inhibitors**

Note

When plating, prepare a cell suspension of roughly 500k cells/mL of neuronal media. Depending on the platform used for the desired experiments ( 140 μ L for 96-well plates,  1 mL of 12-well plates...)



Day 4 - onwards

- 18 **Day 7:** remove half of the media and add equal volume of fresh **neuronal media** (without Rock inhibitors and Doxycyclin from this point)

Note

The cells tend to become **very peely** from this point, proceed with care when handling

- 19 **Day 14:** remove half of the media and add equal volume of fresh **neuronal media**

Note

Cells are considered mature and are ready for experiments from Day 14.

- 20 **Day 21:** remove half of the media and add equal volume of fresh **neuronal media**

- 21 **Day 28:** remove half of the media and add equal volume of fresh **neuronal media**

Protocol references

Chen et al., 2020 dx.doi.org/10.17504/protocols.io.bcrjiv4n