



Feb 18, 2025

ICC Staining of Fixed Organoid Cells for Flow Analysis

DOI

<https://dx.doi.org/10.17504/protocols.io.bp2l6de6dvqe/v1>

Katharina Meyer¹

¹Wyss Institute



Katharina Meyer

Wyss Institute

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bp2l6de6dvqe/v1>

Protocol Citation: Katharina Meyer 2025. ICC Staining of Fixed Organoid Cells for Flow Analysis . **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bp2l6de6dvqe/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 18, 2025



Last Modified: February 18, 2025

Protocol Integer ID: 120848

Keywords: pfa fixed dissociated organoid cells for flow analysis, organoid cells for flow analysis, pfa fixed dissociated organoid cell, organoid cell, cell, flow analysis, indirect icc, staining

Abstract

Indirect ICC Staining of PFA Fixed Dissociated Organoid Cells for Flow Analysis



Dissociation of Organoids

- 1 Wash organoids in each well with PBS 2 times.
- 2 Transfer organoids into a 2-ml Eppendorf tube containing 300-500 μ l of cold Cell Recovery Solution for Matrigel embedded organoids (Corning).
- 3 Incubate on ice for 1 hour.
- 4 Spin down at 1000g for 1 minute.
- 5 Remove supernatant and wash with 1 ml cold PBS once.
- 6 Incubate organoids with 300-500 μ l cold Accutase for 15 minutes.
- 7 Spin down at 1000g for 1 minute.
- 8 Remove supernatant and wash with cold PBS once.
- 9 Dissociate and resuspend the cells in 500 μ l cold 4% PFA to start fixation.
- 10 Incubate on ice for 1 hour.
- 11 Spin down at 3000g for 1 minute.
- 12 Remove supernatant and wash with 3 $\times$ 1 ml PBS. Resuspend pellet each time by vortexing.
- 13 Keep cells in PBS until permeabilization.

ICC Staining

- 14 Recover cells by spinning down at 3000g for 1 minute.
- 15 Count cells.
- 16 Remove supernatant and resuspend the pellet with 300 μ l blocking buffer (3% Normal donkey serum and 0.3% Triton X-100 in PBS).
- 17 Incubate for 30 minutes and briefly vortex the cells at 30% speed.
- 18 Spin down at 3000g for 1 minute.
- 19 Remove supernatant and resuspend the cells with blocking buffer + primary antibody/ies (100-250 μ l and minimum of 100K cells) via vortexing at 70% speed.
- 20 Lay all tubes on a rack holder and secure them with tape.
- 21 Incubate on a shaker at 4°C overnight.
- 22 Wash the cells 3 $\times$ 500 μ l PBS. Spin down at 3000g for 1 minute.
- 23 Resuspend the cells with 100-250 μ l of blocking buffer + secondary antibody/ies.
- 24 Incubate the cells at room temperature for 1 hour and briefly vortex the cells at 30% speed.
- 25 Wash the cells 3 $\times$ 500 μ l PBS. Spin down at 3000g for 1 minute.



- 26 Resuspend the cells in 0.5-1 ml PBS and filter the cells through a mesh of 0.45 μm . Keep the cells from light. On ice.