

Oct 30, 2025

IDEEL- UNC Implementation of DBLa Var-Typing

DOI

<https://dx.doi.org/10.17504/protocols.io.q26g7n781lwz/v1>



Infectious Disease Epidemiology and Ecology Lab¹

¹Infectious Disease Epidemiology and Ecology Lab, University of North Carolina at Chapel Hill, Chapel Hill, NC, USA

IDEEL



Infectious Disease Epidemiology and Ecology Lab

University of North Carolina

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.q26g7n781lwz/v1>

Protocol Citation: Infectious Disease Epidemiology and Ecology Lab 2025. IDEEL- UNC Implementation of DBLa Var-Typing. protocols.io <https://dx.doi.org/10.17504/protocols.io.q26g7n781lwz/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 27, 2025

Last Modified: October 30, 2025

Protocol Integer ID: 231055

Keywords: DBLa PCR, Amplification, Sequencing, Primers, unc implementation of dbla var, dbla var, continued plasmodium falciparum disease transmission, pcr assay, plasmodium falciparum disease transmission after an outbreak, use with the current dbla, molecular epidemiology, existing primer sequence, current dbla, pcr, primer sequence, sequencing approach, adapted primer, reverse primer, forward primer, standardized oligo segment, example publication of this protocol

Abstract

This protocol is designed for use with the current DBLa sequencing approach developed by the Day laboratory. An example publication of this protocol is: Ruybal-Pesántez S, Sáenz FE, Deed SL, Johnson EK, Larremore DB, Vera-Arias CA, Tiedje KE, Day KP. Molecular epidemiology of continued Plasmodium falciparum disease transmission after an outbreak in Ecuador. Front Trop Dis. 2023;4:1085862. doi: 10.3389/fitd.2023.1085862. Epub 2023 Mar 16. PMID: 39525803; PMCID: PMC11546077.

PCR assays are “Nextera-Adapted” by adding standardized oligo segments onto existing primer sequences.

For example:

- Non-Nextera-Adapted Primer:

1. Forward primer: CCATCAGGGAAATGTCCAGT
2. Reverse primer : TTTCCTGCATGTCTTGAACA
3. Forward primer addition: **TCGTCGGCAGCGTCAGATGTGTATAAGAGACAG**
4. Reverse primer addition: **GTCTCGTGGGCTCGGAGATGTGTATAAGAGACAG**

Example:

- Nextera-Adapted Primer:

1. Forward: **TCGTCGGCAGCGTCAGATGTGTATAAGAGACAG**CCATCAGGGAAATGTCCAGT
2. Reverse: **GTCTCGTGGGCTCGGAGATGTGTATAAGAGACAG**TTTCCTGCATGTCTTGAACA

PCR-based library preparation allows for a significant reduction in time to complete libraries, reduces the consumables required, and lessens the impact of library preparation on technicians. Pre-combined 96-well Unique Dual Index plates are utilized in combination with standardized barcoding worksheets to simplify the generation of final library indexing manifests.

Materials

Materials

Consumables:

- P200 pipettor and tips/ single or multichannel
- P20 pipettor and tips/ single or multichannel
- 25mL reagent reservoir
- 96-well PCR plates
- Nuclease-free PCR microtubes/strips/plates.
- Nuclease-free 1.5µL Eppendorf tubes (for master mixes).
- Cold blocks for microplates and 1.5µL tubes.

Reagents:

- 80% Ethanol
- Qubit or PicoGreen supplies for quantifying product
- Speedbeed mix (dx.doi.org/10.17504/protocols.io.x54v9p7b1g3e/v2)
- Roche Fast Start High Fidelity PCR System (Cat #: 03553400001)
- Molecular Grade Water (MGW)
- Pre-Combined UDI 96-well plate at 10µM each index (i7 and i5 Index primers)
- DBLa primers
- KAPA HiFi HotStart ReadyMix 2X (Cat #: 07958927001)

Materials: DBLa Primers

	A	B
	Primer	Sequence
	Nextera DBLaF	TCG TCG GCA GC GTC AGA TGT GTA TAA GAG ACA GG CAC GM AGT

	A	B
		TTY GC
	NexteraDBLaR	GTC TCG TGG GCT CG GAG ATG TGT ATA AGA GAC AGG CCC ATT CST CGA ACC A

Materials: I7 Index primers

	A	B	C	D
	Beginn ing	Index	End	To order
	CAA GCA GAA GAC GG CAT ACG AGA T	CGC TCA GTT C	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG AGA TCG CTC AGT TCG TCT CGT GG GCT CG G
	CAA GCA GAA GAC GG CAT ACG	TAT CTG ACC T	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG



	A	B	C	D
	AGA T			AGA TTA TCT GAC CTG TCT CGT GG GCT CG G
	CAA GCA GAA GAC GG CAT ACG AGA T	TCG GAT GTC G	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG AGA TTC GGA TGT CG GTC TCG TGG GCT CG G
	CAA GCA GAA GAC GG CAT ACG AGA T	CTT ATG GAA T	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG AGA TCT TAT GGA ATG TCT CGT GG GCT CG G
	CAA GCA GAA GAC GG CAT ACG AGA T	TCC TAT TGT G	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG AGA TTC CTA TTG



	A	B	C	D
				TGG TCT CGT GG GCT CG G
	CAA GCA GAA GAC GG CAT ACG AGA T	GC GC GAT GTT	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG AGA TGC GC GAT GTT GTC TCG TGG GCT CG G
	CAA GCA GAA GAC GG CAT ACG AGA T	AGA GCA CTA G	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG AGA TAG AGC ACT AGG TCT CGT GG GCT CG G
	CAA GCA GAA GAC GG CAT ACG AGA T	TGC CTT GAT C	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG AGA TTG CCT TGA TCG TCT CGT GG

	A	B	C	D
				GCT CG G
	CAA GCA GAA GAC GG CAT ACG AGA T	CTA CTC AGT C	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG AGA TCT ACT CAG TCG TCT CGT GG GCT CG G
	CAA GCA GAA GAC GG CAT ACG AGA T	TCG TCT GAC T	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG AGA TTC GTC TGA CTG TCT CGT GG GCT CG G
	CAA GCA GAA GAC GG CAT ACG AGA T	GAA CAT ACG G	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG AGA TGA ACA TAC GG GTC TCG TGG GCT CG G



	A	B	C	D
	CAA GCA GAA GAC GG CAT ACG AGA T	CCT ATG ACT C	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG AGA TCC TAT GAC TCG TCT CGT GG GCT CG G
	CAA GCA GAA GAC GG CAT ACG AGA T	TAA TGG CAA G	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG AGA TTA ATG GCA AGG TCT CGT GG GCT CG G
	CAA GCA GAA GAC GG CAT ACG AGA T	GTG CCG CTT C	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG AGA TGT GCC GCT TCG TCT CGT GG GCT CG G
	CAA GCA GAA	CG GCA	GTC TCG TGG	CAA GCA GAA



	A	B	C	D
	GAC GG CAT ACG AGA T	ATG GA	GCT CG G	GAC GG CAT ACG AGA TCG GCA ATG GAG TCT CGT GG GCT CG G
	CAA GCA GAA GAC GG CAT ACG AGA T	GCC GTA ACC G	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG AGA TGC CGT AAC CG GTC TCG TGG GCT CG G
	CAA GCA GAA GAC GG CAT ACG AGA T	AAC CAT TCT C	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG AGA TAA CCA TTC TCG TCT CGT GG GCT CG G
	CAA GCA GAA GAC GG CAT ACG	GGT TGC CTC T	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG

	A	B	C	D
	AGA T			AGA TGG TTG CCT CTG TCT CGT GG GCT CG G
	CAA GCA GAA GAC GG CAT ACG AGA T	CTA ATG ATG G	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG AGA TCT AAT GAT GG GTC TCG TGG GCT CG G
	CAA GCA GAA GAC GG CAT ACG AGA T	TCG GCC TAT C	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG AGA TTC GG CCT ATC GTC TCG TGG GCT CG G
	CAA GCA GAA GAC GG CAT ACG AGA T	AGT CAA CCA T	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG AGA TAG TCA ACC

	A	B	C	D
				ATG TCT CGT GG GCT CG G
	CAA GCA GAA GAC GG CAT ACG AGA T	GAG CGC AAT A	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG AGA TGA GC GCA ATA GTC TCG TGG GCT CG G
	CAA GCA GAA GAC GG CAT ACG AGA T	AAC AAG GC GT	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG AGA TAA CAA GG CGT GTC TCG TGG GCT CG G
	CAA GCA GAA GAC GG CAT ACG AGA T	GTA TGT AGA A	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG AGA TGT ATG TAG AAG TCT CGT GG

	A	B	C	D
				GCT CG G
	CAA GCA GAA GAC GG CAT ACG AGA T	TTC TAT GGT T	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG AGA TTT CTA TGG TTG TCT CGT GG GCT CG G
	CAA GCA GAA GAC GG CAT ACG AGA T	CCT CGC AAC C	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG AGA TCC TCG CAA CCG TCT CGT GG GCT CG G
	CAA GCA GAA GAC GG CAT ACG AGA T	TGG ATG CTT A	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG AGA TTG GAT GCT TAG TCT CGT GG GCT CG G



	A	B	C	D
	CAA GCA GAA GAC GG CAT ACG AGA T	ATG TCG TGG T	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG AGA TAT GTC GTG GTG TCT CGT GG GCT CG G
	CAA GCA GAA GAC GG CAT ACG AGA T	AGA GTG CG GC	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG AGA TAG AGT GC GG CGT CTC GTG GG CTC GG
	CAA GCA GAA GAC GG CAT ACG AGA T	TGC CTG GTG G	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG AGA TTG CCT GGT GG GTC TCG TGG GCT CG G
	CAA GCA GAA	TGC GTG	GTC TCG TGG	CAA GCA GAA



	A	B	C	D
	GAC GG CAT ACG AGA T	TCA C	GCT CG G	GAC GG CAT ACG AGA TTG CGT GTC ACG TCT CGT GG GCT CG G
	CAA GCA GAA GAC GG CAT ACG AGA T	CAT ACA CTG T	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG AGA TCA TAC ACT GTG TCT CGT GG GCT CG G
	CAA GCA GAA GAC GG CAT ACG AGA T	CGT ATA ATC A	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG AGA TCG TAT AAT CAG TCT CGT GG GCT CG G
	CAA GCA GAA GAC GG CAT ACG	TAC GC GG CTG	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG



	A	B	C	D
	AGA T			AGA TTA CGC GG CTG GTC TCG TGG GCT CG G
	CAA GCA GAA GAC GG CAT ACG AGA T	GC GAG TTA CC	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG AGA TGC GAG TTA CCG TCT CGT GG GCT CG G
	CAA GCA GAA GAC GG CAT ACG AGA T	TAC GG CCG GT	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG AGA TTA CG GCC GGT GTC TCG TGG GCT CG G
	CAA GCA GAA GAC GG CAT ACG AGA T	GTC GAT TAC A	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG AGA TGT CGA TTA



	A	B	C	D
				CAG TCT CGT GG GCT CG G
	CAA GCA GAA GAC GG CAT ACG AGA T	CTG TCT GCA C	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG AGA TCT GTC TGC ACG TCT CGT GG GCT CG G
	CAA GCA GAA GAC GG CAT ACG AGA T	CAG CCG ATT G	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG AGA TCA GCC GAT TGG TCT CGT GG GCT CG G
	CAA GCA GAA GAC GG CAT ACG AGA T	TGA CTA CAT A	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG AGA TTG ACT ACA TAG TCT CGT GG

	A	B	C	D
				GCT CG G
	CAA GCA GAA GAC GG CAT ACG AGA T	ATT GCC GAG T	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG AGA TAT TGC CGA GTG TCT CGT GG GCT CG G
	CAA GCA GAA GAC GG CAT ACG AGA T	GCC ATT AGA C	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG AGA TGC CAT TAG ACG TCT CGT GG GCT CG G
	CAA GCA GAA GAC GG CAT ACG AGA T	GG CGA GAT GG	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG AGA TGG CGA GAT GG GTC TCG TGG GCT CG G



	A	B	C	D
	CAA GCA GAA GAC GG CAT ACG AGA T	TGG CTC GCA G	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG AGA TTG GCT CGC AGG TCT CGT GG GCT CG G
	CAA GCA GAA GAC GG CAT ACG AGA T	TAG AAT AAC G	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG AGA TTA GAA TAA CG GTC TCG TGG GCT CG G
	CAA GCA GAA GAC GG CAT ACG AGA T	TCC ATG TTG C	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG AGA TTC CAT GTT GC GTC TCG TGG GCT CG G
	CAA GCA GAA	TAT CCA	GTC TCG TGG	CAA GCA GAA

	A	B	C	D
	GAC GG CAT ACG AGA T	GGA C	GCT CG G	GAC GG CAT ACG AGA TTA TCC AGG ACG TCT CGT GG GCT CG G
	CAA GCA GAA GAC GG CAT ACG AGA T	AGT GCC ACT G	GTC TCG TGG GCT CG G	CAA GCA GAA GAC GG CAT ACG AGA TAG TGC CAC TGG TCT CGT GG GCT CG G

Materials: 15 Index Primers

	A	B	C	D
	Beginn ing	Index	End	To order
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	TCG TGG AGC G	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC TCG TGG AGC GTC



	A	B	C	D
				GTC GG CAG CGT C
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	CTA CAA GAT A	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC CTA CAA GAT ATC GTC GG CAG CGT C
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	TAC GTT CAT T	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC TAC GTT CAT TTC GTC GG CAG CGT C
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	TGC CTG GTG G	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC TGC CTG GTG GTC GTC



	A	B	C	D
				GG CAG CGT C
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	TCC ATC CGA G	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC TCC ATC CGA GTC GTC GG CAG CGT C
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	GTC CAC TTG T	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC GTC CAC TTG TTC GTC GG CAG CGT C
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	TGG AAC AGT A	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC TGG AAC AGT ATC GTC GG



	A	B	C	D
				CAG CGT C
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	CCT TGT TAA T	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC CCT TGT TAA TTC GTC GG CAG CGT C
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	GTT GAT AGT G	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC GTT GAT AGT GTC GTC GG CAG CGT C
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	ACC AGC GAC A	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC ACC AGC GAC ATC GTC GG CAG



	A	B	C	D
				CGT C
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	CAT ACA CTG T	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC CAT ACA CTG TTC GTC GG CAG CGT C
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	GTG TGG CGC T	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC GTG TGG CGC TTC GTC GG CAG CGT C
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	ATC ACG AAG G	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC ATC ACG AAG GTC GTC GG CAG



	A	B	C	D
				CGT C
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	CG GCT CTA CT	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC CG GCT CTA CTT CGT CG GCA GC GTC
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	GAA TGC ACG A	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC GAA TGC ACG ATC GTC GG CAG CGT C
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	AAG ACT ATA G	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC AAG ACT ATA GTC GTC GG CAG



	A	B	C	D
				CGT C
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	TCG GCA GCA A	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC TCG GCA GCA ATC GTC GG CAG CGT C
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	CTA ATG ATG G	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC CTA ATG ATG GTC GTC GG CAG CGT C
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	GGT TGC CTC T	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC GGT TGC CTC TTC GTC GG CAG



	A	B	C	D
				CGT C
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	CGC ACA TGG C	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC CGC ACA TGG CTC GTC GG CAG CGT C
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	GG CCT GTC CT	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC GG CCT GTC CTT CGT CG GCA GC GTC
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	CTG TGT TAG G	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC CTG TGT TAG GTC GG CAG



	A	B	C	D
				CGT C
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	TAA GGA ACG T	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC TAA GGA ACG TTC GTC GG CAG CGT C
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	CTA ACT GTA A	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC CTA ACT GTA ATC GTC GG CAG CGT C
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	GG CGA GAT GG	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC GG CGA GAT GGT CGT CG GCA



	A	B	C	D
				GC GTC
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	AAT AGA GCA A	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC AAT AGA GCA ATC GTC GG CAG CGT C
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	TCA ATC CAT T	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC TCA ATC CAT TTC GTC GG CAG CGT C
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	TCG TAT GC GG	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC TCG TAT GC GGT CGT CG GCA

	A	B	C	D
				GC GTC
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	TCC GAC CTC G	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC TCC GAC CTC GTC GTC GG CAG CGT C
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	CTT ATG GAA T	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC CTT ATG GAA TTC GTC GG CAG CGT C
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	GCT TAC GGA C	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC GCT TAC GGA CTC GTC GG CAG



	A	B	C	D
				CGT C
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	GAA CAT ACG G	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC GAA CAT ACG GTC GTC GG CAG CGT C
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	GTC GAT TAC A	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC GTC GAT TAC ATC GTC GG CAG CGT C
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	ACT AGC CGT G	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC ACT AGC CGT GTC GTC GG CAG



	A	B	C	D
				CGT C
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	AAG TTG GTG A	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC AAG TTG GTG ATC GTC GG CAG CGT C
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	TGG CAA TAT T	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC TGG CAA TAT TTC GTC GG CAG CGT C
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	GAT CAC CGC G	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC GAT CAC CGC GTC GG CAG



	A	B	C	D
				CGT C
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	TAC CAT CCG T	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC TAC CAT CCG TTC GTC GG CAG CGT C
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	GCT GTA GGA A	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC GCT GTA GGA ATC GTC GG CAG CGT C
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	CGC ACT AAT G	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC CGC ACT AAT GTC GTC GG CAG



	A	B	C	D
				CGT C
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	GAC AAC TGA A	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC GAC AAC TGA ATC GTC GG CAG CGT C
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	AGT GGT CAG G	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC AGT GGT CAG GTC GTC GG CAG CGT C
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	CCA TCT CGC C	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC CCA TCT CGC CTC GTC GG CAG



	A	B	C	D
				CGT C
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	CTG CGA GCC A	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC CTG CGA GCC ATC GTC GG CAG CGT C
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	CGT TAT TCT A	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC CGT TAT TCT ATC GTC GG CAG CGT C
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	GCA ACA TGG A	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC GCA ACA TGG ATC GTC GG CAG

	A	B	C	D
				CGT C
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	GTC CTG GAT A	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC GTC CTG GAT ATC GTC GG CAG CGT C
	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC	CAG TGG CAC T	TCG TCG GCA GC GTC	AAT GAT ACG GC GAC CAC CGA GAT CTA CAC CAG TGG CAC TTC GTC GG CAG CGT C

Step 1: General Lab measures

- 1 Ensure all work areas and pipettors are sterilized per laboratory protocol.

Note

- Clean pipettes and working bench tops with 70% ethanol.
- UV-treat the working surface, all equipment, and tubes/strips.

- 2 DBLa PCR (PCR1) should be set up in a clean room, free of PCR amplification.
- 3 Indexing PCR (PCR2), clean-up, and quantification should be done in an area where PCR products can be worked with (i.e., a post-PCR space).
- 4 Thaw reagents and DNA samples  On ice .
- 5 Briefly vortex and centrifuge each reagent before use, with the exception of enzymes. 

Note

Gently mix enzymes by inverting the tube a few times, then centrifuge.

Step 2: Setting up Primers

- 6 Dilute the  100 micromolar (μM) stock primers down to  20 micromolar (μM) working solutions for DBLa forward and reverse primers.

- 7 Generation of UDI Plates

- 7.1 The final concentration of the i5/i7 primer mix should be  10 micromolar (μM) . This can be done any number of ways, including diluting the  100 micromolar (μM) stock 

of each primer to [M] 10 micromolar (μM) and then mixing.

- 7.2 UDI Combinatorial Plates can be made ahead of time and stored at $-20\text{ }^{\circ}\text{C}$ for use.
- 7.3 Typically, we recommend making plates for a maximum of 10 indexing PCR reactions; thus, if contamination should occur, a new plate can be made without wasting reagents.




Figure legend: Generation of UDI (Unique Dual Index) plates. By ordering i5 (yellow) and i7 (blue) on plates in opposite orientation, UDI plates of combined i5 and i7 indexes (yellow + blue = green) can be made using a multichannel pipette to mix primers in a LoBind DNA plate. The column of i5 primers should be pipetted across all columns of the UDI plate and the i7 primer row pipetted across all rows of the plate.


Step 3: Initial DBLa PCR

- 8 Calculate the amount of reagents needed with a 10% excess to account for pipette loss.
- 9 Prepare the PCR master mix based on the table below.

A	B
Reagent	Volume/rxn (μL)
Roche Fast Start High Fidelity PCR System Buffer (2X)	2.5
Nucleotide Mix	0.5
Fwd Primer (20uM)	2.5
Rev Primer (20uM)	2.5
Enzyme Blend	0.2
Water	14.8

Note

- All reagents, PCR master mix, primer pools, and samples should be held in freezer blocks (tube or 96-well) unless stated otherwise.
- All assays are carried out with positive and negative controls (preferably 2 of each). For a positive control, we use a mocked dried blood spot containing an estimated parasite density of 100 Pf/μL. The positive control you use should be extracted the same way as the samples to be tested. For non-template control, we recommend using human DNA at a clinically relevant concentration (e.g. 16.77 ng/μL). If human DNA is not available, water can be used as a non-template control.

- 10 Combine the master mix and  2 μL sample DNA (25μl total reaction volume), seal caps/strips/plates, and place the reaction in a thermal cycler using the following conditions:



A	B	C	D
Step	Duration	Cycles	Temperature
Activation	5 min	1 cycle	95° C
	1 min	1 cycle	42° C
	1 min	1 cycle	60° C
Cycling	1 min	29 cycles	95° C
	1 min		42° C
	1 min		60° C



Step 4: PCR cleanup and library prep

45m

11 Bring Speedbeads to  Room temperature for  00:30:00 .

30m

12 Mix fresh 80% Ethanol using 200 proof ethanol (100%, anhydrous) and molecular grade water.



13 Vortex Speed-Beads thoroughly, then add  50 μ L to each well. Mix gently by pipetting.



14 Allow the mixture to incubate at  Room temperature for  00:10:00 .

10m



Note

If there are droplets on the sides of the wells following incubation, cover the plate with a foil seal and briefly spin the plate to avoid losing product.

15 Place the plate on a 96-well magnetic block and wait  00:02:00 -  00:03:00 for the solution to clear.

3m

16 Aspirate the clear supernatant and discard.

17 Dispense  120 μ L 80% Ethanol into each well, incubate  00:00:30 , then aspirate to wash.

30s



18 Repeat for a total of 2 washes.



19 Aspirate final drops of 80% Ethanol with P20 pipette tips. Allow bead pellets to dry momentarily for  00:00:30 to  00:01:30 .

1m 30s



20 Remove plate from the magnet.

- 21 Dispense  25 µL of Elution Buffer onto each bead pellet. Mix by pipetting to resuspend DNA. 
- 22 Cover plate with foil seal and briefly centrifuge to collect product in the bottom of wells before storing or proceeding to library preparation. *The beads can remain in the eluent and may be added to the library preparation reaction. The eluent is not transferred to a new plate to limit cross-contamination of concentrated amplicon products.* 

Step 5: Library Preparation (PCR2)

- 23 Thaw Pre-Combined UDI 96-well plate, purified PCR reaction, and KAPA Master Mix.
- 24 Calculate the amount of reagents needed with a 10% excess to account for pipette loss.
- 25 Mix master mix as described in the following table: 

A	B
Reagent	Volume/rxn
KAPA HiFi HotStart ReadyMix 2X	12.5 µL
Molecular Grade Water	4.5 µL
Total	17.0 µL

- 26 Combine KAPA HiFi HotStart ReadyMix 2X and water in a tube, then aliquot equal volumes into 0.2mL microtube strips. Pipette  17.0 µL Master Mix into each well with a multichannel pipettor. 
- 27 Spin UDI plate to ensure there is no liquid on the plate sealing foil.
- 28 Remove the foil seal carefully to avoid contamination.

29 Add  2.5 µL of the appropriate UDI to each well. No need to mix yet.



30 Seal stock UDI plate before opening and adding samples to prevent contamination.

31 Add  5.0 µL of PCR1. Gently mix by pipetting.



- Sample and UDI plate should now have the following volumes in each well:

A	B
Reagent	Volume/rxn
KAPA HiFi HotStart ReadyMix 2X	12.5 µL
Molecular Grade Water	4.5 µL
UDI Primer	2.5 µL
Purified amplicons	5.0 µL
Total	24.5 µL

32 Program the thermal cycler to the following conditions for 8-12 cycles.

A	B	C	D
Step	Duration	Cycles	Temperature
Initial Activation Step	3 min	1	95° C
Denaturation	30 sec	8-12	95° C
Annealing	30 sec		55° C
Extension	30 sec		72° C
Final Extension	5 min	1	72° C
Hold	Indefinite	1	4° C

Step 6: Contamination check



- 33 Prior to proceeding to sample pooling into a final library, the libraries generated from positive and non-template controls should be checked for contamination.
- 34 It is possible that negative controls will have a detectable band on gel below the bands consistent with amplification products. Non-template control bands typically represent primer dimer that should not yield data when processed through bioinformatic analysis.

Step 7: Pooling based on PicoGreen Concentration

54m 30s

- 35 Bring SpeedBeads to Room temperature ~ 00:30:00 . 30m
- 36 Mix fresh 80% Ethanol using 200 proof ethanol (100%, anhydrous) and molecular grade water.
- 37 Vortex Speed-Beads thoroughly before use.
- 38 Combine libraries with Speed-Beads at a 1:1 ratio.
- 39 Allow the mixture to incubate at Room temperature for 00:10:00 . 10m
- 40 Gently spin the plate to collect contents at the bottom of the wells.
- 41 Place the plate on a plate magnet and wait 00:02:00 - 00:03:00 for the solution to clear. 3m
- 42 Aspirate the clear supernatant and discard.
- 43 Keeping the plate on the magnet, dispense 80% Ethanol into the wells to wash. Ethanol wash volume should be higher than the combined library and Speed-Bead mixture.
- 44 Incubate 00:00:30 , then aspirate to wash. 30s



45 Repeat for a total of 2 washes.



46 Aspirate final drops of ethanol with P20 pipette tips. Allow bead pellets to dry on the magnet for 00:00:30 to 00:01:00 .

1m

Note

Dry beads to a satin-matte finish, but not to the point of cracking. This happens quickly if ethanol has been aspirated completely.

47 Remove the plate from the magnet and dispense 25 μL Elution Buffer onto each bead pellet. Mix by pipetting to resuspend DNA.



48 Allow the mixture to incubate at Room temperature for 00:05:00 .

5m



49 Gently spin the plate to collect product in the bottom of wells, then place on a plate magnet and wait 00:02:00 - 00:03:00

3m



50 Remove eluent to a new LoBind plate.

51 Quantify 2 μL of amplified libraries in 48 μL of 1X PicoGreen reagent in a black, flat bottom microplate appropriate for reading fluorescence.

Note

- If the reagent is not already provided at 1X, combine dye and buffer in a 1:200 ratio to create a working solution.
- If you are quantifying an entire 96-well plate of samples, a second plate will be necessary for assay standards.

52 Add 2 μL of each sample to a black, flat-bottomed polystyrene assay plate.



- The 90-degree angle where the well wall meets the flat bottom holds this droplet nicely.
- Add  2 μL of each standard to the appropriate wells.
- Standards 1-8 in duplicate = 16 wells required.

53 Fill each sample and standard wells with  48 μL 1X PicoGreen Reagent. 

Note

Pipette to mix thoroughly, careful not to introduce bubbles.

54 Cover plate with optical film and incubate for  00:02:00 at  Room temperature before assaying. Protect from light during incubation. 

2m

55 Place the prepared quantitation plate into the multimode plate reader.

Note

The plate reader should be set up to analyze Fluorescence Intensity (FI) with an excitation wavelength around 480nm and an emission wavelength around 580nm.

56 Calculate each sample's concentration by creating a standard curve and plotting unknown fluorescence values. Divide the resulting value by two to report nanograms per microliter.

57 Negative controls can have detectable DNA by PicoGreen, typically in the form of dimers. To determine the concentration of the library in order to pool equimolar volumes of each sample, one should use a baseline (averaged NTC concentration) subtracted from the concentration for each sample to determine pooling volumes.

58 The library pool is ready to prepare for sequencing.