

Jan 14, 2026

Version 1

IDEEL- UNC Implementation of Targeted csp Sequencing- MIP to Nanopore V.1

DOI

<https://dx.doi.org/10.17504/protocols.io.8epv55xz5v1b/v1>

Jonathan Juliano¹

¹Infectious Disease Epidemiology and Ecology Lab, University of North Carolina

IDEEL



Infectious Disease Epidemiology and Ecology Lab

University of North Carolina

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.8epv55xz5v1b/v1>

Protocol Citation: Jonathan Juliano 2026. IDEEL- UNC Implementation of Targeted csp Sequencing- MIP to Nanopore. protocols.io <https://dx.doi.org/10.17504/protocols.io.8epv55xz5v1b/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: November 10, 2025

Last Modified: January 14, 2026

Protocol Integer ID: 238553


Keywords: PCR2 primer preparation, Pooling and Cleanup, Nanopore library preparation, Adaptor ligation, Nanopore sequencing, nanopore this protocol detail, amplicon generation through pcr2 primer preparation, nanopore library preparation, pcr2 primer preparation, nanopore, sequencing, amplicon generation, derived amplicon generation, unc implementation, mip, csp

Abstract

This protocol details the IDEEL-UNC targeted csp sequencing workflow, from MIP-derived amplicon generation through PCR2 primer preparation, pooling, cleanup, Nanopore library preparation, adaptor ligation, and sequencing.

Guidelines

Location of amplicon relative to *csp*:



Alignment of reverse primer



XM_001351086.1 reverse	ATGATGAGAAAATTAGCTATTTTATCTGTTTCTTCCTTTTATTTGTTGAGGCCTTATTC -----CCTTATTC *****
XM_001351086.1 reverse	CAGGAATACCACTGCTATGGAAGTTCGTCAACACAAGGGTTCTAAATGAATTAATTAT CAGGAATACCACTGC----- *****
XM_001351086.1 reverse	GATAATGCAGGCACTAATTTATATAATGAATTAGAAATGAATTATTATGGGAAACAGGAA -----
XM_001351086.1 reverse	AATTGGTATAGTCTTAAAAAATAGTAGATCACTTGAGAGAAATGATGATGGAAATAAC -----
XM_001351086.1 reverse	GAAGACAACGAGAAATTAAGGAAACCAAAACATAAAAAATTAAAGCAACCAGCGGATGGT -----
XM_001351086.1 reverse	AATCCTGATCCAAATGCAAACCCAAATGTAGATCCAATGCCAACCCAAATGTAGATCCA -----
XM_001351086.1 reverse	AATGCAAACCCAAATGTAGATCCAAATGCAAACCCAAATGCAAACCCAAATGCAAACCCA -----
XM_001351086.1 reverse	AATGCAAACCCAAATGCAAACCCAAATGCAAACCCAAATGCAAACCCAAATGCAAACCCA -----
XM_001351086.1 reverse	AATGCAAACCCAAATGCAAACCCAAATGCAAACCCAAATGCAAACCCAAATGCAAACCCA -----
XM_001351086.1 reverse	AATGCAAACCCCAATGCAATCCTAATGCAAACCCCAATGCAAACCCCAACGTAGATCCT -----
XM_001351086.1 reverse	AATGCAATCCAAATGCAAACCCCAACGCAACCCCAATGCAATCCTAATGCAAACCCC -----
XM_001351086.1 reverse	AATGCAATCCTAATGCAATCCTAATGCCAATCCAAATGCAATCCAAATGCAAACCCA -----
XM_001351086.1 reverse	AACGCAAACCCCAATGCAATCCTAATGCCAATCCAAATGCAATCCAAATGCAAACCCA -----
XM_001351086.1 reverse	AATGCAAACCCCAATGCAAACCCCAATGCAATCCTAATAAAACAATCAAGGTAATGGA -----
XM_001351086.1 reverse	CAAGGTCACAATATGCCAAATGACCCAAACCGAAATGTAGATGAAATGCTAATGCCAAC -----
XM_001351086.1 reverse	AGTGCTGTAAAAATAATAATAACGAAGAACCAAGTGATAAGCACATAAAGAATATTTA -----
XM_001351086.1 reverse	AACAAAATACAAAATTCTTTCAACTGAATGGTCCCCATGTAGTGAACCTGTGGAAAT -----
XM_001351086.1 reverse	GGTATTCAAGTTAGAATAAAGCCTGGCTCTGCTAATAAACCTAAAGACGAATTAGATTAT -----
XM_001351086.1 reverse	GCAAAATGATATTGAAAAAAAAATTTGTAATGGAAAAATGTTCCAGTGTGTTAATGTC -----
XM_001351086.1 reverse	GTAAATAGTTCAATAGGATTAATAATGGTATTATCCTTCTTGTTCTTAATTAG -----

Materials

- Molecular Grade Water (MGW)
- Primers
 - sp-rtss-repeat_v4_F_MIP GACTCGCCAAGCTGAAGNNNTCGCAAACGTAATTAAATATTACAAAA
 - csp-rtss-repeat_v4_R_MIP ACGTGTGCTCTTCCGATCTNNNNCCTTATTCCAGGAATACCAAGTGC
-  GoTaq® Long PCR Master Mix, 2X **Promega Catalog #M402A**
- Illumina FWD and REV barcode primers
- Gel electrophoresis reagents and equipment.
- 1X TE Low EDTA Buffer
- PromethION 2 solo and device specific flow cells
-  Ligation Sequencing Kit V14 **Oxford Nanopore Technologies Catalog #SQK-LSK114**

PCR2 Forward Primers

	A	B	C
	Well Position	Name	Sequence
	A01	hybrid_dual_fw_001	AATGATAC GGCGACCA CCGAGATC TACACGTT AAGACACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	A02	hybrid_dual_fw_002	AATGATAC GGCGACCA CCGAGATC TACACTCTA AGTTACAC TCTTTCCC TACACGAC TCGCCAAG CTGA
	A03	hybrid_dual_fw_003	AATGATAC GGCGACCA CCGAGATC TACACTGTA CATTACACT CTTTCCCT ACACGACT CGCCAAGC TGA
	A04	hybrid_dual_fw_004	AATGATAC GGCGACCA CCGAGATC



	A	B	C
			TACACAATA CCGCACAC TCTTTCCC TACACGAC TCGCCAAG CTGA
	A05	hybrid_dual_ fw_005	AATGATAC GGCGACCA CCGAGATC TACACTGT CTATGACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	A06	hybrid_dual_ fw_006	AATGATAC GGCGACCA CCGAGATC TACACTCA CTGTTACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	A07	hybrid_dual_ fw_007	AATGATAC GGCGACCA CCGAGATC TACACTTG TATTCACAC TCTTTCCC TACACGAC TCGCCAAG CTGA
	A08	hybrid_dual_ fw_008	AATGATAC GGCGACCA CCGAGATC TACACAGA CTAACACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	A09	hybrid_dual_ fw_009	AATGATAC GGCGACCA CCGAGATC TACACATTG TCAAACAC TCTTTCCC TACACGAC TCGCCAAG CTGA
	A10	hybrid_dual_ fw_010	AATGATAC GGCGACCA CCGAGATC



	A	B	C
			TACACGTT CTACAACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	A11	hybrid_dual_ fw_011	AATGATAC GGCGACCA CCGAGATC TACACTTC ATTGACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	A12	hybrid_dual_ fw_012	AATGATAC GGCGACCA CCGAGATC TACACAAG CGTCAACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	B01	hybrid_dual_ fw_013	AATGATAC GGCGACCA CCGAGATC TACACCGT AGGATACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	B02	hybrid_dual_ fw_014	AATGATAC GGCGACCA CCGAGATC TACACAGT TAGCGACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	B03	hybrid_dual_ fw_015	AATGATAC GGCGACCA CCGAGATC TACACGAA CGTCTACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	B04	hybrid_dual_ fw_016	AATGATAC GGCGACCA CCGAGATC



	A	B	C
			TACACAGA ACACTACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	B05	hybrid_dual_ fw_017	AATGATAC GGCGACCA CCGAGATC TAACTGT TGGCAACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	B06	hybrid_dual_ fw_018	AATGATAC GGCGACCA CCGAGATC TACACGGC TAATAACAC TCTTTCCC TACACGAC TCGCCAAG CTGA
	B07	hybrid_dual_ fw_019	AATGATAC GGCGACCA CCGAGATC TACACGCT AGAATACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	B08	hybrid_dual_ fw_020	AATGATAC GGCGACCA CCGAGATC TAACTGTA CTAAACAC TCTTTCCC TACACGAC TCGCCAAG CTGA
	B09	hybrid_dual_ fw_021	AATGATAC GGCGACCA CCGAGATC TACACAGT TACACACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	B10	hybrid_dual_ fw_022	AATGATAC GGCGACCA CCGAGATC



	A	B	C
			TACACGTA CACGAACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	B11	hybrid_dual_ fw_023	AATGATAC GGCGACCA CCGAGATC TACACACT GTAAGACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	B12	hybrid_dual_ fw_024	AATGATAC GGCGACCA CCGAGATC TACACGAT CATGAACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	C01	hybrid_dual_ fw_025	AATGATAC GGCGACCA CCGAGATC TACACTGC TAACGACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	C02	hybrid_dual_ fw_026	AATGATAC GGCGACCA CCGAGATC TACACTAAT GCTCACAC TCTTTCCC TACACGAC TCGCCAAG CTGA
	C03	hybrid_dual_ fw_027	AATGATAC GGCGACCA CCGAGATC TACACCTT GTTGAACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	C04	hybrid_dual_ fw_028	AATGATAC GGCGACCA CCGAGATC



	A	B	C
			TACACGATA CAAGACAC TCTTTCCC TACACGAC TCGCCAAG CTGA
	C05	hybrid_dual_ fw_029	AATGATAC GGCGACCA CCGAGATC TACACATAC AGCGACAC TCTTTCCC TACACGAC TCGCCAAG CTGA
	C06	hybrid_dual_ fw_030	AATGATAC GGCGACCA CCGAGATC TACACCTG GATAGACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	C07	hybrid_dual_ fw_031	AATGATAC GGCGACCA CCGAGATC TACACTCAT CAGCACAC TCTTTCCC TACACGAC TCGCCAAG CTGA
	C08	hybrid_dual_ fw_032	AATGATAC GGCGACCA CCGAGATC TACACTTAT ACCGACAC TCTTTCCC TACACGAC TCGCCAAG CTGA
	C09	hybrid_dual_ fw_033	AATGATAC GGCGACCA CCGAGATC TACACGCT GATTCACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	C10	hybrid_dual_ fw_034	AATGATAC GGCGACCA CCGAGATC



	A	B	C
			TACACATAA GGCCACAC TCTTTCCC TACACGAC TCGCCAAG CTGA
	C11	hybrid_dual_ fw_035	AATGATAC GGCGACCA CCGAGATC TACACAAG ACTAGACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	C12	hybrid_dual_ fw_036	AATGATAC GGCGACCA CCGAGATC TACACTAG CAACAACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	D01	hybrid_dual_ fw_037	AATGATAC GGCGACCA CCGAGATC TACACCAA GTAATACAC TCTTTCCC TACACGAC TCGCCAAG CTGA
	D02	hybrid_dual_ fw_038	AATGATAC GGCGACCA CCGAGATC TACACTCAT GATGACAC TCTTTCCC TACACGAC TCGCCAAG CTGA
	D03	hybrid_dual_ fw_039	AATGATAC GGCGACCA CCGAGATC TACACTAG CTCTTACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	D04	hybrid_dual_ fw_040	AATGATAC GGCGACCA CCGAGATC



	A	B	C
			TACACACT ACATGACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	D05	hybrid_dual_ fw_041	AATGATAC GGCGACCA CCGAGATC TAACTGT GCCATACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	D06	hybrid_dual_ fw_042	AATGATAC GGCGACCA CCGAGATC TACACGAC TTCCAACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	D07	hybrid_dual_ fw_043	AATGATAC GGCGACCA CCGAGATC TAACTATT CTCGACAC TCTTTCCC TACACGAC TCGCCAAG CTGA
	D08	hybrid_dual_ fw_044	AATGATAC GGCGACCA CCGAGATC TACACCAC TAGATACAC TCTTTCCC TACACGAC TCGCCAAG CTGA
	D09	hybrid_dual_ fw_045	AATGATAC GGCGACCA CCGAGATC TACACCTT CGCATACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	D10	hybrid_dual_ fw_046	AATGATAC GGCGACCA CCGAGATC



	A	B	C
			TACACCTA GCGATACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	D11	hybrid_dual_ fw_047	AATGATAC GGCGACCA CCGAGATC TACACTAC GTAAACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	D12	hybrid_dual_ fw_048	AATGATAC GGCGACCA CCGAGATC TACACAAG TAGTCACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	E01	hybrid_dual_ fw_049	AATGATAC GGCGACCA CCGAGATC TACACGAC GTATCACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	E02	hybrid_dual_ fw_050	AATGATAC GGCGACCA CCGAGATC TACACACG TTCAGACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	E03	hybrid_dual_ fw_051	AATGATAC GGCGACCA CCGAGATC TACACAGC TATATACAC TCTTTCCC TACACGAC TCGCCAAG CTGA
	E04	hybrid_dual_ fw_052	AATGATAC GGCGACCA CCGAGATC



	A	B	C
			TACACAGT GCGAAACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	E05	hybrid_dual_ fw_053	AATGATAC GGCGACCA CCGAGATC TACACCAA TGGTAACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	E06	hybrid_dual_ fw_054	AATGATAC GGCGACCA CCGAGATC TACACATTC CGATACAC TCTTTCCC TACACGAC TCGCCAAG CTGA
	E07	hybrid_dual_ fw_055	AATGATAC GGCGACCA CCGAGATC TACACAGC TGTTGACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	E08	hybrid_dual_ fw_056	AATGATAC GGCGACCA CCGAGATC TACACTTC GACTCACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	E09	hybrid_dual_ fw_057	AATGATAC GGCGACCA CCGAGATC TACACCGA ATGGTACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	E10	hybrid_dual_ fw_058	AATGATAC GGCGACCA CCGAGATC



	A	B	C
			TACACTAGT GTACACAC TCTTTCCC TACACGAC TCGCCAAG CTGA
	E11	hybrid_dual_ fw_059	AATGATAC GGCGACCA CCGAGATC TACACACG CATGAACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	E12	hybrid_dual_ fw_060	AATGATAC GGCGACCA CCGAGATC TACACTTC GTACAACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	F01	hybrid_dual_ fw_061	AATGATAC GGCGACCA CCGAGATC TACACTTC TGCTAACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	F02	hybrid_dual_ fw_062	AATGATAC GGCGACCA CCGAGATC TACACTGG AAGCAACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	F03	hybrid_dual_ fw_063	AATGATAC GGCGACCA CCGAGATC TACACGAT GCTCAACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	F04	hybrid_dual_ fw_064	AATGATAC GGCGACCA CCGAGATC



	A	B	C
			TACACAGC CTGAAACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	F05	hybrid_dual_ fw_065	AATGATAC GGCGACCA CCGAGATC TACACCAC CAGTAACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	F06	hybrid_dual_ fw_066	AATGATAC GGCGACCA CCGAGATC TACACGTTA CACCACAC TCTTTCCC TACACGAC TCGCCAAG CTGA
	F07	hybrid_dual_ fw_067	AATGATAC GGCGACCA CCGAGATC TACACGCG TTAGTACAC TCTTTCCC TACACGAC TCGCCAAG CTGA
	F08	hybrid_dual_ fw_068	AATGATAC GGCGACCA CCGAGATC TACACTGA CGTTCACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	F09	hybrid_dual_ fw_069	AATGATAC GGCGACCA CCGAGATC TACACTAG GCGTAACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	F10	hybrid_dual_ fw_070	AATGATAC GGCGACCA CCGAGATC



	A	B	C
			TACACAAG TGAGCACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	F11	hybrid_dual_ fw_071	AATGATAC GGCGACCA CCGAGATC TACACCCG GAATAACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	F12	hybrid_dual_ fw_072	AATGATAC GGCGACCA CCGAGATC TACACAATA GGCAACAC TCTTTCCC TACACGAC TCGCCAAG CTGA
	G01	hybrid_dual_ fw_073	AATGATAC GGCGACCA CCGAGATC TACACTTC AGGTCACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	G02	hybrid_dual_ fw_074	AATGATAC GGCGACCA CCGAGATC TACACTAAT CGCTACAC TCTTTCCC TACACGAC TCGCCAAG CTGA
	G03	hybrid_dual_ fw_075	AATGATAC GGCGACCA CCGAGATC TACACCTT AGGTGACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	G04	hybrid_dual_ fw_076	AATGATAC GGCGACCA CCGAGATC



	A	B	C
			TACACTTG ATCTCACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	G05	hybrid_dual_ fw_077	AATGATAC GGCGACCA CCGAGATC TACACTCT TACTGACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	G06	hybrid_dual_ fw_078	AATGATAC GGCGACCA CCGAGATC TACACGAG TACACACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	G07	hybrid_dual_ fw_079	AATGATAC GGCGACCA CCGAGATC TACACACG TTACAACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	G08	hybrid_dual_ fw_080	AATGATAC GGCGACCA CCGAGATC TACACTGG TAAGCACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	G09	hybrid_dual_ fw_081	AATGATAC GGCGACCA CCGAGATC TACACGATA GCACACAC TCTTTCCC TACACGAC TCGCCAAG CTGA
	G10	hybrid_dual_ fw_082	AATGATAC GGCGACCA CCGAGATC



	A	B	C
			TACACGAC TTAGGACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	G11	hybrid_dual_ fw_083	AATGATAC GGCGACCA CCGAGATC TACACAGT CGTTAACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	G12	hybrid_dual_ fw_084	AATGATAC GGCGACCA CCGAGATC TACACCGT GATTAACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	H01	hybrid_dual_ fw_085	AATGATAC GGCGACCA CCGAGATC TACACTCA GTAGAACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	H02	hybrid_dual_ fw_086	AATGATAC GGCGACCA CCGAGATC TACACATGT GTCTACAC TCTTTCCC TACACGAC TCGCCAAG CTGA
	H03	hybrid_dual_ fw_087	AATGATAC GGCGACCA CCGAGATC TACACCAG AATATACAC TCTTTCCC TACACGAC TCGCCAAG CTGA
	H04	hybrid_dual_ fw_088	AATGATAC GGCGACCA CCGAGATC



	A	B	C
			TACACTAA GTCTGACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	H05	hybrid_dual_ fw_089	AATGATAC GGCGACCA CCGAGATC TACACATC GACATACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	H06	hybrid_dual_ fw_090	AATGATAC GGCGACCA CCGAGATC TACACACA GCTGTACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	H07	hybrid_dual_ fw_091	AATGATAC GGCGACCA CCGAGATC TACACATCA CAAGACAC TCTTTCCC TACACGAC TCGCCAAG CTGA
	H08	hybrid_dual_ fw_092	AATGATAC GGCGACCA CCGAGATC TACACGTC GAATTACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	H09	hybrid_dual_ fw_093	AATGATAC GGCGACCA CCGAGATC TACACTGC CGATTACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	H10	hybrid_dual_ fw_094	AATGATAC GGCGACCA CCGAGATC

	A	B	C
			TACACACT AATGCACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	H11	hybrid_dual_ fw_095	AATGATAC GGCGACCA CCGAGATC TACACCTA GTAGGACA CTCTTTCC CTACACGA CTCGCCAA GCTGA
	H12	hybrid_dual_ fw_096	AATGATAC GGCGACCA CCGAGATC TACACTTAA GCCAACAC TCTTTCCC TACACGAC TCGCCAAG CTGA

PCR2 Reverse Primers

	Well Position	Name	Sequence
	A01	hybrid_dual_rev_001	CAAGCAGAAGACGGCATAACGAGATGTAAAGACGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	A02	hybrid_dual_rev_002	CAAGCAGAAGACGGCATAACGAGATTCTAAGTTGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	A03	hybrid_dual_rev_003	CAAGCAGAAGACGGCATAACGAGATTGTACATTGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	A04	hybrid_dual_rev_004	CAAGCAGAAGACGGCATAACGAGATAATACCGCGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	A05	hybrid_dual_rev_005	CAAGCAGAAGACGGCATAACGAGATTGTCTATGGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	A06	hybrid_dual_rev_006	CAAGCAGAAGACGGCATAACGAGATTCACTGTTGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	A07	hybrid_dual_rev_007	CAAGCAGAAGACGGCATAACGAGATTTGTATTTCGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT



	A08	hybrid_dual_rev_008	CAAGCAGAAGACGGCATAACGAGATAGACTAACGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	A09	hybrid_dual_rev_009	CAAGCAGAAGACGGCATAACGAGATATTGTCAAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	A10	hybrid_dual_rev_010	CAAGCAGAAGACGGCATAACGAGATGTTCTACAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	A11	hybrid_dual_rev_011	CAAGCAGAAGACGGCATAACGAGATTTTCATTCCGGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	A12	hybrid_dual_rev_012	CAAGCAGAAGACGGCATAACGAGATAAGCGTCAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	B01	hybrid_dual_rev_013	CAAGCAGAAGACGGCATAACGAGATCGTAGGATGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	B02	hybrid_dual_rev_014	CAAGCAGAAGACGGCATAACGAGATAGTTAGCGGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	B03	hybrid_dual_rev_015	CAAGCAGAAGACGGCATAACGAGATGAACGTCTGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	B04	hybrid_dual_rev_016	CAAGCAGAAGACGGCATAACGAGATAGAACACTGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	B05	hybrid_dual_rev_017	CAAGCAGAAGACGGCATAACGAGATTGTTGGCAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	B06	hybrid_dual_rev_018	CAAGCAGAAGACGGCATAACGAGATGGCTAATAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	B07	hybrid_dual_rev_019	CAAGCAGAAGACGGCATAACGAGATGCTAGAATGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	B08	hybrid_dual_rev_020	CAAGCAGAAGACGGCATAACGAGATTGTACTAAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	B09	hybrid_dual_rev_021	CAAGCAGAAGACGGCATAACGAGATAGTTACACGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	B10	hybrid_dual_rev_022	CAAGCAGAAGACGGCATAACGAGATGTACACGAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	B11	hybrid_dual_rev_023	CAAGCAGAAGACGGCATAACGAGATACTGTAAGGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	B12	hybrid_dual_rev_024	CAAGCAGAAGACGGCATAACGAGATGATCATGAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	C01	hybrid_dual_rev_025	CAAGCAGAAGACGGCATAACGAGATTGCTAACGGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	C02	hybrid_dual_rev_026	CAAGCAGAAGACGGCATAACGAGATTAATGCTCGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT



	C03	hybrid_dual_rev_027	CAAGCAGAAGACGGCATAACGAGATCTTGTTGAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	C04	hybrid_dual_rev_028	CAAGCAGAAGACGGCATAACGAGATGATACAAGGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	C05	hybrid_dual_rev_029	CAAGCAGAAGACGGCATAACGAGATATACAGCGGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	C06	hybrid_dual_rev_030	CAAGCAGAAGACGGCATAACGAGATCTGGATAGGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	C07	hybrid_dual_rev_031	CAAGCAGAAGACGGCATAACGAGATTCATCAGCGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	C08	hybrid_dual_rev_032	CAAGCAGAAGACGGCATAACGAGATTTATACCGGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	C09	hybrid_dual_rev_033	CAAGCAGAAGACGGCATAACGAGATGCTGATTCGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	C10	hybrid_dual_rev_034	CAAGCAGAAGACGGCATAACGAGATATAAGGCCGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	C11	hybrid_dual_rev_035	CAAGCAGAAGACGGCATAACGAGATAAGACTAGGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	C12	hybrid_dual_rev_036	CAAGCAGAAGACGGCATAACGAGATTAGCAACAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	D01	hybrid_dual_rev_037	CAAGCAGAAGACGGCATAACGAGATCAAGTAATGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	D02	hybrid_dual_rev_038	CAAGCAGAAGACGGCATAACGAGATTCATGATGGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	D03	hybrid_dual_rev_039	CAAGCAGAAGACGGCATAACGAGATTAGCTCTTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	D04	hybrid_dual_rev_040	CAAGCAGAAGACGGCATAACGAGATACTACATGGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	D05	hybrid_dual_rev_041	CAAGCAGAAGACGGCATAACGAGATTGTGCCATGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	D06	hybrid_dual_rev_042	CAAGCAGAAGACGGCATAACGAGATGACTTCCAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	D07	hybrid_dual_rev_043	CAAGCAGAAGACGGCATAACGAGATTATTCTCGGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	D08	hybrid_dual_rev_044	CAAGCAGAAGACGGCATAACGAGATCACTAGATGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	D09	hybrid_dual_rev_045	CAAGCAGAAGACGGCATAACGAGATCTTCGCATGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT



	D10	hybrid_dual_rev_046	CAAGCAGAAGACGGCATAACGAGATCTAGCGATGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	D11	hybrid_dual_rev_047	CAAGCAGAAGACGGCATAACGAGATTACGTTAAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	D12	hybrid_dual_rev_048	CAAGCAGAAGACGGCATAACGAGATAAGTAGTCGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	E01	hybrid_dual_rev_049	CAAGCAGAAGACGGCATAACGAGATGACGTATCGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	E02	hybrid_dual_rev_050	CAAGCAGAAGACGGCATAACGAGATACGTTTCAGGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	E03	hybrid_dual_rev_051	CAAGCAGAAGACGGCATAACGAGATAGCTATATGTGACTGGAGTTCAGACGTG GCTCTTCCGAT
	E04	hybrid_dual_rev_052	CAAGCAGAAGACGGCATAACGAGATAGTGCGAAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	E05	hybrid_dual_rev_053	CAAGCAGAAGACGGCATAACGAGATCAATGGTAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	E06	hybrid_dual_rev_054	CAAGCAGAAGACGGCATAACGAGATATTCCGATGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	E07	hybrid_dual_rev_055	CAAGCAGAAGACGGCATAACGAGATAGCTGTTGGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	E08	hybrid_dual_rev_056	CAAGCAGAAGACGGCATAACGAGATTTCGACTCGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	E09	hybrid_dual_rev_057	CAAGCAGAAGACGGCATAACGAGATCGAATGGTGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	E10	hybrid_dual_rev_058	CAAGCAGAAGACGGCATAACGAGATTAGTGTACGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	E11	hybrid_dual_rev_059	CAAGCAGAAGACGGCATAACGAGATACGCATGAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	E12	hybrid_dual_rev_060	CAAGCAGAAGACGGCATAACGAGATTTCTGACAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	F01	hybrid_dual_rev_061	CAAGCAGAAGACGGCATAACGAGATTTCTGCTAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	F02	hybrid_dual_rev_062	CAAGCAGAAGACGGCATAACGAGATTGGAAGCAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	F03	hybrid_dual_rev_063	CAAGCAGAAGACGGCATAACGAGATGATGCTCAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	F04	hybrid_dual_rev_064	CAAGCAGAAGACGGCATAACGAGATAGCCTGAAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT




	F05	hybrid_dual_rev_065	CAAGCAGAAGACGGCATAACGAGATCACCAGTAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	F06	hybrid_dual_rev_066	CAAGCAGAAGACGGCATAACGAGATGTTACACCGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	F07	hybrid_dual_rev_067	CAAGCAGAAGACGGCATAACGAGATGCGTTAGTGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	F08	hybrid_dual_rev_068	CAAGCAGAAGACGGCATAACGAGATTGACGTTTCGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	F09	hybrid_dual_rev_069	CAAGCAGAAGACGGCATAACGAGATTAGGCGTAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	F10	hybrid_dual_rev_070	CAAGCAGAAGACGGCATAACGAGATAAGTGAGCGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	F11	hybrid_dual_rev_071	CAAGCAGAAGACGGCATAACGAGATCCGGAATAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	F12	hybrid_dual_rev_072	CAAGCAGAAGACGGCATAACGAGATAATAGGCAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	G01	hybrid_dual_rev_073	CAAGCAGAAGACGGCATAACGAGATTTAGGTCGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	G02	hybrid_dual_rev_074	CAAGCAGAAGACGGCATAACGAGATTAATCGCTGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	G03	hybrid_dual_rev_075	CAAGCAGAAGACGGCATAACGAGATCTTAGGTGGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	G04	hybrid_dual_rev_076	CAAGCAGAAGACGGCATAACGAGATTTGATCTCGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	G05	hybrid_dual_rev_077	CAAGCAGAAGACGGCATAACGAGATTCTTACTGGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	G06	hybrid_dual_rev_078	CAAGCAGAAGACGGCATAACGAGATGAGTACACGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	G07	hybrid_dual_rev_079	CAAGCAGAAGACGGCATAACGAGATACGTTACAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	G08	hybrid_dual_rev_080	CAAGCAGAAGACGGCATAACGAGATTGGTAAGCGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	G09	hybrid_dual_rev_081	CAAGCAGAAGACGGCATAACGAGATGATAGCACGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	G10	hybrid_dual_rev_082	CAAGCAGAAGACGGCATAACGAGATGACTTAGGGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	G11	hybrid_dual_rev_083	CAAGCAGAAGACGGCATAACGAGATAGTCGTTAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT



	G12	hybrid_dual_rev_084	CAAGCAGAAGACGGCATAACGAGATCGTGATTAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	H01	hybrid_dual_rev_085	CAAGCAGAAGACGGCATAACGAGATTCAGTAGAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	H02	hybrid_dual_rev_086	CAAGCAGAAGACGGCATAACGAGATATGTGTCTGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	H03	hybrid_dual_rev_087	CAAGCAGAAGACGGCATAACGAGATCAGAATATGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	H04	hybrid_dual_rev_088	CAAGCAGAAGACGGCATAACGAGATTAAGTCTGGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	H05	hybrid_dual_rev_089	CAAGCAGAAGACGGCATAACGAGATATCGACATGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	H06	hybrid_dual_rev_090	CAAGCAGAAGACGGCATAACGAGATACAGCTGTGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	H07	hybrid_dual_rev_091	CAAGCAGAAGACGGCATAACGAGATATCACAAGGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	H08	hybrid_dual_rev_092	CAAGCAGAAGACGGCATAACGAGATGTGCAATTGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	H09	hybrid_dual_rev_093	CAAGCAGAAGACGGCATAACGAGATTGCCGATTGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	H10	hybrid_dual_rev_094	CAAGCAGAAGACGGCATAACGAGATACTAATGCGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	H11	hybrid_dual_rev_095	CAAGCAGAAGACGGCATAACGAGATCTAGTAGGGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	H12	hybrid_dual_rev_096	CAAGCAGAAGACGGCATAACGAGATTTAAGCCAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT

Preparation of PCR2 Primer plates

- 1 To allow for high levels of combinations of indexing primers, we generate multiple forward primer plates (12) from a single 96-well plate.
- 2 By using these 12 plates in combination with a single 96-well plate of reverse primers, we can index 1,152 samples.
- 3 Additional primers can be designed to increase this multiplexing. Here we present one plate of forward and one plate of reverse primers to allow for this level of combination.
- 4 To generate forward primer plates, we shift the column to the right one column when replating to a new plate (e.g. column 1 becomes column 2 and column 12 becomes column 1). See the figure below. Be sure to uniquely label the plate so you know which primers are used in library preparation.



Note

We do NOT do this for reverse primers.

Generate Primers mix

- 5 Spin lyophilized Primer.
- 6 Resuspend in MGW using the provided volume for a 100 micromolar (μM) stock.
- 7 Vortex and spin down briefly.
- 8 Let the mixture sit for about 00:10:00 for complete dissolution.
- 9 Prepare the stock primer pool (100 micromolar (μM)) by mixing 1 volume of sp-rtss-repeat_v4_F_MIP with 1 volume of csp-rtss-repeat_v4_R_MIP.
- 10 Dilute primer pool with MGW to obtain 20 micromolar (μM) working solution: 1 volume of stock + 4 volumes of MGW.

PCR1

2m

- 11 Pre-set the thermocycler before starting master mix preparation using the programme below.

A	B	C	D
Step	Temperature (°C)	Time (minute)	Number of cycles
Step1 (preheat)	95	2	1
Step2 (Activation)	95	0.5	35
Step3 (Denaturation)	94	0.5	
Step4 (Annealing)	59	1	

	A	B	C	D
	Step5 (final Extension)	72	10	1
	Step6	4	hold	NA

12 All sample plates and controls are to be kept on cold blocks - switch out cold blocks if they are getting warm.

13 All reagents are to be kept  On ice and the tube used to prepare the master mix should also remain  On ice .



14 Centrifuge sample plates at  4000 rpm, 00:01:00 .

1m



15 Prepare the master mix using the template (here 96 reactions) below and adjust it to the number of reactions.

These volumes could be halved but not tested at that volume.

	A	B	C
	Master mix	x1	x96* 1.1
	Gotaq	25	2,640
	Molecular grade water	19	2,6006.4
	Primer pool (P2_Px1Block3_v2_F and P2_Px1Frag3_v8_R)	1	105.6
	Template DNA	5	

16 Divide the total master mix volume by 12 and aliquot the designated volume of into 12 tube strip on cold block to allow for multichannel pipetting.

17 Using a multichannel, pipette  16 µL of master mix into each well of 96 well plate..



- 18 Using a multichannel, plate  4 µL of sample into each well on the PCR1 plate.
- Add at least one positive control (lab strain,  4 µL) and one PCR1 negative control (MGW,  4 µL).



- 19 Centrifuge the plate at  4000 rpm, 00:01:00 .

1m



- 20 Place the plate in the thermocycler and start running the pre-set PCR program.



PCR2 or barcoding PCR

3m

- 21 Once the first PCR is complete, pre-set the thermocycler before starting the master mix preparation for PCR2 using the program below.



	A	B	C	D
	Step	Temperature (°C)	Time (minute)	Number of cycles
	Step1	98	hold	NA
	Step2	98	0:30	NA
	Step3	95	0:10	x10
	Step4	62	8:0	x10
	Step5	72	0:20	x10
	Step6	72	1:0	NA
	Step7	4	hold	NA

- 22 Prepare the master mix for PCR2 using the template (here 96 reactions) below and adjust it to the number of reactions.



	A	B	C
	Master mix	x1	x96* 1.1
	Gotaq	10	1056
	Molecular grade water	6	633.6

23 Divide the total master mix volume by 12 and aliquot the designated volume ( 140.8 μL for 96 sample run) into 12 tube strip on cold block to allow for multichannel pipetting. 

24 Using a multichannel, pipette  16 μL of master mix into each well. 

25 Spin down the forward and reverse PCR2 primer plates needed for the experiment at  4000 rpm, 00:01:00 and add  1 μL of forward primer and  1 μL of reverse primer into each well using the multichannel.  

- On your sample sheet, please denote the name of the forward and reverse primer plate used and the well of the sample to be sure you know the correct index for the sample you are genotyping.

26 In the post-PCR lab, add  2 μL of the PCR1 product to the master mix using a multichannel. 

27 **Optional:** shake the plate gently on a shaker for  00:01:00 to mix the master mix added to the samples.  

28 Centrifuge the plate at  4000 rpm, 00:01:00 .  

29 Place the plate in the thermocycler and start running the pre-set PCR program. 

Check a random set of 10 wells per plate for amplification by gel electrophoresis based on your lab protocol (including the NTC and positive control).



PCR Product Pooling and Cleanup

10m

- 30 Pool  5 μL of PCR2 product into a single Eppendorf tube. If the pool is very large,  3 μL is sufficient. 
- 31 Bring Ampure Beads to  Room temperature . This can also be done with Spherotech beads (DIY). 
- 32 Add an equal volume of beads to the pool (1X bead clean). Let it incubate for at least  00:10:00 at  Room temperature . For a higher yield, place the tube on a Hula Mixer or invert it in your hand.  
- 33 Place the pool on a magnet and wait until the beads have pelleted on the magnet and the solution is clear. Prepare enough fresh 80% ethanol to cover the beads for 2 washes.
- 34 Remove the supernatant and reserve it in another tube.
- 35 Wash with 80% ethanol for at least  00:00:30 . Discard the ethanol. 
- 36 Repeat the previous step for a second wash. 
- 37 Remove the ethanol and let the beads dry for around  00:01:00 (do not let the beads dry until cracking).
- 38 Elute in  51 μL of 1X low EDTA TE Buffer. If the pool is very large, use your discretion and elute in a higher volume (we eluted in  150 μL with 6 96-well plates of sample). Incubate at  Room temperature for  00:10:00 .    10m
- 39 Pellet the beads on a magnet and remove the elution into another tube.
- 40 Quantify the elution with a Qubit fluorometer.



Nanopore Library preparation

- 41 We use the SQK-LSK114 Kit (Ligation Sequencing) and have adapted this protocol from the Ligation Sequencing Amplicons V14 protocol from Nanopore. We use the P2 sequencing machine, but this protocol could be adapted for Minion.

Nanopore Library preparation - End Prep

26m 30s

- 42 Transfer 200 fmol of cleaned DNA into a PCR tube or plate. Calculation for fmol can be done on this website: <https://www.promega.com/resources/tools/biomath/> (1 pmol = 1,000 fmol). Adjust the volume to  50 μL with nuclease-free water. 

- 43 Add the following to a PCR tube: 

A	B
Reagent	Volume
DNA Template	50 μL
Ultra II End Prep Reaction Buffer	7 μL
Ultra II End Prep Enzyme Mix	3 μL
Total	60μL

- 44 Mix thoroughly by pipetting.  

- 45 Incubate at  20 $^{\circ}\text{C}$ for  00:05:00 and  65 $^{\circ}\text{C}$ for  00:05:00 .  

10m

- 46 Spin down and transfer the sample to a clean 1.5ml Eppendorf tube. 

- 47 Bring beads to  Room temperature and add  60 μL to the DNA sample.  



48 Incubate for 00:10:00 at Room temperature . For increased yield, place the tube on a Hula Mixer or invert in your hand. Meanwhile, prepare 1 mL of 80% ethanol.

10m



49 Pellet on a magnet and remove the supernatant - keep the supernatant for extra security.

50 Wash with 500 μ L of 80% ethanol. Remove after 00:00:30 .

30s



51 Repeat the previous step for a total of 2 washes.

51.1 Wash with 500 μ L of 80% ethanol. Remove after 00:00:30 (1/2).

30s



51.2 Wash with 500 μ L of 80% ethanol. Remove after 00:00:30 (2/2).

30s



52 Let the beads dry for around 00:01:00 . Do not let the beads crack on the side of the tube.

53 Elute in 61 μ L of nuclease-free water. Let it incubate at Room temperature for 00:05:00 .

5m



54 Pellet on the magnet and remove the elution.

55 Quantify the End-prepped DNA with a qubit fluorometer. This library can be stored at 4 °C for short-term storage, or at -20 °C or -80 °C for long-term storage.



Nanopore Library preparation - Adapter Ligation

21m 30s

56 Thaw the Short Fragment Buffer (SFB) and the Elution Buffer (EB) at  Room temperature . Place the Ligation Adapter (LA), Ligation Buffer (LNB), and Quick T4 Ligase  On ice .



57 Mix the following in a 1.5 ml Eppendorf tube, mixing 10-20 times after each addition:



A	B
Reagent	Volume
End-prepped DNA	60 μ L
Ligation Buffer (LNB)	25 μ L
NEBNext Quick T4 Ligase	10 μ L
Ligation Adapter (LA)	5 μ L
Total	100 μL

58 Mix well by pipetting. Incubate at  Room temperature for at least  00:10:00 . A longer incubation of  00:20:00 will increase the amount of ligated adapters.



59 Bring the beads to  Room temperature and resuspend.



60 Add  60 μ L of beads and mix by inverting. Incubate for  00:10:00 at  Room temperature . For higher yield, invert in your hand or place on a Hula mixer

10m



61 Pellet on a magnet and remove the supernatant.

62 Wash with  250 μ L of the Short Fragment buffer. Remove after around  00:00:30 .

30s





63 Repeat the previous step for a total of 2 washes.



63.1 Wash with  250 μL of the Short Fragment buffer. Remove after around  00:00:30 (1/2).

30s

63.2 Wash with  250 μL of the Short Fragment buffer. Remove after around  00:00:30 (2/2).

30s

64 Allow beads to dry on the magnet, and pipette out any residual buffer. Do not let beads dry to the point of cracking.

65 Elute in  25 μL of the Elution Buffer (EB). Incubate for  00:10:00 at  Room temperature .

10m



66 Pellet on a magnet and pipette the elution into a tube (ideally a screw cap tube for best storage).

67 Quantify the library with a qubit fluorometer. Be sure that you have at least 120 fmol of DNA, which is around  25 ng of DNA for a 4CAST library. This library can be stored at  4 $^{\circ}\text{C}$ for short-term storage, or at  -20 $^{\circ}\text{C}$ or  -80 $^{\circ}\text{C}$ for long-term storage.



Nanopore sequencing - Loading the Flow Cell

20m

68 Take the Promethion Flow Cell you plan to use out of the  4 $^{\circ}\text{C}$ refrigerator at least  00:20:00 before use. It is highly recommended to do a flow cell check before any sequencing.



69 Determine 100 - 120 fmols of library and add to an Eppendorf tube (Nanopore recommends 20 fmol for sequencing; however, this leads to underloading of the flow cell in most cases. If you have less than 100 fmols, add the entire library).  25 ng of 4CAST Library is around 120 fmols.





70 Bring up the 100 - 120 fmols of library to  32 μL with the elution buffer (EB). 

71 After  00:20:00 at  Room temperature , gently insert the flow cell into one of the ports of the P2. The lights will turn blue if inserted correctly, and red if inserted incorrectly. 

72 In a 1.5 Eppendorf tube, mix the following to create the flow cell priming mix: 

A	B
Reagent	Volume
Flow Cell Flush (FCF)	1170 μL
Flow Cell Tether (FCT)	30 μL
Total	1200 μL

73 Set a P1000 pipette to 200. Open the inlet port and place the tip inside the port. Turn the dial until it shows 220-230. A small amount of buffer should enter the tip. Discard the tip and buffer. 

74 Load  500 μL of the flow cell priming mix into the inlet port. Load by slowly twisting the dial in the counterclockwise direction. Wait  00:05:00 .  

75 While the priming mix is incubating, prepare the library: 

A	B
Reagent	Volume
Sequencing Buffer (SB)	100 μL
Library Beads (LIB)	68 μL



A	B
DNA Library	32 µL
Total	200 µL

76 After the  00:05:00 incubation, add another  500 µL of flow cell priming mix to the inlet port.

5m



77 Immediately add the  200 µL of the library slowly into the inlet port.



78 Close the inlet port and wait for the library to incubate for  00:10:00 before starting any sequencing.

10m



Nanopore sequencing - Flow cell wash

1h

79 A flow cell wash should be done as soon as possible after the sequencing run is complete.

80 Place the Wash Mix (WMX)  On ice , and thaw the Wash Diluent (DIL) at  Room temperature .



81 Vortex the Wash Diluent thoroughly.



82 Make the Flow Cell Wash Mix by adding the following to a 1.5ml Eppendorf tube:



A	B
Reagent	Volume
Wash Diluent (DIL)	398 µL
Wash Mix (WMX)	2 µL



83 Pipette to mix and place  On ice .



84 Stop the sequencing run on Minknow, and leave the flow cell in the port.

85 Make sure the inlet port is closed, and remove waste from port 2 or port 3.

86 Open the inlet port and set the P1000 to 200. Remove the air bubble by turning the wheel until it reaches 220-230, and a small volume of buffer enters the pipette tip.

87 Load  400 μL of the Flow Cell Wash Mix into the inlet port. Incubate for  01:00:00 .

1h



88 After  01:00:00 , remove the waste again from ports 2 or 3. Make sure the inlet port is closed.

89 Remove the air bubble if it is present after removing the waste. Then add  500 μL of Storage Buffer (S). The flow cell can now be placed back into a  4 $^{\circ}\text{C}$ refrigerator for long-term storage. Incubate the flow cell at least  Overnight before doing another flow cell check.

