

Oct 24, 2025

Version 1

IDEEL: UNC Implementation of Targeted PX1 Sequencing-MIP to Nanopore V.1



DOI

<https://dx.doi.org/10.17504/protocols.io.yxmvmmbbcbg3p/v1>

Infectious Disease Epidemiology and Ecology Lab¹

¹Infectious Disease Epidemiology and Ecology Lab, University of North Carolina at Chapel Hill, Chapel Hill, NC, USA

IDEEL



Infectious Disease Epidemiology and Ecology Lab

University of North Carolina

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.yxmvmmbbcbg3p/v1>

Protocol Citation: Infectious Disease Epidemiology and Ecology Lab 2025. IDEEL: UNC Implementation of Targeted PX1 Sequencing- MIP to Nanopore. protocols.io <https://dx.doi.org/10.17504/protocols.io.yxmvmmbbcbg3p/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 16, 2025

Last Modified: October 24, 2025

Protocol Integer ID: 230584

Keywords: PCR, Barcoding, Nanopore Library preparation, Amplification, Adapter Ligation, targeted sequencing of the px1 locus, targeted sequencing, integrating molecular inversion probe, molecular inversion probe, nanopore the protocol, px1 locus, oxford nanopore, read sequencing, nanopore

Abstract

The protocol describes a novel workflow integrating Molecular Inversion Probes (MIPs) with Oxford Nanopore long-read sequencing for targeted sequencing of the PX1 locus.

Materials

Required Materials:

- Molecular Grade Water (MGW)
- Primers

P2_Px1Block3_v2_F: GACTCGCCAAGCTGAAGNNNNAGAGAACAGGAAAAAGTTTGTACGT
P2_Px1Frag3_v8_R: ACGTGTGCTCTTCCGATCTNNNNTCATCACTGTCGCTACTTTTGTGTC

-  GoTaq® Long PCR Master Mix, 2X **Promega Catalog #M402A**
- Illumina FWD and REV barcode primers
- Gel electrophoresis reagents and equipment.
- 1X TE Low EDTA Buffer
- PromethION 2 solo and device specific flow cells
- ONT Ligation sequencing kit (SQK-LSK114 Kit)

PCR Round 2 Primers

Forward Primers:

	Well Position	Name	Sequence
	A01	hybrid_dual_fw_001	AATGATACGGCGACCACCGAGATCTACACGTTAAGACACACTCTTTCCCTACAC GACTCGCCAAGCTGA
	A02	hybrid_dual_fw_002	AATGATACGGCGACCACCGAGATCTACACTCTAAGTTACACTCTTTCCCTACAC GACTCGCCAAGCTGA
	A03	hybrid_dual_fw_003	AATGATACGGCGACCACCGAGATCTACACTGTACATTACACTCTTTCCCTACAC GACTCGCCAAGCTGA
	A04	hybrid_dual_fw_004	AATGATACGGCGACCACCGAGATCTACACAATACCGCACACTCTTTCCCTACAC GACTCGCCAAGCTGA
	A05	hybrid_dual_fw_005	AATGATACGGCGACCACCGAGATCTACACTGTCTATGACACTCTTTCCCTACAC GACTCGCCAAGCTGA
	A06	hybrid_dual_fw_006	AATGATACGGCGACCACCGAGATCTACACTCACTGTTACACTCTTTCCCTACAC GACTCGCCAAGCTGA
	A07	hybrid_dual_fw_007	AATGATACGGCGACCACCGAGATCTACACTTGTATTACACTCTTTCCCTACAC GACTCGCCAAGCTGA
	A08	hybrid_dual_fw_008	AATGATACGGCGACCACCGAGATCTACACAGACTAACACACTCTTTCCCTACAC GACTCGCCAAGCTGA



	A09	hybrid_dual_fw_009	AATGATACGGCGACCACCGAGATCTACACATTGTCAAACACTCTTTCCCTACACGACTCGCCAAGCTGA
	A10	hybrid_dual_fw_010	AATGATACGGCGACCACCGAGATCTACACGTTCTACAACACTCTTTCCCTACACGACTCGCCAAGCTGA
	A11	hybrid_dual_fw_011	AATGATACGGCGACCACCGAGATCTACACTTCATTGACACTCTTTCCCTACACGACTCGCCAAGCTGA
	A12	hybrid_dual_fw_012	AATGATACGGCGACCACCGAGATCTACACAAGCGTCAACACTCTTTCCCTACACGACTCGCCAAGCTGA
	B01	hybrid_dual_fw_013	AATGATACGGCGACCACCGAGATCTACACCGTAGGATACACTCTTTCCCTACACGACTCGCCAAGCTGA
	B02	hybrid_dual_fw_014	AATGATACGGCGACCACCGAGATCTACACAGTTAGCGACACTCTTTCCCTACACGACTCGCCAAGCTGA
	B03	hybrid_dual_fw_015	AATGATACGGCGACCACCGAGATCTACACGAACGTCTACACTCTTTCCCTACACGACTCGCCAAGCTGA
	B04	hybrid_dual_fw_016	AATGATACGGCGACCACCGAGATCTACACAGAACACTACACTCTTTCCCTACACGACTCGCCAAGCTGA
	B05	hybrid_dual_fw_017	AATGATACGGCGACCACCGAGATCTACACTGTTGGCAACACTCTTTCCCTACACGACTCGCCAAGCTGA
	B06	hybrid_dual_fw_018	AATGATACGGCGACCACCGAGATCTACACGGCTAATAACACTCTTTCCCTACACGACTCGCCAAGCTGA
	B07	hybrid_dual_fw_019	AATGATACGGCGACCACCGAGATCTACACGCTAGAATACACTCTTTCCCTACACGACTCGCCAAGCTGA
	B08	hybrid_dual_fw_020	AATGATACGGCGACCACCGAGATCTACACTGTACTAAACACTCTTTCCCTACACGACTCGCCAAGCTGA
	B09	hybrid_dual_fw_021	AATGATACGGCGACCACCGAGATCTACACAGTTACACACACTCTTTCCCTACACGACTCGCCAAGCTGA
	B10	hybrid_dual_fw_022	AATGATACGGCGACCACCGAGATCTACACGTACACGAACACTCTTTCCCTACACGACTCGCCAAGCTGA
	B11	hybrid_dual_fw_023	AATGATACGGCGACCACCGAGATCTACACACTGTAAGACACTCTTTCCCTACACGACTCGCCAAGCTGA
	B12	hybrid_dual_fw_024	AATGATACGGCGACCACCGAGATCTACACGATCATGAACACTCTTTCCCTACACGACTCGCCAAGCTGA
	C01	hybrid_dual_fw_025	AATGATACGGCGACCACCGAGATCTACACTGCTAACGACACTCTTTCCCTACACGACTCGCCAAGCTGA
	C02	hybrid_dual_fw_026	AATGATACGGCGACCACCGAGATCTACACTAATGCTCACACTCTTTCCCTACACGACTCGCCAAGCTGA
	C03	hybrid_dual_fw_027	AATGATACGGCGACCACCGAGATCTACACCTTGTTGAACACTCTTTCCCTACACGACTCGCCAAGCTGA



	C04	hybrid_dual_fw_028	AATGATACGGCGACCACCGAGATCTACACGATACAAGACACTCTTTCCCTACACGACTCGCCAAGCTGA
	C05	hybrid_dual_fw_029	AATGATACGGCGACCACCGAGATCTACACATACAGCGACACTCTTTCCCTACACGACTCGCCAAGCTGA
	C06	hybrid_dual_fw_030	AATGATACGGCGACCACCGAGATCTACACCTGGATAGACACTCTTTCCCTACACGACTCGCCAAGCTGA
	C07	hybrid_dual_fw_031	AATGATACGGCGACCACCGAGATCTACACTCATCAGCACACTCTTTCCCTACACGACTCGCCAAGCTGA
	C08	hybrid_dual_fw_032	AATGATACGGCGACCACCGAGATCTACACTTATACCGACACTCTTTCCCTACACGACTCGCCAAGCTGA
	C09	hybrid_dual_fw_033	AATGATACGGCGACCACCGAGATCTACACGCTGATTCACACTCTTTCCCTACACGACTCGCCAAGCTGA
	C10	hybrid_dual_fw_034	AATGATACGGCGACCACCGAGATCTACACATAAGGCCACACTCTTTCCCTACACGACTCGCCAAGCTGA
	C11	hybrid_dual_fw_035	AATGATACGGCGACCACCGAGATCTACACAAGACTAGACACTCTTTCCCTACACGACTCGCCAAGCTGA
	C12	hybrid_dual_fw_036	AATGATACGGCGACCACCGAGATCTACACTAGCAACAACACTCTTTCCCTACACGACTCGCCAAGCTGA
	D01	hybrid_dual_fw_037	AATGATACGGCGACCACCGAGATCTACACCAAGTAATACACTCTTTCCCTACACGACTCGCCAAGCTGA
	D02	hybrid_dual_fw_038	AATGATACGGCGACCACCGAGATCTACACTCATGATGACACTCTTTCCCTACACGACTCGCCAAGCTGA
	D03	hybrid_dual_fw_039	AATGATACGGCGACCACCGAGATCTACACTAGCTCTTACACTCTTTCCCTACACGACTCGCCAAGCTGA
	D04	hybrid_dual_fw_040	AATGATACGGCGACCACCGAGATCTACACACTACATGACACTCTTTCCCTACACGACTCGCCAAGCTGA
	D05	hybrid_dual_fw_041	AATGATACGGCGACCACCGAGATCTACACTGTGCCATACACTCTTTCCCTACACGACTCGCCAAGCTGA
	D06	hybrid_dual_fw_042	AATGATACGGCGACCACCGAGATCTACACGACTTCCAACACTCTTTCCCTACACGACTCGCCAAGCTGA
	D07	hybrid_dual_fw_043	AATGATACGGCGACCACCGAGATCTACACTATTCTCGACACTCTTTCCCTACACGACTCGCCAAGCTGA
	D08	hybrid_dual_fw_044	AATGATACGGCGACCACCGAGATCTACACCACTAGATACACTCTTTCCCTACACGACTCGCCAAGCTGA
	D09	hybrid_dual_fw_045	AATGATACGGCGACCACCGAGATCTACACCTTCGCATACACTCTTTCCCTACACGACTCGCCAAGCTGA
	D10	hybrid_dual_fw_046	AATGATACGGCGACCACCGAGATCTACACCTAGCGATACACTCTTTCCCTACACGACTCGCCAAGCTGA



	D11	hybrid_dual_fw_047	AATGATACGGCGACCACCGAGATCTACACTACGTTAAACACTCTTTCCCTACACGACTCGCCAAGCTGA
	D12	hybrid_dual_fw_048	AATGATACGGCGACCACCGAGATCTACACAAGTAGTCACACTCTTTCCCTACACGACTCGCCAAGCTGA
	E01	hybrid_dual_fw_049	AATGATACGGCGACCACCGAGATCTACACGACGTATCACACTCTTTCCCTACACGACTCGCCAAGCTGA
	E02	hybrid_dual_fw_050	AATGATACGGCGACCACCGAGATCTACACACGTTTCAGACACTCTTTCCCTACACGACTCGCCAAGCTGA
	E03	hybrid_dual_fw_051	AATGATACGGCGACCACCGAGATCTACACAGCTATATACACTCTTTCCCTACACGACTCGCCAAGCTGA
	E04	hybrid_dual_fw_052	AATGATACGGCGACCACCGAGATCTACACAGTGCGAAACACTCTTTCCCTACACGACTCGCCAAGCTGA
	E05	hybrid_dual_fw_053	AATGATACGGCGACCACCGAGATCTACACCAATGGTAACACTCTTTCCCTACACGACTCGCCAAGCTGA
	E06	hybrid_dual_fw_054	AATGATACGGCGACCACCGAGATCTACACATTCCGATACACTCTTTCCCTACACGACTCGCCAAGCTGA
	E07	hybrid_dual_fw_055	AATGATACGGCGACCACCGAGATCTACACAGCTGTTGACACTCTTTCCCTACACGACTCGCCAAGCTGA
	E08	hybrid_dual_fw_056	AATGATACGGCGACCACCGAGATCTACACTTCGACTCACACTCTTTCCCTACACGACTCGCCAAGCTGA
	E09	hybrid_dual_fw_057	AATGATACGGCGACCACCGAGATCTACACCGAATGGTACACTCTTTCCCTACACGACTCGCCAAGCTGA
	E10	hybrid_dual_fw_058	AATGATACGGCGACCACCGAGATCTACACTAGTGTACACACTCTTTCCCTACACGACTCGCCAAGCTGA
	E11	hybrid_dual_fw_059	AATGATACGGCGACCACCGAGATCTACACACGCATGAACACTCTTTCCCTACACGACTCGCCAAGCTGA
	E12	hybrid_dual_fw_060	AATGATACGGCGACCACCGAGATCTACACTTCGTACAACACTCTTTCCCTACACGACTCGCCAAGCTGA
	F01	hybrid_dual_fw_061	AATGATACGGCGACCACCGAGATCTACACTTCTGCTAACACTCTTTCCCTACACGACTCGCCAAGCTGA
	F02	hybrid_dual_fw_062	AATGATACGGCGACCACCGAGATCTACACTGGAAGCAACACTCTTTCCCTACACGACTCGCCAAGCTGA
	F03	hybrid_dual_fw_063	AATGATACGGCGACCACCGAGATCTACACGATGCTCAACACTCTTTCCCTACACGACTCGCCAAGCTGA
	F04	hybrid_dual_fw_064	AATGATACGGCGACCACCGAGATCTACACAGCCTGAAACACTCTTTCCCTACACGACTCGCCAAGCTGA
	F05	hybrid_dual_fw_065	AATGATACGGCGACCACCGAGATCTACACCACCAGTAACACTCTTTCCCTACACGACTCGCCAAGCTGA



	F06	hybrid_dual_fw_066	AATGATACGGCGACCACCGAGATCTACACGTTACACCACACTCTTTCCCTACAC GACTCGCCAAGCTGA
	F07	hybrid_dual_fw_067	AATGATACGGCGACCACCGAGATCTACACGCGTTAGTACACTCTTTCCCTACAC GACTCGCCAAGCTGA
	F08	hybrid_dual_fw_068	AATGATACGGCGACCACCGAGATCTACACTGACGTTACACTCTTTCCCTACAC GACTCGCCAAGCTGA
	F09	hybrid_dual_fw_069	AATGATACGGCGACCACCGAGATCTACACTAGGCGTAACACTCTTTCCCTACAC GACTCGCCAAGCTGA
	F10	hybrid_dual_fw_070	AATGATACGGCGACCACCGAGATCTACACAAGTGAGCACACTCTTTCCCTACA CGACTCGCCAAGCTGA
	F11	hybrid_dual_fw_071	AATGATACGGCGACCACCGAGATCTACACCCGGAATAACACTCTTTCCCTACA CGACTCGCCAAGCTGA
	F12	hybrid_dual_fw_072	AATGATACGGCGACCACCGAGATCTACACAATAGGCAACACTCTTTCCCTACAC GACTCGCCAAGCTGA
	G01	hybrid_dual_fw_073	AATGATACGGCGACCACCGAGATCTACACTTCAGGTCACACTCTTTCCCTACAC GACTCGCCAAGCTGA
	G02	hybrid_dual_fw_074	AATGATACGGCGACCACCGAGATCTACACTAATCGCTACACTCTTTCCCTACAC GACTCGCCAAGCTGA
	G03	hybrid_dual_fw_075	AATGATACGGCGACCACCGAGATCTACACCTTAGGTGACACTCTTTCCCTACAC GACTCGCCAAGCTGA
	G04	hybrid_dual_fw_076	AATGATACGGCGACCACCGAGATCTACACTTGATCTCACACTCTTTCCCTACAC GACTCGCCAAGCTGA
	G05	hybrid_dual_fw_077	AATGATACGGCGACCACCGAGATCTACACTCTTACTGACACTCTTTCCCTACAC GACTCGCCAAGCTGA
	G06	hybrid_dual_fw_078	AATGATACGGCGACCACCGAGATCTACACGAGTACACACACTCTTTCCCTACAC GACTCGCCAAGCTGA
	G07	hybrid_dual_fw_079	AATGATACGGCGACCACCGAGATCTACACACGTTACAACACTCTTTCCCTACAC GACTCGCCAAGCTGA
	G08	hybrid_dual_fw_080	AATGATACGGCGACCACCGAGATCTACACTGGTAAGCACACTCTTTCCCTACAC GACTCGCCAAGCTGA
	G09	hybrid_dual_fw_081	AATGATACGGCGACCACCGAGATCTACACGATAGCACACACTCTTTCCCTACAC GACTCGCCAAGCTGA
	G10	hybrid_dual_fw_082	AATGATACGGCGACCACCGAGATCTACACGACTTAGGACACTCTTTCCCTACAC GACTCGCCAAGCTGA
	G11	hybrid_dual_fw_083	AATGATACGGCGACCACCGAGATCTACACAGTCGTTAACACTCTTTCCCTACAC GACTCGCCAAGCTGA
	G12	hybrid_dual_fw_084	AATGATACGGCGACCACCGAGATCTACACCGTGATTAACACTCTTTCCCTACAC GACTCGCCAAGCTGA



	H01	hybrid_dual_fw_085	AATGATACGGCGACCACCGAGATCTACACTCAGTAGAACACTCTTTCCCTACAC GACTCGCCAAGCTGA
	H02	hybrid_dual_fw_086	AATGATACGGCGACCACCGAGATCTACACATGTGTCTACACTCTTTCCCTACAC GACTCGCCAAGCTGA
	H03	hybrid_dual_fw_087	AATGATACGGCGACCACCGAGATCTACACCAGAATATACACTCTTTCCCTACAC GACTCGCCAAGCTGA
	H04	hybrid_dual_fw_088	AATGATACGGCGACCACCGAGATCTACACTAAGTCTGACACTCTTTCCCTACAC GACTCGCCAAGCTGA
	H05	hybrid_dual_fw_089	AATGATACGGCGACCACCGAGATCTACACATCGACATACACTCTTTCCCTACAC GACTCGCCAAGCTGA
	H06	hybrid_dual_fw_090	AATGATACGGCGACCACCGAGATCTACACACAGCTGTACACTCTTTCCCTACAC GACTCGCCAAGCTGA
	H07	hybrid_dual_fw_091	AATGATACGGCGACCACCGAGATCTACACATCACAAGACACTCTTTCCCTACAC GACTCGCCAAGCTGA
	H08	hybrid_dual_fw_092	AATGATACGGCGACCACCGAGATCTACACGTGGAATTACACTCTTTCCCTACAC GACTCGCCAAGCTGA
	H09	hybrid_dual_fw_093	AATGATACGGCGACCACCGAGATCTACACTGCCGATTACACTCTTTCCCTACAC GACTCGCCAAGCTGA
	H10	hybrid_dual_fw_094	AATGATACGGCGACCACCGAGATCTACACACTAATGCACACTCTTTCCCTACAC GACTCGCCAAGCTGA
	H11	hybrid_dual_fw_095	AATGATACGGCGACCACCGAGATCTACACCTAGTAGGACACTCTTTCCCTACAC GACTCGCCAAGCTGA
	H12	hybrid_dual_fw_096	AATGATACGGCGACCACCGAGATCTACACTTAAGCCAACACTCTTTCCCTACAC GACTCGCCAAGCTGA

Reverse Primers

	Well Position	Name	Sequence
	A01	hybrid_dual_rev_001	CAAGCAGAAGACGGCATAACGAGATGTTAAGACGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	A02	hybrid_dual_rev_002	CAAGCAGAAGACGGCATAACGAGATTCTAAGTTGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	A03	hybrid_dual_rev_003	CAAGCAGAAGACGGCATAACGAGATTGTACATTGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	A04	hybrid_dual_rev_004	CAAGCAGAAGACGGCATAACGAGATAATACCGCGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT



	A05	hybrid_dual_rev_005	CAAGCAGAAGACGGCATAACGAGATTGTCTATGGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	A06	hybrid_dual_rev_006	CAAGCAGAAGACGGCATAACGAGATTCACTGTTGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	A07	hybrid_dual_rev_007	CAAGCAGAAGACGGCATAACGAGATTTGTATTTCGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	A08	hybrid_dual_rev_008	CAAGCAGAAGACGGCATAACGAGATAGACTAACGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	A09	hybrid_dual_rev_009	CAAGCAGAAGACGGCATAACGAGATATTGTCAAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	A10	hybrid_dual_rev_010	CAAGCAGAAGACGGCATAACGAGATGTTCTACAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	A11	hybrid_dual_rev_011	CAAGCAGAAGACGGCATAACGAGATTTTCATTTCGGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	A12	hybrid_dual_rev_012	CAAGCAGAAGACGGCATAACGAGATAAGCGTCAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	B01	hybrid_dual_rev_013	CAAGCAGAAGACGGCATAACGAGATCGTAGGATGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	B02	hybrid_dual_rev_014	CAAGCAGAAGACGGCATAACGAGATAGTTAGCGGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	B03	hybrid_dual_rev_015	CAAGCAGAAGACGGCATAACGAGATGAACGTCTGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	B04	hybrid_dual_rev_016	CAAGCAGAAGACGGCATAACGAGATAGAACACTGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	B05	hybrid_dual_rev_017	CAAGCAGAAGACGGCATAACGAGATTGTTGGCAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	B06	hybrid_dual_rev_018	CAAGCAGAAGACGGCATAACGAGATGGCTAATAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	B07	hybrid_dual_rev_019	CAAGCAGAAGACGGCATAACGAGATGCTAGAATGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	B08	hybrid_dual_rev_020	CAAGCAGAAGACGGCATAACGAGATTGTACTAAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	B09	hybrid_dual_rev_021	CAAGCAGAAGACGGCATAACGAGATAGTTACACGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	B10	hybrid_dual_rev_022	CAAGCAGAAGACGGCATAACGAGATGTACACGAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	B11	hybrid_dual_rev_023	CAAGCAGAAGACGGCATAACGAGATACTGTAAGGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT



	B12	hybrid_dual_rev_024	CAAGCAGAAGACGGCATAACGAGATGATCATGAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	C01	hybrid_dual_rev_025	CAAGCAGAAGACGGCATAACGAGATTGCTAACGGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	C02	hybrid_dual_rev_026	CAAGCAGAAGACGGCATAACGAGATTAATGCTCGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	C03	hybrid_dual_rev_027	CAAGCAGAAGACGGCATAACGAGATCTTGTTGAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	C04	hybrid_dual_rev_028	CAAGCAGAAGACGGCATAACGAGATGATACAAGGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	C05	hybrid_dual_rev_029	CAAGCAGAAGACGGCATAACGAGATATACAGCGGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	C06	hybrid_dual_rev_030	CAAGCAGAAGACGGCATAACGAGATCTGGATAGGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	C07	hybrid_dual_rev_031	CAAGCAGAAGACGGCATAACGAGATTCATCAGCGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	C08	hybrid_dual_rev_032	CAAGCAGAAGACGGCATAACGAGATTTATACCGGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	C09	hybrid_dual_rev_033	CAAGCAGAAGACGGCATAACGAGATGCTGATTCTGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	C10	hybrid_dual_rev_034	CAAGCAGAAGACGGCATAACGAGATATAAGGCCGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	C11	hybrid_dual_rev_035	CAAGCAGAAGACGGCATAACGAGATAAGACTAGGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	C12	hybrid_dual_rev_036	CAAGCAGAAGACGGCATAACGAGATTAGCAACAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	D01	hybrid_dual_rev_037	CAAGCAGAAGACGGCATAACGAGATCAAGTAATGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	D02	hybrid_dual_rev_038	CAAGCAGAAGACGGCATAACGAGATTCATGATGGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	D03	hybrid_dual_rev_039	CAAGCAGAAGACGGCATAACGAGATTAGCTCTTGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	D04	hybrid_dual_rev_040	CAAGCAGAAGACGGCATAACGAGATACTACATGGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	D05	hybrid_dual_rev_041	CAAGCAGAAGACGGCATAACGAGATTGTGCCATGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	D06	hybrid_dual_rev_042	CAAGCAGAAGACGGCATAACGAGATGACTTCCAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT



	D07	hybrid_dual_rev_043	CAAGCAGAAGACGGCATAACGAGATTATTCTCGGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	D08	hybrid_dual_rev_044	CAAGCAGAAGACGGCATAACGAGATCACTAGATGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	D09	hybrid_dual_rev_045	CAAGCAGAAGACGGCATAACGAGATCTTCGCATGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	D10	hybrid_dual_rev_046	CAAGCAGAAGACGGCATAACGAGATCTAGCGATGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	D11	hybrid_dual_rev_047	CAAGCAGAAGACGGCATAACGAGATTACGTTAAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	D12	hybrid_dual_rev_048	CAAGCAGAAGACGGCATAACGAGATAAGTAGTCGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	E01	hybrid_dual_rev_049	CAAGCAGAAGACGGCATAACGAGATGACGTATCGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	E02	hybrid_dual_rev_050	CAAGCAGAAGACGGCATAACGAGATACGTTTCAGGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	E03	hybrid_dual_rev_051	CAAGCAGAAGACGGCATAACGAGATAGCTATATGTGACTGGAGTTCAGACGTG GCTCTTCCGAT
	E04	hybrid_dual_rev_052	CAAGCAGAAGACGGCATAACGAGATAGTGCGAAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	E05	hybrid_dual_rev_053	CAAGCAGAAGACGGCATAACGAGATCAATGGTAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	E06	hybrid_dual_rev_054	CAAGCAGAAGACGGCATAACGAGATATTCCGATGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	E07	hybrid_dual_rev_055	CAAGCAGAAGACGGCATAACGAGATAGCTGTTGGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	E08	hybrid_dual_rev_056	CAAGCAGAAGACGGCATAACGAGATTTCTGACTCGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	E09	hybrid_dual_rev_057	CAAGCAGAAGACGGCATAACGAGATCGAATGGTGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	E10	hybrid_dual_rev_058	CAAGCAGAAGACGGCATAACGAGATTAGTGTACGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	E11	hybrid_dual_rev_059	CAAGCAGAAGACGGCATAACGAGATACGCATGAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	E12	hybrid_dual_rev_060	CAAGCAGAAGACGGCATAACGAGATTTCTGTACAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	F01	hybrid_dual_rev_061	CAAGCAGAAGACGGCATAACGAGATTTCTGCTAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT




	F02	hybrid_dual_rev_062	CAAGCAGAAGACGGCATAACGAGATTGGAAGCAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	F03	hybrid_dual_rev_063	CAAGCAGAAGACGGCATAACGAGATGATGCTCAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	F04	hybrid_dual_rev_064	CAAGCAGAAGACGGCATAACGAGATAGCCTGAAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	F05	hybrid_dual_rev_065	CAAGCAGAAGACGGCATAACGAGATCACCAGTAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	F06	hybrid_dual_rev_066	CAAGCAGAAGACGGCATAACGAGATGTTACACCGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	F07	hybrid_dual_rev_067	CAAGCAGAAGACGGCATAACGAGATGCGTTAGTGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	F08	hybrid_dual_rev_068	CAAGCAGAAGACGGCATAACGAGATTGACGTTTCGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	F09	hybrid_dual_rev_069	CAAGCAGAAGACGGCATAACGAGATTAGGCGTAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	F10	hybrid_dual_rev_070	CAAGCAGAAGACGGCATAACGAGATAAGTGAGCGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	F11	hybrid_dual_rev_071	CAAGCAGAAGACGGCATAACGAGATCCGGAATAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	F12	hybrid_dual_rev_072	CAAGCAGAAGACGGCATAACGAGATAATAGGCAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	G01	hybrid_dual_rev_073	CAAGCAGAAGACGGCATAACGAGATTTTCAGGTCGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	G02	hybrid_dual_rev_074	CAAGCAGAAGACGGCATAACGAGATTAATCGCTGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	G03	hybrid_dual_rev_075	CAAGCAGAAGACGGCATAACGAGATCTTAGGTGGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	G04	hybrid_dual_rev_076	CAAGCAGAAGACGGCATAACGAGATTTGATCTCGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	G05	hybrid_dual_rev_077	CAAGCAGAAGACGGCATAACGAGATTCTTACTGGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	G06	hybrid_dual_rev_078	CAAGCAGAAGACGGCATAACGAGATGAGTACACGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	G07	hybrid_dual_rev_079	CAAGCAGAAGACGGCATAACGAGATACGTTACAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	G08	hybrid_dual_rev_080	CAAGCAGAAGACGGCATAACGAGATTGGTAAGCGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT



	G09	hybrid_dual_rev_081	CAAGCAGAAGACGGCATAACGAGATGATAGCACGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	G10	hybrid_dual_rev_082	CAAGCAGAAGACGGCATAACGAGATGACTTAGGGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	G11	hybrid_dual_rev_083	CAAGCAGAAGACGGCATAACGAGATAGTCGTTAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	G12	hybrid_dual_rev_084	CAAGCAGAAGACGGCATAACGAGATCGTGATTAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	H01	hybrid_dual_rev_085	CAAGCAGAAGACGGCATAACGAGATTCAGTAGAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	H02	hybrid_dual_rev_086	CAAGCAGAAGACGGCATAACGAGATATGTGTCTGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	H03	hybrid_dual_rev_087	CAAGCAGAAGACGGCATAACGAGATCAGAATATGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	H04	hybrid_dual_rev_088	CAAGCAGAAGACGGCATAACGAGATTAAGTCTGGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	H05	hybrid_dual_rev_089	CAAGCAGAAGACGGCATAACGAGATATCGACATGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	H06	hybrid_dual_rev_090	CAAGCAGAAGACGGCATAACGAGATACAGCTGTGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	H07	hybrid_dual_rev_091	CAAGCAGAAGACGGCATAACGAGATATCACAAGGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	H08	hybrid_dual_rev_092	CAAGCAGAAGACGGCATAACGAGATGTCTGAATTGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	H09	hybrid_dual_rev_093	CAAGCAGAAGACGGCATAACGAGATTGCCGATTGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	H10	hybrid_dual_rev_094	CAAGCAGAAGACGGCATAACGAGATACTAATGCGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	H11	hybrid_dual_rev_095	CAAGCAGAAGACGGCATAACGAGATCTAGTAGGGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT
	H12	hybrid_dual_rev_096	CAAGCAGAAGACGGCATAACGAGATTTAAGCCAGTGACTGGAGTTCAGACGTG TGCTCTTCCGAT

Section 1: Preparation of PCR2 Primer plates

- 1 To allow for high levels of combinations of indexing primers, we generate multiple forward primer plates (12) from a single 96-well plate. By using these 12 plates in combination with a single 96-well plate of reverse primers, we can index 1,152 samples. Additional primers can be designed to increase this multiplexing. Here we present one plate of forward and one plate of reverse primers to allow for this level of combination.
- 2 To generate forward primer plates, we shift the column to the right one column when replating to a new plate (e.g. column 1 becomes column 2 and column 12 becomes column 1). See the figure below. Be sure to uniquely label the plate so you know which primers are used in library preparation.



Primer plate

Note

We do NOT do this for reverse primers.

Section 2: Generate Primers mix

10m



- 3 Spin lyophilized Primer.
- 4 Resuspend in MGW using the provided volume for a 100 micromolar (μM) stock.
- 5 Vortex and spin down briefly.
- 6 Let the mixture sit for about 00:10:00 for complete dissolution. 10m
- 7 Prepare the stock primer pool (100 micromolar (μM)) by mixing 1 volume of P2_Px1Block3_v2_F with 1 volume of P2_Px1Frag3_v8_R.
- 8 Dilute primer pool with MGW to obtain 25 micromolar (μM) working solution: 1 volume of stock + 3 volumes of MGW.

Section 3: PCR1

- 9 Pre-set the thermocycler before starting master mix preparation using the programme below.

	A	B	C	D
	Step	Temperature ($^{\circ}\text{C}$)	Time (minute)	Number of cycles
	Step1 (preheat)	95	hold	NA
	Step2 (Activation)	95	02:00	NA
	Step3 (Denaturation)	94	00:30	x32
	Step4 (Annealing)	65	08:00	x32^ go to step 3
	Step5 (final Extension)	72	10:00	NA
	Step6	4	hold	NA



- 10 All sample plates and controls are to be kept on cold blocks - switch out cold blocks if they are getting warm.
- 11 All reagents are to be kept  On ice and the tube used to prepare the master mix should also remain  On ice .

- 12 Centrifuge sample plates at  4000 rpm, 00:01:00 .

1m



- 13 Prepare the master mix using the template (here 96 reactions) below and adjust it to the number of reactions.

	A	B	C
	Master Mix	x1	x96*1.1
	Gotaq	10	1056
	Molecular grade water	5.2	549.12
	Primer pool (P2_Px1Block3_v2_F and P2_Px1Frag3_v8_R)	0.8	84.48
	template DNA	4	

- 14 Divide the total master mix volume by 12 and aliquot the designated volume of into 12 tube strip on cold block to allow for multichannel pipetting.

- 15 Using a multichannel, pipette  16 μ L of master mix into each well of 96 well plate..



- 16 Using a multichannel, plate  4 μ L of sample into each well on the PCR1 plate.



Note

Add at least one positive control (lab strain, 4uL) and one PCR1 negative control (MGW, 4uL).

17 Centrifuge the plate at  4000 rpm, 00:01:00 .

1m



18 Place the plate in the thermocycler and start running the pre-set PCR program.



Section 4: PCR2 or barcoding PCR

19 Once the first PCR is complete, pre-set the thermocycler before starting the master mix preparation for PCR2 using the programme below.

	A	B	C	D
	Step	Temperature (°C)	Time (minute)	Number of cycles
	Step1	98	hold	NA
	Step2	98	0:30	NA
	Step3	95	0:10	x10
	Step4	62	8:0	x10
	Step5	72	0:20	x10
	Step6	72	1:0	NA
	Step7	4	hold	NA

20 Prepare the master mix for PCR2 using the template (here 96 reactions) below and adjust it to the number of reactions.

	A	B	C
	Master mix	x1	x96*1.1
	Gotaq	10	1056
	Molecular grade water	6	633.6

21 Divide the total master mix volume by 12 and aliquot the designated volume ( 140.8 µL for 96 sample run) into 12 tube strip on cold block to allow for multichannel



pipetting.

22 Using a multichannel, pipette  16 μL of master mix into each well. 

23 Spin down the forward and reverse PCR2 primer plates needed for the experiment at  4000 rpm, 00:01:00 and add  1 μL of forward primer and  1 μL of reverse primer into each well using the multichannel.   1m

Note

On your sample sheet, please denote the name of the forward and reverse primer plate used and the well of the sample to be sure you know the correct index for the sample you are genotyping.

24 In the post-PCR lab, add  2 μL of the PCR1 product to the master mix using a multichannel. 

25 Optional: shake the plate gently on a shaker for  00:01:00 to mix the master mix added to the samples.  1m

26 Centrifuge the plate at  4000 rpm, 00:01:00 . 

27 Place the plate in the thermocycler and start running the pre-set PCR program. 

Note

Check a random set of 10 wells per plate for amplification by gel electrophoresis based on your lab protocol (including the NTC and positive control).

Section 5: PCR Product Pooling and Cleanup

28 Pool  5 μL of PCR2 product into a single Eppendorf tube. If the pool is very large,  3 μL is sufficient. 



- 29 Bring Ampure Beads to Room temperature . The Bailey Lab makes their own bead solution with Spherotech beads (DIY).
- 30 Add an equal volume of beads to the pool (1X bead clean). Let it incubate for at least 00:10:00 at Room temperature . For a higher yield, place the tube on a Hula Mixer or invert it in your hand.
- 31 Place the pool on a magnet and wait until the beads have pelleted on the magnet and the solution is clear. Prepare enough fresh 80% ethanol to cover the beads for 2 washes.
- 32 Remove the supernatant and reserve it in another tube.
- 33 Wash with 80% ethanol for at least 00:00:30 . Discard the ethanol.
- 34 Repeat the above step for a second wash.
- 35 Remove the ethanol and let the beads dry for around 00:01:00 (do not let the beads dry until cracking). 1m
- 36 Elute in 51 μL of 1X low EDTA TE Buffer. If the pool is very large, use your discretion and elute in a higher volume (we eluted in 150 μL with 6 96-well plates of sample). Incubate at Room temperature for 00:10:00 . 10m
- 37 Pellet the beads on a magnet and remove the elution into another tube.
- 38 Quantify the elution with a Qubit fluorometer.

Section 6: Nanopore Library preparation

36m 30s

39

Note

We use the SQK-LSK114 Kit (Ligation Sequencing) and have adapted this protocol from the Ligation Sequencing Amplicons V14 protocol from Nanopore. We use the P2 sequencing machine, but this protocol could be adapted for Minion.

40 End Prep

Transfer 200 fmol of cleaned DNA into a PCR tube or plate. Calculation for fmol can be done on this website: <https://www.promega.com/resources/tools/biomath/> (1 pmol = 1,000 fmol). Adjust the volume to  50 µL with nuclease-free water.

40.1 Add the following to a PCR tube:

A	B
Reagent	Volume
DNA Template	50 uL
Ultra II End Prep Reaction Buffer	7 uL
Ultra II End Prep Enzyme Mix	3 uL
Total	60uL

40.2 Mix thoroughly by pipetting.

40.3 Incubate at 20 °C for 00:05:00 and 65 °C for 00:05:00 .

10m

40.4 Spin down and transfer the sample to a clean 1.5ml Eppendorf tube

40.5 Bring beads to Room temperature and add 60 µL to the DNA sample.



- 40.6 Incubate for 00:10:00 at Room temperature . For increased yield, place the tube on a Hula Mixer or invert in your hand. Meanwhile, prepare 1 mL of 80% ethanol. 10m
- 40.7 Pellet on a magnet and remove the supernatant - keep the supernatant for extra security.
- 40.8 Wash with 500 μ L of 80% ethanol. Remove after 00:00:30 . 30s
- 40.9 Repeat the above step for a total of 2 washes.
- 40.10 Let the beads dry for around 00:01:00 . Do not let the beads crack on the side of the tube. 1m
- 40.11 Elute in 61 μ L of nuclease-free water. Let it incubate at Room temperature for 00:05:00 . 5m
- 40.12 Pellet on the magnet and remove the elution.
- 40.13 Quantify the End-prepped DNA with a qubit fluorometer. This library can be stored at 4 °C for short-term storage, or at -20 °C or -80 °C for long-term storage.

41 Adapter Ligation

Thaw the Short Fragment Buffer (SFB) and the Elution Buffer (EB) at

Room temperature . Place the Ligation Adapter (LA), Ligation Buffer (LNB), and Quick T4 Ligase On ice .

- 41.1 Mix the following in a 1.5 ml Eppendorf tube, mixing 10-20 times after each addition:



A	B
Reagent	Volume
End-prepped DNA	60 uL
Ligation Buffer (LNB)	25 uL
NEBNext Quick T4 Ligase	10 uL
Ligation Adapter (LA)	5 uL
Total	100 uL

41.2 Mix well by pipetting. Incubate at Room temperature for at least 00:10:00 . A longer incubation of 00:20:00 will increase the amount of ligated adapters.



41.3 Bring the beads to Room temperature and resuspend.



41.4 Add 60 μL of beads and mix by inverting. Incubate for 00:10:00 at Room temperature . For higher yield, invert in your hand or place on a Hula mixer.

10m



41.5 Pellet on a magnet and remove the supernatant.

41.6 Wash with 250 μL of the Short Fragment buffer. Remove after around 00:00:30 .



41.7 Repeat the above step for a total of 2 washes.



41.8 Allow beads to dry on the magnet, and pipette out any residual buffer. Do not let beads dry to the point of cracking.

41.9 Elute in 25 μL of the Elution Buffer (EB). Incubate for 00:10:00 at Room temperature .



- 41.10 Pellet on a magnet and pipette the elution into a tube (ideally a screw cap tube for best storage).
- 41.11 Quantify the library with a qubit fluorometer. Be sure that you have at least 120 fmol of DNA, which is around  25 ng of DNA for a 4CAST library. This library can be stored at  4 °C for short-term storage, or at  -20 °C or  -80 °C for long-term storage.



Section 7: Nanopore sequencing

1h 15m

42 Loading the Flow Cell

- 42.1 Take the Promethion Flow Cell you plan to use out of the  4 °C refrigerator at least  00:20:00 before use. It is highly recommended to do a flow cell check before any sequencing.
- 42.2 Determine 100 - 120 fmols of library and add to an Eppendorf tube (Nanopore recommends 20 fmol for sequencing; however, this leads to underloading of the flow cell in most cases. If you have less than 100 fmols, add the entire library).  25 ng of 4CAST Library is around 120 fmols.
- 42.3 Bring up the 100 - 120 fmols of library to  32 µL with the elution buffer (EB).
- 42.4 After  00:20:00 at  Room temperature , gently insert the flow cell into one of the ports of the P2. The lights will turn blue if inserted correctly, and red if inserted incorrectly.
- 42.5 In a 1.5 Eppendorf tube, mix the following to create the flow cell priming mix:



A	B
Reagent	Volume
Flow Cell Flush (FCF)	1170 uL

A	B
Flow Cell Tether (FCT)	30 uL
Total	1200 uL

42.6 Set a P1000 pipette to 200. Open the inlet port and place the tip inside the port. Turn the dial until it shows 220-230. A small amount of buffer should enter the tip. Discard the tip and buffer.

42.7 Load  500 µL of the flow cell priming mix into the inlet port. Load by slowly twisting the dial in the counterclockwise direction. Wait  00:05:00 .

5m



42.8 While the priming mix is incubating, prepare the library:



A	B
Reagent	Volume
Sequencing Buffer (SB)	100 uL
Library Beads (LIB)	68 uL
DNA Library	32 uL
Total	200 uL

42.9 After the 5-minute incubation, add another  500 µL of flow cell priming mix to the inlet port.



42.10 Immediately add the  200 µL of the library slowly into the inlet port.



42.11 Close the inlet port and wait for the library to incubate for  00:10:00 before starting any sequencing.

10m



43 Flow cell wash

43.1 A flow cell wash should be done as soon as possible after the sequencing run is complete. 

43.2 Place the Wash Mix (WMX)  On ice , and thaw the Wash Diluent (DIL) at  Room temperature .

43.3 Vortex the Wash Diluent thoroughly. 

43.4 Make the Flow Cell Wash Mix by adding the following to a 1.5ml Eppendorf tube:

A	B
Reagent	Volume
Wash Diluent (DIL)	398 uL
Wash Mix (WMX)	2 uL

43.5 Pipette to mix and place  On ice . 

43.6 Stop the sequencing run on Minknow, and leave the flow cell in the port.

43.7 Make sure the inlet port is closed, and remove waste from port 2 or port 3.

43.8 Open the inlet port and set the P1000 to 200. Remove the air bubble by turning the wheel until it reaches 220-230, and a small volume of buffer enters the pipette tip.

43.9 Load  400 µL of the Flow Cell Wash Mix into the inlet port. Incubate for  01:00:00 . 




43.10 After 1 hour, remove the waste again from ports 2 or 3. Make sure the inlet port is closed.

43.11 Remove the air bubble if it is present after removing the waste. Then add  500 μL of Storage Buffer (S). The flow cell can now be placed back into a  4 °C refrigerator for long-term storage. Incubate the flow cell at least  Overnight before doing another flow cell check.

