



Oct 16, 2025

Version 1

IDEEL- UNC Implemented 4CAST Modification (Nextera Library Preparation Protocol) V.1

 [Frontiers in Parasitology](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.eq2ly4qdm1x9/v1>

Jonathan Juliano¹

¹University of North Carolina

IDEEL



Infectious Disease Epidemiology and Ecology Lab

University of North Carolina



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.eq2ly4qdm1x9/v1>

External link: <https://doi.org/10.3389/fpara.2026.1749169>



Protocol Citation: Jonathan Juliano 2025. IDEEL- UNC Implemented 4CAST Modification (Nextera Library Preparation Protocol). protocols.io <https://dx.doi.org/10.17504/protocols.io.eq2ly4qdm1x9/v1>

Manuscript citation:

Ali IM, Martinez BV, Tchuente VPK, Sadler JM, Gorman CC, Nsango SE, Junior VV, Ateba JM, Wanda GY, Bailey JA, Dinglasan RR, Lin JT, Juliano JJ Haplotype diversity of Plasmodium falciparum circumsporozoite protein prior to malaria vaccine implementation in Kaelé Health District, Cameroon, 2022–2023. Frontiers in Parasitology 5(). doi: [10.3389/fpara.2026.1749169](https://doi.org/10.3389/fpara.2026.1749169)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 12, 2025

Last Modified: October 16, 2025

Protocol Integer ID: 229871

Keywords: Infectious Disease Epidemiology and Ecology Lab, IDEEL, 4CAST amplicon panel, Multiplexed amplicon sequencing, Nextera Library Preparation, nextera library preparation protocol, 4cast modification of nextera library preparation, nextera library preparation, 4cast modification, protocol detail, protocol

Abstract

This protocol details the 4CAST modification of Nextera Library Preparation.

Guidelines

Overview:

This protocol is designed for use with the current 4CAST amplicon panel. The original publication of this protocol is: LaVerriere, E., Schwabl, P., Carrasquilla, M., Taylor, A. R., Johnson, Z. M., Shieh, M., et al. (2022). Design and implementation of multiplexed amplicon sequencing panels to serve genomic epidemiology of infectious disease: A malaria case study. *Mol Ecol Resour* 22, 2285–2303.

Portions of this protocol are taken verbatim from Supplementary Protocol S1 from this manuscript.

PCR assays are “Nextera-Adapted” by adding standardized, identical oligo segments onto existing primer sequences.

For example:

Non-Nextera-Adapted Primer→

Forward primer: CCATCAGGGAAATGTCCAGT

Reverse primer : TTCCTGCATGTCTTGAACA

Forward primer addition: TCGTCGGCAGCGTCAGATGTGTATAAGAGACAG

Reverse primer addition: GTCTCGTGGGCTCGGAGATGTGTATAAGAGACAG

Nextera-Adapted Primer:

Forward: TCGTCGGCAGCGTCAGATGTGTATAAGAGACAGCCATCAGGGAAATGTCCAGT

Reverse: GTCTCGTGGGCTCGGAGATGTGTATAAGAGACAGTTTCCTGCATGTCTTGAACA

PCR-based library preparation allows for a significant reduction in time to complete libraries, reduces the consumables required, and lessens the impact of library preparation on technicians. Pre-combined 96-well Unique Dual Index plates are utilized in combination with standardized barcoding worksheets to simplify the generation of final library indexing manifests.

4CAST Primers

Primer	To Order
4-CSP-F	TCGTCGGCAGCGTCAGATGTGTATAAGAGACAGTTAAGGAACAAGAAGGATAATACCA
4-CSP-R	GTCTCGTGGGCTCGGAGATGTGTATAAGAGACAGAAATGACCCAAACCGAAATG
4-TRAP-F	TCGTCGGCAGCGTCAGATGTGTATAAGAGACAGTCCAGCACATGCGAGTAAAG
4-TRAP-R	GTCTCGTGGGCTCGGAGATGTGTATAAGAGACAGAAACCCGAAAATAAGCACGA
4-SERA-F	TCGTCGGCAGCGTCAGATGTGTATAAGAGACAGTACTTTCCCTTGCCCTTGTG
4-SERA-R	GTCTCGTGGGCTCGGAGATGTGTATAAGAGACAGCACTACAGATGAATCTGCTACAGGA
4-AMA-F	TCGTCGGCAGCGTCAGATGTGTATAAGAGACAGCCATCAGGGAAATGTCCAGT
4-AMA-R	GTCTCGTGGGCTCGGAGATGTGTATAAGAGACAGTTTCCTGCATGTCTTGAACA

17 Index primers

	Beginning	Index	End	To order
	CAAGCAGAAGACGGCA TACGAGAT	CGCTCAG TTC	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATCGCTCAGTTC GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	TATCTGA CCT	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATTATCTGACCT GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	TCGGATG TCG	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATTCGGATGTCTG GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	CTTATGG AAT	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATCTTATGGAAT GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	TCCTATT GTG	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATTCCTATTGTG GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	GCGCGAT GTT	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATGCGCGATGTT GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	AGAGCAC TAG	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATAGAGCACTAG GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	TGCCTTG ATC	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATTGCCTTGATC GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	CTACTCA GTC	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATCTACTCAGTC GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	TCGTCTG ACT	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATTCGTCTGACT GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	GAACATA CGG	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATGAACATACGG GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	CCTATGA CTC	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATCCTATGACTC GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	TAATGGC AAG	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATTAATGGCAAG GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	GTGCCGC TTC	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATGTGCCGCTTC GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	CGGCAAT GGA	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATCGGCAATGGA GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	GCCGTAA CCG	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATGCCGTAACCG GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	AACCATT CTC	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATAACCATTCTC GTCTCGTGGGCTCGG



	CAAGCAGAAGACGGCA TACGAGAT	GGTTGCC TCT	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATGGTTGCCTCT GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	CTAATGA TGG	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATCTAATGATGG GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	TCGGCCT ATC	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATTCGGCCTATC GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	AGTCAAC CAT	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATAGTCAACCAT GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	GAGCGCA ATA	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATGAGCGCAATA GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	AACAAGG CGT	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATAACAAGGCGT GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	GTATGTA GAA	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATGTATGTAGAA GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	TTCTATG GTT	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATTTCTATGGTT GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	CCTCGCA ACC	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATCCTCGCAACC GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	TGGATGC TTA	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATTGGATGCTTA GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	ATGTCGT GGT	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATATGTCGTGGT GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	AGAGTGC GGC	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATAGAGTGCGGC GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	TGCCTGG TGG	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATTGCCTGGTGG GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	TGCGTGT CAC	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATTGCGTGTAC GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	CATACAC TGT	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATCATACACTGT GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	CGTATAA TCA	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATCGTATAATCAG TCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	TACGCGG CTG	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATTACGCGGCTG GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	GCGAGTT ACC	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATGCGAGTTACC GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	TACGGCC GGT	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATTACGGCCGGT GTCTCGTGGGCTCGG

	CAAGCAGAAGACGGCA TACGAGAT	GTCGATT ACA	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATGTCTGATTACA GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	CTGTCTG CAC	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATCTGTCTGCAC GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	CAGCCGA TTG	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATCAGCCGATTG GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	TGACTAC ATA	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATTGACTACATA GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	ATTGCCG AGT	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATATTGCCGAGT GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	GCCATTA GAC	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATGCCATTAGAC GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	GGCGAGA TGG	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATGGCGAGATGG GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	TGGCTCG CAG	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATTGGCTCGCAG GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	TAGAATA ACG	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATTAGAATAACG GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	TCCATGT TGC	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATTCCATGTTGC GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	TATCCAG GAC	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATTATCCAGGAC GTCTCGTGGGCTCGG
	CAAGCAGAAGACGGCA TACGAGAT	AGTGCCA CTG	GTCTCGTGGG CTCGG	CAAGCAGAAGACGGCATAACGAGATAGTGCCACTG GTCTCGTGGGCTCGG

15 Index Primers

	A	B	C	D
	Beginning	Index	End	To order
	AATGATAC GGCGACCA CCGAGATC TACAC	TCGTGGAGCG	TCGTGGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACTCG TGGAGCGT CGTCGGCA GCGTC
	AATGATAC GGCGACCA CCGAGATC TACAC	CTACAAGATA	TCGTGGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACCTA CAAGATATC



	A	B	C	D
				GTCGGCAG CGTC
	AATGATAC GGCGACCA CCGAGATC TACAC	TACGTTTCATT	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACTAC GTTCATTT CGTCGGCA GCGTC
	AATGATAC GGCGACCA CCGAGATC TACAC	TGCCTGGTGG	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACTGC CTGGTGGT CGTCGGCA GCGTC
	AATGATAC GGCGACCA CCGAGATC TACAC	TCCATCCGAG	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACTCC ATCCGAGT CGTCGGCA GCGTC
	AATGATAC GGCGACCA CCGAGATC TACAC	GTCCACTTGT	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACGTC CACTTGTT CGTCGGCA GCGTC
	AATGATAC GGCGACCA CCGAGATC TACAC	TGGAACAGTA	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACTGG AACAGTAT CGTCGGCA GCGTC
	AATGATAC GGCGACCA CCGAGATC TACAC	CCTTGTTAAT	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACCCT TGTTAATTC GTCGGCAG CGTC
	AATGATAC GGCGACCA CCGAGATC TACAC	GTTGATAGTG	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACGTT GATAGTGT CGTCGGCA GCGTC



	A	B	C	D
	AATGATAC GGCGACCA CCGAGATC TACAC	ACCAGCGACA	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACACC AGCGACAT CGTCGGCA GCGTC
	AATGATAC GGCGACCA CCGAGATC TACAC	CATACACTGT	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACCATA CACTGTTC GTCGGCAG CGTC
	AATGATAC GGCGACCA CCGAGATC TACAC	GTGTGGCGCT	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACGTG TGGCGCTT CGTCGGCA GCGTC
	AATGATAC GGCGACCA CCGAGATC TACAC	ATCACGAAGG	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACATCA CGAAGGTC GTCGGCAG CGTC
	AATGATAC GGCGACCA CCGAGATC TACAC	CGGCTCTACT	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACCGG CTCTACTT CGTCGGCA GCGTC
	AATGATAC GGCGACCA CCGAGATC TACAC	GAATGCACGA	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACGAA TGCACGAT CGTCGGCA GCGTC
	AATGATAC GGCGACCA CCGAGATC TACAC	AAGACTATAG	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACAAG ACTATAGTC GTCGGCAG CGTC
	AATGATAC GGCGACCA	TCGGCAGCAA	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC



	A	B	C	D
	CCGAGATC TACAC			TACACTCG GCAGCAAT CGTCGGCA GCGTC
	AATGATAC GGCGACCA CCGAGATC TACAC	CTAATGATGG	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACCTA ATGATGGT CGTCGGCA GCGTC
	AATGATAC GGCGACCA CCGAGATC TACAC	GGTTGCCTCT	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACGGT TGCCTCTT CGTCGGCA GCGTC
	AATGATAC GGCGACCA CCGAGATC TACAC	CGCACATGGC	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACCGC ACATGGCT CGTCGGCA GCGTC
	AATGATAC GGCGACCA CCGAGATC TACAC	GGCCTGTCCT	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACGGC CTGTCCTT CGTCGGCA GCGTC
	AATGATAC GGCGACCA CCGAGATC TACAC	CTGTGTTAGG	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACCTG TGTTAGGT CGTCGGCA GCGTC
	AATGATAC GGCGACCA CCGAGATC TACAC	TAAGGAACGT	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACTAA GGAACGTT CGTCGGCA GCGTC



	A	B	C	D
	AATGATAC GGCGACCA CCGAGATC TACAC	CTAACTGTAA	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACCTA ACTGTAAT CGTCGGCA GCGTC
	AATGATAC GGCGACCA CCGAGATC TACAC	GGCGAGATGG	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACGGC GAGATGGT CGTCGGCA GCGTC
	AATGATAC GGCGACCA CCGAGATC TACAC	AATAGAGCAA	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACAATA GAGCAATC GTCGGCAG CGTC
	AATGATAC GGCGACCA CCGAGATC TACAC	TCAATCCATT	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACTCA ATCCATTTC GTCGGCAG CGTC
	AATGATAC GGCGACCA CCGAGATC TACAC	TCGTATGCGG	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACTCG TATGCGGT CGTCGGCA GCGTC
	AATGATAC GGCGACCA CCGAGATC TACAC	TCCGACCTCG	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACTCC GACCTCGT CGTCGGCA GCGTC
	AATGATAC GGCGACCA CCGAGATC TACAC	CTTATGGAAT	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACCTT ATGGAATT CGTCGGCA GCGTC
	AATGATAC GGCGACCA	GCTTACGGAC	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC



	A	B	C	D
	CCGAGATC TACAC			TACACGCT TACGGACT CGTCGGCA GCGTC
	AATGATAC GGCGACCA CCGAGATC TACAC	GAACATACGG	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACGAA CATACGGT CGTCGGCA GCGTC
	AATGATAC GGCGACCA CCGAGATC TACAC	GTCGATTACA	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACGTC GATTACATC GTCGGCAG CGTC
	AATGATAC GGCGACCA CCGAGATC TACAC	ACTAGCCGTG	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACACT AGCCGTGT CGTCGGCA GCGTC
	AATGATAC GGCGACCA CCGAGATC TACAC	AAGTTGGTGA	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACAAG TTGGTGAT CGTCGGCA GCGTC
	AATGATAC GGCGACCA CCGAGATC TACAC	TGGCAATATT	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACTGG CAATATTTTC GTCGGCAG CGTC
	AATGATAC GGCGACCA CCGAGATC TACAC	GATCACCGCG	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACGAT CACCGCGT CGTCGGCA GCGTC



	A	B	C	D
	AATGATAC GGCGACCA CCGAGATC TACAC	TACCATCCGT	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACTAC CATCCGTT CGTCGGCA GCGTC
	AATGATAC GGCGACCA CCGAGATC TACAC	GCTGTAGGAA	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACGCT GTAGGAAT CGTCGGCA GCGTC
	AATGATAC GGCGACCA CCGAGATC TACAC	CGCACTAATG	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACCGC ACTAATGT CGTCGGCA GCGTC
	AATGATAC GGCGACCA CCGAGATC TACAC	GACAACTGAA	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACGAC AACTGAAT CGTCGGCA GCGTC
	AATGATAC GGCGACCA CCGAGATC TACAC	AGTGGTCAGG	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACAGT GGTCAGGT CGTCGGCA GCGTC
	AATGATAC GGCGACCA CCGAGATC TACAC	CCATCTCGCC	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACCCA TCTCGCCT CGTCGGCA GCGTC
	AATGATAC GGCGACCA CCGAGATC TACAC	CTGCGAGCCA	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACCTG CGAGCCAT CGTCGGCA GCGTC
	AATGATAC GGCGACCA	CGTTATTCTA	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC

	A	B	C	D
	CCGAGATC TACAC			TACACCGT TATTCTATC GTCGGCAG CGTC
	AATGATAC GGCGACCA CCGAGATC TACAC	GCAACATGGA	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACGCA ACATGGAT CGTCGGCA GCGTC
	AATGATAC GGCGACCA CCGAGATC TACAC	GTCCTGGATA	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACGTC CTGGATAT CGTCGGCA GCGTC
	AATGATAC GGCGACCA CCGAGATC TACAC	CAGTGGCACT	TCGTCGGC AGCGTC	AATGATAC GGCGACCA CCGAGATC TACACCAG TGGCACTT CGTCGGCA GCGTC



Materials

Consumables:

1. P200 pipettor and tips
2. P20 pipettor and tips
3. 25 μ L reagent reservoir (3)
4. 96-well PCR plates

Reagents:

1. 80% Ethanol
2. AmpureXP beads
3. KAPA HiFi HotStart ReadyMix 2X
4. Molecular Grade Water (MGW)
5. Pre-Combined UDI 96-well plate @ 10 μ M each index (i7 and i5 Index primers)
6. 4CAST primer

4CAST Primer Mix Preparation

- 1 Combine  10 μL of each primer (forward and reverse for CSP, AMA1, SERA2, and TRAP) at  100 micromolar (μM) into 1.5mL Eppendorf tube. 
- 1.1 Store primer cocktail at  $-20\text{ }^{\circ}\text{C}$. 
- 1.2 Volume will need to be modified if more than 96 samples are being prepared at once.
- 2 Perform a 1:1 dilution with MGW to achieve a working concentration of  6.250 micromolar (μM) of each primer. 
- 2.1 For primer mix with  10 μL of each primer, use  80 μL of primer mix and  80 μL dH_2O . 

PCR1

- 3 Calculate the amount of reagents needed with a 10% excess to account for pipette loss.
- 4 Prepare the PCR1 cocktail using the following quantities for each sample being prepared and pipet  6.5 μL of mix into each well of a 96-well PCR plate. 

A	B
REAGENT	VOLUME/REACTION
KAPA HiFi HotStart ReadyMix 2X	5.0 μL
Primer Cocktail	1.5 μL
Total MM Volume	6.5 μL

- 5 Add  4.0 μL sample genomic DNA (or negative/positive templates) to each well. 



5.1 Recommended to place controls randomly throughout the plate.

6 Place the plate on a thermocycler and program with the following conditions:



<i>Step</i>	<i># cycles</i>	<i>Temp C</i>	<i>Time</i>
Activation	1	95	3 min.
-Denaturation	30	98	20 sec.
-Annealing		57	15 sec.
-Extension		62	30 sec.
Final Extension	1	72	1 min.
hold	1	4	forever

Setup for Library Preparation PCR: Preparation of Materials

7 Generation of UDI plates is described at Step 19 at dx.doi.org/10.17504/protocols.io.e6nvwbwy7vmk/v2.

8 Thaw Pre-Combined UDI 96-well plate, Purified PCR reaction, and KAPA Master Mix.

9 Mix fresh 80% Ethanol using 200 proof ethanol (100%, anhydrous) molecular grade water.



10 **Record the combinatorial index plate used (e.g. UDI-A1, UDI-C3, UDI-E5).**

Setup for Library Preparation PCR: Preparing UDI Plate

11 **Spin UDI plate** to ensure there is no liquid on the plate sealing foil.

12 Remove foil from desired UDI wells row-by-row or column-by-column, using a razor to carefully score between wells.

- 12.1 Use the edge of the razor to lift the foil before gently peeling the foil off.
- 12.2 This method allows for easy tracking of what indexes have been added and helps limit contamination of index plates.
- 13 Dilute UDI in plate as necessary to reach a concentration of **10 micromolar (μM)** .
- 13.1 Add **40 μL** Tris with no EDTA (if possible) for dilution from UDI plates.
- 13.2 Perform dilution in the UDI plate.
- 14 Add **2.2 μL** of the appropriate UDI to each well. **No need to mix yet.**
- 14.1 Record which well of the UDI plate was used for each sample.
- 15 Seal UDI plate with a fresh piece of foil immediately after use.

Library Preparation PCR

- 16 Calculate the amount of master mix reagents needed with a 10% excess to account for pipette loss. (Use tool [here](#) for calculations)
- 17 Prepare LP-PCR plate for indexing samples with the following reaction volumes.

A	B
REAGENT	VOLUME/REACTION
KAPA HiFi HotStart ReadyMix 2X	5.0 μL
Molecular H2O	2.0 μL
Unique Dual Index	2.2 μL

A	B
Total MM Volume	9.2 µL

18 Add  3.0 µL of PCR 1 product to each well.



A	B	C	D
STEP	# CYCLES	TEMP C	TIME
Activation	1	95	1 min.
-Denaturation	15	95	15 sec.
-Annealing		55	15 sec
-Extension		72	30 sec
Final Extension	1	72	1 min.
hold	1	4	forever

19 Place the plate on a thermocycler and program with the following conditions:



Note

Sample can be frozen at  -20 °C at this point.

Sample pooling and double-sided clean-up/size selection

20m

20



Note

Equilibrate beads to  Room temperature for  00:30:00 prior to use; take Tapestation reagents out of 4°C to come to room temperature prior to use (if using).



Combine approximately  8 μL of each sample into a 2.0mL Eppendorf tube to create the sample pool.

20.1 If completing with less than 13 samples, can adjust volume of each sample accordingly (ex: use  10 μL of each sample) to achieve an adequate sample pool volume.

20.2 Mix gently and spin down.



21 Perform bead clean-up using AmpureXP beads for size selection.

21.1 Transfer  100 μL combined sample into a PCR tube and add  55 μL beads, mixing well.



21.2 Incubate at  Room temperature for  00:05:00 .

5m



21.3 Place on a magnetic rack for  00:05:00 .

5m

21.4 **KEEP** and remove supernatant from beads and transfer into a new PCR tube.

21.5 Add  20 μL of beads, mixing well and incubating at  Room temperature for  00:05:00 .

5m



21.6 Place on the magnetic rack for  00:03:00 .

3m

21.7 Discard supernatant and wash twice with 80% ethanol while still on rack.



- Use approximately  150 μL for each wash.

21.8 Elute in  30 μL  10 millimolar (mM) Tris-HCl pH 8.0, incubate for  00:02:00 , place on the magnetic rack, and allow the pellet to form (about 3 minutes).

2m





- If using stock Tris-HCl, concentration is typically 1M; therefore, dilute down to 1M 10 millimolar (mM) with 1 μ L of Tris-HCl into 99 μ L molecular grade water.

21.9 **Keep** and remove supernatant and transfer to a 1.5mL Eppendorf tube.

21.10 Add 2.5 μ L 10.0 millimolar (mM) Tris-HCl containing 0.05% Tween-20 for final suspension.



QC/library selection for sequencing

22 A full qPCR quantification is not necessary for 4CAST.

23 We recommend:

23.1 Qubit

23.2 Tapestation

24 QC can be done by the core or sequencing facility.