

Jan 15, 2019

iDisco immunolabeling in brown adipose tissue (BAT)

 In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.wqmfd6>

Seoeun Lee¹, Lori Zeltser¹

¹Columbia University

Optical Clearing of Tissue

SPARC



Seoeun Lee

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.wqmfd6>

External link: <https://idisco.info/>

Protocol Citation: Seoeun Lee, Lori Zeltser 2019. iDisco immunolabeling in brown adipose tissue (BAT). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.wqmfd6>

**Manuscript citation:**

[https://www.cell.com/cell/fulltext/S0092-8674\(14\)01297-5](https://www.cell.com/cell/fulltext/S0092-8674(14)01297-5)

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 02, 2019

Last Modified: January 15, 2019

Protocol Integer ID: 18925

Keywords: iDisco, BAT, immunostaining, confocal, brown adipose tissue, idisco protocol, using idisco protocol, adipose tissue, idisco, bat, result via confocal microscope, microscope

Abstract

This protocol describes how to stain Brown Adipose Tissue using iDisco protocol and visualize the result via confocal microscope.

Attachments



[iDisco.pdf](#)

474KB

Materials

	PBS 10X	Ambion AM9624
	Triton-X100	Sigma X100-500ML
	Tween-20	Sigma P9416-100ML
	DMSO	Fisher D128-4
	Sodium Azide	58032-100G
	Donkey Serum	Jackson ImmunoResearch 017-000-121
	Glycine	Sigma G7126-500G
	Heparin	Sigma H3393-50KU
	Methanol	Fisher A412SK-4
	Hydrogen Peroxide 30%	Sigma 216763-100ML
	DiChloroMethane (DCM)	Sigma 270997-12X100ML
	(Di)BenzylEhter (DBE)	Sigma 108014-1KG
	ParaFormAldehyde 16%	EMS 15710-S

Reagents and references

	Tubes (small samples)	Eppendorf 2ml or USA scientific SealRite 2ml (Cat. 1620-2707)
	Tubes (large samples)	Eppendorf 5ml transparent or amber
	Tube Revolver/Rotator	Thermofisher 88881001
	Rotating incubator	Scientific Industries Incubator-Genie (SKU: SI-1400)

Consumables and hardware

Before start

Starting the process (sample collection) on Tuesday is helpful to avoid doing the longer steps during weekends



Sample Collection

- 1 Anesthetize the mouse, perfuse with 20ml Saline and 4% PFA in 0.1M PB
- 2 Dissect BAT, carefully remove white adipose tissue surrounding the BAT. The two bilateral depots can be left connected or separated, as desired

Sample Collection

- 3 Fix in 4% PFA in 0.1M PB at 4°C, overnight with shaking, then RT 1h
- 4 Wash in PBS with shaking: RT 30min x 3times

Sample Pretreatment with Methanol

- 5 Dehydrate with methanol/H₂O series: 20%, 40%, 60%, 80%, 100%; at RT for 1h each with shaking
- 6 Wash further with 100% methanol for 1h at RT and then chill the sample at 4°C with shaking
- 7 Overnight incubation, with shaking, in 66% DCM / 33% methanol / 1% deionized water at RT
- 8 Wash twice in 100% methanol at RT with shaking, and then chill the sample at 4°C
- 9 Bleach in chilled fresh 5% H₂O₂ in methanol (1 volume 30% H₂O₂ to 5 volumes MeOH), overnight at 4°C with shaking
- 10 Rehydrate with methanol/H₂O series: 80%, 60%, 40%, 20%, PBS; 1h each at RT with shaking
- 11 Wash in PTx.2 (0.2% TritonX-100 in PBS) for 1h x2 at RT with shaking



Immunolabeling

- 12 Incubate samples in Permeabilization Solution (40ml PTx.2 + 10ml DMSO + 1.15g glycine), at 37°C with shaking for 2 days
- 13 Block in Blocking Solution (42ml PT2.x + 3ml Donkey Serum + 5ml DMSO) at 37°C with shaking for 2 days
- 14 Incubate with primary antibody in PTwH (0.2% Tween-20 in PBS + 0.1% of 10mg/ml Heparin stock solution in PBS)/5%DMSO/3% Donkey Serum, at 37°C with shaking for 7 days
 - *Validated primary antibodies
 - Beta 3 tubulin (Abcam ab18207, 1:400 dilution)
 - Tyrosine Hydroxylase (Millipore AB152, 1:400 dilution)
 - Pgp9.5 (Abcam ab108986, 1:400 dilution)
- 15 Wash in PTwH for 4-5 times (1-2 hours each) until the next day at RT with shaking
- 16 Incubate with secondary antibody in PTwH/3% Donkey Serum, at 37°C with shaking for 7 days
 - Centrifuge secondary antibody at 20000g for 10 minutes before use
 - Use Amber tube from this step forward
 - Wavelength longer than 555 (ex. AF555, AF647) is recommended
- 17 Wash in PTwH for 4-5 times (1-2 hours each) until the next day at RT with shaking

Clearing

- 18 Dehydrate with methanol/H₂O series: 20%, 40%, 60%, 80%, 100%; at RT for 1h each with shaking. Option to leave overnight at this point, is desired
- 19 3h incubation with shaking, in 66% DCM / 33% methanol / 1% deionized water at RT
- 20 Incubate in 100% DCM for 15 minutes twice at RT with shaking to wash the methanol
- 21 Incubate in DiBenzyl Ether (no shaking). The tube should be filled almost completely with DBE to prevent the air from oxidizing the sample. Before imaging, invert the tube a couple of times to finish mixing the solution



Imaging

- 22 The cleared sample should be transferred to a chamber of an appropriate size (we create one with a 3D printer), sealed on a coverglass with Kwik-sil epoxy (VWR). This chamber is then mounted on Attolfluor cell chamber (Thermofisher A7816) and filled with DBE
- 23 We typically use a Nikon Ti Eclipse inverted microscope for scanning confocal microscopy, but a resonance scanner is recommended for faster imaging