


May 08, 2019

## iGEM Calibration Protocol - Flow Cytometry Cell Size

DOI

[dx.doi.org/10.17504/protocols.io.2pegdje](https://dx.doi.org/10.17504/protocols.io.2pegdje)



Jacob Beal<sup>1</sup>, Cheryl Telmer<sup>1</sup>, Richard Tennant<sup>1</sup>, Paul Rutten<sup>1</sup>

<sup>1</sup>iGEM Measurement Committee

iGEM Measurement

Tech. support email: [pauljrutten@gmail.com](mailto:pauljrutten@gmail.com)



**Paul Rutten**

The University of Oxford

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.2pegdje>

**Protocol Citation:** Jacob Beal, Cheryl Telmer, Richard Tennant, Paul Rutten 2019. iGEM Calibration Protocol - Flow Cytometry Cell Size. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.2pegdje>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** In development

**Created:** May 07, 2019

**Last Modified:** August 07, 2019

**Protocol Integer ID:** 22982

**Keywords:** flow cytometry cell size conversion from forward scatter, flow cytometry cell size conversion, fluorescence calibration, based fluorescence calibration, igem calibration protocol, data with size calibration, size calibration, calibration protocol, eum, forward scatter

## Abstract

Conversion from forward scatter to E $\mu$ m is not a linear function, so data cannot be converted simply by multiplying with a scaling factor. We thus recommend use of software tools for processing data with size calibration.

This protocol can be combined with bead-based fluorescence calibration.

## Materials

### MATERIALS

 Size calibration beads **SpheroTech Catalog #PPS-6K**


- Experimental samples: transformed per your desired experimental procedure
- Flow cytometer with a forward scatter channel (typically a 488nm excitation and a 488nm/10nm emission filter). Measurements on this channel will be converted to equivalent  $\mu$ m diameter (E $\mu$ m).
- [TASBE Flow Analytics](#) - Analysis software

## Troubleshooting


### Before start

Read through this entire protocol carefully before you start your experiment and prepare any materials you may need. See the "Results" section for an example of a completed data analysis spreadsheet.

- 1 Prepare and culture experimental samples according to your desired experimental procedure.
- 2 Prepare experimental samples as needed for running through your flow cytometer.
- 3 Follow SpheroTech directions for preparation for PPS-6K bead sample.
- 4 Measure all samples in flow cytometer
  - Using the bead sample, adjust FSC channel voltage so that the top bead peak is below detector saturation.
  - Using the negative process control sample, adjust forward-scatter and side-scatter voltages to place the strong cell peak as close to the center of the detector range as possible (see example illustrated below).



- Instrument gating should be set to ensure that no cell events are discarded (see example illustrated below).



- Collect at least 10,000 events per sample.

- 5 Compute the  $E_{\mu\text{m/a.u.}}$  scaling function from bead sample using **TASBE Flow Analytics** to create a TASBE color model.
- 6 Convert experimental data to  $E_{\mu\text{m}}$  using **TASBE Flow Analytics** analysis functions.