



Oct 08, 2021

IHC AD neuropathology protocol

DOI

dx.doi.org/10.17504/protocols.io.btbmnik6

Christiana Bjorkli¹

¹The Norwegian University of Science and Technology

Sandvig lab



Christiana.bjorkli

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.btbmnik6>

Protocol Citation: Christiana Bjorkli 2021. IHC AD neuropathology protocol . **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.btbmnik6>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: March 13, 2021

Last Modified: October 08, 2021

Protocol Integer ID: 48205

Keywords: ihc ad neuropathology protocol, implanted transgenic mice, immunohistochemical processing, ad neuropathology, transgenic mice for characterization, implanted animal



Abstract

Immunohistochemical processing was conducted on tissue from chronically implanted animals, as well as non-implanted transgenic mice for characterization of AD neuropathology.



Troubleshooting



Day 1 fluorescent IHC

- 1 Heat-induced antigen retrieval on all tissue at  60 °C for 2 hours in phosphate buffer (PB).
- 2 Wash sections 3 × 10min in PB containing 0.2 % Triton X-100 (PBT+).
- 3 Block sections using 5 % normal goat serum in PBT+ for 1 hour.
- 4 Incubate sections with the primary antibodies Iba1 (1:1000; ab15690, Abcam, Cambridge, UK), McSA1 (1:1000; MediMabs, Montreal, Canada; characterization), and Aβ42 (1:1000; IBL America, Minnesota, USA; characterization) in PBT+ for 4 hours at  4 °C .

Day 1 DAB

- 5 Heat-induced antigen retrieval on all tissue at  60 °C for 2 hours in phosphate buffer (PB).
- 6 Wash in PB for 2× 10 minutes.
- 7 Permeabilize with 0.5 % Triton-X-100 in Tris-buffered saline (TBS-Tx; 50 mm Tris, 150 mm NaCl, pH 8.0) for 10 minutes.
- 8 Block with 10 % normal goat serum in TBS-Tx for 30 minutes.
- 9 Incubate with the primary antibody (AT8, 1:1000) in TBS-Tx overnight at  4 °C

Day 2 fluorescent IHC

- 10 Wash sections for 3 × 10 minutes with PBT.




- 11 Incubate sections with Alexa Fluor 546 to visualize Iba1 and McSA1, and Alexa Fluor 488 to visualize A β 42, for 2 hours at room temperature, protected from light.
- 12 Wash sections for 10 minutes with 4 ' , 6-diamidino-2-fenylindol (DAPI; 1:10 000; Sigma-Aldrich, Saint-Louis, MO, USA) and PB, followed by 3 washes for 10 minutes with PB.

Day 2 DAB

- 13 Wash sections with TBS-Tx for 3 \times 10 minutes.
- 14 Incubate with a biotinylated goat anti-mouse secondary antibody (1:500; Sigma-Aldrich, St Louis, MO, USA) in TBS-Tx for 90 minutes.
- 15 Wash sections for 3 \times 10 minutes with TBS-Tx.
- 16 Incubate with ABC (Vectastain ABC kit, Vector Laboratories, Burlingame, CA, USA) for 90 minutes.
- 17 Wash with TBS-Tx for 3 \times 10 minutes.
- 18 Wash with Tris-HCl for 2 \times 5 minutes.
- 19 Incubate tissue with 0.67 % diaminobenzidine (DAB) and 0.024 % H₂O₂ for 10 minutes before a final wash with Tris-HCl (2 \times 5 minutes).

Mounting sections

- 20 Mount tissue on cut edge frosted glass slides (VWR International, Radnor, PA, USA) with Tris-HCl and leave to dry for at least 4 hours on a  38 °C heating plate, protected from light.
- 21 Place mounted tissue in Xylene (VWR International, Radnor, PA, USA) for at least 5 minutes for defatting and to remove excess water from the tissue, and then coverslip with Entellan (VWR International, Radnor, PA, USA) containing Xylene.



22 Leave the mounted tissue to dry overnight, protected from light.