

Apr 25, 2023

IHC_cFOS+Parvalbumin

DOI

<https://dx.doi.org/10.17504/protocols.io.81wgbyb9yvpk/v1>

Bilge Buyukdemirtas¹

¹Wesleyan University



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.81wgbyb9yvpk/v1>

Protocol Citation: Bilge Buyukdemirtas 2023. IHC_cFOS+Parvalbumin. protocols.io
<https://dx.doi.org/10.17504/protocols.io.81wgbyb9yvpk/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 25, 2023



Last Modified: June 01, 2023

Protocol Integer ID: 80981

Keywords: staining of mouse brain tissue slice, mouse brain tissue slice, parvalbumin co, immunohistochemistry

Abstract

Immunohistochemistry protocol for c-FOS and Parvalbumin co-staining of mouse brain tissue slices (30 microns thick).

Guidelines

All shaker incubations are done at 200 rpm speed.



Staining

1 Antigen Retrieval

Remove PBS from wells. Add 500ul 0.3% citrate buffer into each well.

1.1 Incubate at RT on shaker for 30 min

Incubate at  65 °C for 45 min

Cool down on the bench for 20 min

2 3X PBS Washes

2.1 Remove citrate buffer from wells and replace with 500ul 1X PBS (or use a dropper pipette and fill the well halfway)

2.2 Incubate at RT on shaker for 15 min

 remove PBS and replace with new 1X PBS, repeat two more times

3 Block Non-Specific Signal

Make enough: 10% NGS (Vector Labs, #S-1000-20) in 0.3% triton-X PBS for all samples (500ul per well)

3.1 Remove PBS and replace with blocking reagent 500ul per well.

3.2 Incubate on shaker RT for 60min.

4 Primary Antibody

Make enough: mouse anti-parv (1:1000, Sigma-Aldrich, #P3088) and rabbi anti-cfos (1:1000, cell Signaling, #2250S) in 1% NGS (Vector Labs, #S-1000-20) in 0.3% triton-X PBS

4.1 Replace the blocking reagent with primary antibody mix, 500 ul per well, and incubate at

 4 °C

 Overnight

5 DAY 2



3X PBS Washes

5.1 Remove primary Ab mix from wells and replace with 500ul 1X PBS (or use a dropper pipette and fill the well halfway)

5.2 Incubate at RT on shaker for 15 min  repeat two more times

6 Secondary Antibody

Make enough: goat anti-mouse AF-555 (1:500, Invitrogen, #A-21422) and biotin goat anti-rabbit (1:1000, Vector Labs, #PK-6101) in 1% NGS (Vector Labs, #S-1000-20) in 0.3% triton-X PBS

6.1 Remove PBS from wells and replace with 500ul Secondary antibody mix. Cover with aluminum foil. Note: secondary antibodies are light sensitive, so from this point on the plate MUST be covered from light as much as possible.

Incubate at RT on shaker for 120 min.

7 3X PBS washes as described in Step 5 and 5.1 with the plate covered 

8 Tertiary Antibody (for cFOS only)

make enough: AF-488 streptavidin conjugated (1:1000, Thermo Fisher, #S11223) in 1% NGS (Vector Labs, #S-1000-20) in 0.3% triton-X PBS

8.1 Replace PBS in wells with 500uL of tertiary antibody solution, cover the plate with aluminum foil, and incubate at RT on shaker for 120 min.

9 3X PBS washes as described in step 5 and 5.1 with the plated covered 

10 Mounting with medium with and without DAPI