

Apr 10, 2023

Illumina denature and dilute

DOI

dx.doi.org/10.17504/protocols.io.6qpvr4ex2gmk/v1

Angela R.S. Kruse¹, Morad C Malek¹, Jamie Allen¹, Melissa Farrow¹, Jeff Spraggins¹

¹Vanderbilt University

VU Biomolecular Multim...

Human BioMolecular Atl...



Morad C Malek

Vanderbilt University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.6qpvr4ex2gmk/v1>

Protocol Citation: Angela R.S. Kruse, Morad C Malek, Jamie Allen, Melissa Farrow, Jeff Spraggins 2023. Illumina denature and dilute. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.6qpvr4ex2gmk/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: September 07, 2022

Last Modified: October 18, 2023

Protocol Integer ID: 69710

Keywords: illumina denature, illumina ng, preparation, ng

Abstract

This protocol describes the preparation of a library for Illumina NGS sequencing.

Materials

Consumables:

The following consumables are required to prepare DNA libraries for sequencing on the MiSeq.

Consumable	Supplier
HT1 (Hybridization Buffer), thawed and prechilled	Illumina, Provided in the MiSeq Reagent Kit
[Optional] Illumina PhiX Control	Illumina, catalog # FC-110-3001
1.0 N NaOH, molecular biology grade	General lab supplier
Tris-Cl 10 mM, pH 8.5 with 0.1% Tween 20	General lab supplier
Tris-HCl, pH 7.0	General lab supplier
[Protocol C] Low TE	Illumina, Provided in the AmpliSeq Library PLUS kit

Troubleshooting



Standard Normalization Method

- 1 Use the following steps to denature and dilute libraries that have been normalized using standard library quantification and quality control procedures recommended in the library prep documentation.

Follow the steps most appropriate for your library and the version of Illumina Reagent Kit you are using.

Loading concentration can also vary depending on library type and quantification methods.

Chemistry	Compatible Denature and Dilute Steps
MiSeq Reagent Kit v3	4 nM library—Results in a 6–20 pM loading concentration.
MiSeq Reagent Kit v2	4 nM library—Results in a 6–20 pM loading concentration. 2 nM library—Results in a 6–10 pM loading concentration.

The denaturation steps described in this guide make sure that the concentration of NaOH is not more than 0.001 (1 mM) in the final solution after diluting with HT1. Higher concentrations of NaOH in the library inhibit library hybridization to the flow cell and decrease cluster density.

Prepare Reagents

2 Prepare a Fresh Dilution of NaOH

1: Combine the following volumes in a microcentrifuge tube

- Laboratory-grade water (800 µl)
- Stock 1.0 N NaOH (200 µl)

The result is 1 ml of 0.2 N NaOH.

2: Invert the tube several times to mix.

NOTE

Use the fresh dilution within 12 hours.

3 Prepare HT1

1 Remove HT1 from -25°C to -15°C storage and thaw at room temperature.



2 Store at 2°C to 8°C until you are ready to dilute denatured libraries.

Denature a 4 nM Library

- 4 1 Combine the following volumes in a microcentrifuge tube.
 - 4 nM library (5 μ l)
 - 0.2 N NaOH (5 μ l)
- 5 Vortex briefly and then centrifuge at 280 \times g for 1 minute.
- 6 Incubate at room temperature for 5 minutes.
- 7 Add 990 μ l prechilled HT1 to the tube containing denatured library.
The result is 1 ml of a 20 pM denatured library.

Note that the denatured library concentration may vary by Illumina platform.