

Sep 11, 2025

Illumina Dual Index Amplicon Sequencing Sample Preparation 18sV4 (Piredda primer pair)

DOI

<https://dx.doi.org/10.17504/protocols.io.4r3l21423g1y/v1>

Eivind Stensrud¹, Alexander Eiler²

¹eDNA solutions AB; ²Alex private

eDNAsolutions

Tech. support phone: +46 (704) 139-481 email: eivind@ednasolutions.se



Eivind Stensrud

eDNA solutions AB

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.4r3l21423g1y/v1>

Protocol Citation: Eivind Stensrud, Alexander Eiler 2025. Illumina Dual Index Amplicon Sequencing Sample Preparation 18sV4 (Piredda primer pair) . **protocols.io** <https://dx.doi.org/10.17504/protocols.io.4r3l21423g1y/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: July 23, 2025

Last Modified: September 11, 2025

Protocol Integer ID: 223067

Keywords: illumina dual index amplicon sequencing sample preparation, pireda primer pair, pireda et al 2017 primer set, preparation of pcr product, pcr product, primer pair, amplicon, primer, pireda et al, primer set

Abstract

Preparation of PCR products for amplicon sequencing using 18sV4 Piredda et al 2017 primer set.

Guidelines

	Index 2 (i5)	Index 2 (i5) Sequence	Index 1 (i7)	Index 1 (i7) Sequence
	IIIu_N50_1F	TAG ATC GC	IIIu_N70_1R	TCG CCT TA
	IIIu_N50_2F	CTC TCT AT	IIIu_N70_2R	CTA GTA CG
	IIIu_N50_3F	TAT CCT CT	IIIu_N70_3R	TTC TGC CT
	IIIu_N50_4F	AGA GTA GA	IIIu_N70_4R	GCT CAG GA
	IIIu_N50_5F	GTA AGG AG	IIIu_N70_5R	AGG AGT CC
	IIIu_N50_6F	ACT GCA TA	IIIu_N70_6R	CAT GCC TA
	IIIu_N50_7F	AAG GAG TA	IIIu_N70_7R	GTA GAG AG
	IIIu_N50_8F	CTA AGC CT	IIIu_N70_8R	CCT CTC TG
	IIIu_N52_1F	CTT GCT TT	IIIu_N70_9R	AGC GTA GC
	IIIu_N52_2F	GG CTT CAA	IIIu_N71_0R	CAG CCT CG
	IIIu_N52_3F	AAT CG GCA	IIIu_N71_1R	TGC CTC TT
	IIIu_N52_4F	GGT TCA AA	IIIu_N71_2R	TCC TCT AC



	Illu_ N52 5F	ACT TCG AC	Illu_ N73 3R	GGT ATA AG
	Illu_ N52 6F	TGA CTT GC	Illu_ N73 5R	CAG CTA GA
	Illu_ N52 7F	TAG GAC CT	Illu_ N73 6R	CCA TAG CA
	Illu_ N52 8F	GGA GAC TT	Illu_ N73 8R	GGT ATA GC
	Illu_ N52 9F	AGG TTA CG	Illu_ N73 9R	GGT TAT GC
	Illu_ N53 0F	AAT TCG CT	Illu_ N74 0R	TAG GCA AG
	Illu_ N53 1F	TCA GCT AA	Illu_ N74 1R	TTG TCC AT
	Illu_ N53 2F	GC GAT ATG	Illu_ N74 3R	TCT AGG CA

Materials

Thermo Scientific, dNTP Mix (2 mM each): R0242

New England Biolabs, Q5 High-Fidelity DNA Polymerase: M0491LVIAL

New England Biolabs, Q5 Reaction Buffer: B9027SVIAL

Eurofins, NGSgrade Oligos (HPLC purified) Piredda Forward and reverse primers with illumina handles.

Thermo Fisher Scientific, Invitrogen, TE, pH 8.0, RNase-free: AM9858

Thermo Fisher Scientific, Invitrogen, UltraPure DNase/RNase-Free Distilled Water: 10977035

Abstract

- 1 Preparation of PCR products for amplicon sequencing using Piredda et al. 2017, 18sV4 primer set.

Introduction

- 2 Standard procedure for preparing libraries for metabarcoding using the 18sV4 Piredda primer set (Piredda et al., 2017) <https://doi.org/10.1093/femsec/fiw200>.

The primer set is regarded as near universal for protist (phytoplankton) and is the modified version of the Tareuk primer pair (**10.1111/j.1365-294X.2009.04480.x**), and amplifies the ecological important phytoplankton group of haptophytes better than the original primer set <https://doi.org/10.1093/femsec/fiw200>.

Safety warnings

- 3 The lab work should be conducted in designated rooms for DNA-extractions, pre-PCR and post-PCR.
Perform plate set-up in the pre-PCR lab, and all PCR products should only be worked with in the post-PCR.
These rooms should have a strict one-way flow, meaning no equipment or samples should be moved from post-PCR to pre-PCR.
Always include (a) positive control(s) with a known mock-community, and add (a) negative plate control(s) in addition to the negative control(s) from the DNA extraction.

Before starting

- 4 Ensure that metadata files are finalised before starting library preparations in the lab.

Necessary equipment:

Pipettes P10, P100/P200 + P1000

Pipette tips (w + w/o filters)

Eppendorf tubes

96-well plate / PCR strips w/ lids

96-well plate seal

Seal applicator

Microplate rack / PCR strip rack

Racks

Microcentrifuge

Vortex



Clean bench and equipment (ex. pipettes and racks) with 10% chlorine, wait 10 minutes, then proceed to clean with 70% Ethanol and wait additional 10 minutes.

If working in an UV-hood, to sterilize the hood, turn on the UV for 30 minutes.
Avoid using chlorine on metal, and ethanol on hood plastics.

DNA extraction

- 5 Conduct the DNA extraction, we recommend using the DNeasy PowerWater Sterivex Kit (remember to add negative DNA extraction(s)) for water samples. Follow the manufacturer recommendation with the exception of a dual elution of 100 µl volume.
(<https://www.qiagen.com/us/products/discovery-and-translational-research/dna-rna-purification/dna-purification/microbial-dna/dneasy-powerwater-sterivex-kit>)

Quantify the extracted DNA using PicoGreen.

DOI: [dx.doi.org/10.17504/protocols.io.dm6gpdb1gzp/v1](https://doi.org/10.17504/protocols.io.dm6gpdb1gzp/v1)

Primers

- 6 **Piredda primer set (V4 18S Next.For and V4 18S Next.Rev primers)**
NGSgrade Oligos (HPLC purified) (Eurofins), the primer set has attached Illumina adapters from Sinclair et al. 2015. <https://doi.org/10.1371/journal.pone.0116955> and Juttonen et al. <https://doi.org/10.1111/1462-2920.15058> incorporate random nucleotides to reduce amplification bias in the forward primer.

Resuspend the primers with TE Buffer or Nuclease-Free water to obtain a 100 µM stock solution (store in -20 °C).

To reduce risk of cross-contamination and primer degradation, we recommend to make several aliquots of 10 µM primers by diluting with either TE-buffer or Nuclease-Free Water.

Forward primer (59mer):

V4 18S Next.For (5'-3'):

ACA CTC TTT CCC TAC ACG ACG CTC TTC CGA TCT NNN NNN CCA GCA SCY GCG GT
A ATT CC

Reverse primer (48mer):

V4 18S Next.Rev (5'-3'):

AGA CGT GTG CTC TTC CGA TCT NNN NNN ACT TTC GTT CTT GAT YRA TGA

Illumina adapters and spacers are marked in bold, while the taxa specific primers are in roman.

First PCR reactions

7 Work in pre-PCR room

Calculate the volumes of each reagent needed to create the Mastermix for the number of samples you are running (+10% extra to account for pipetting lost and potential errors). Remember to run the samples in duplicates/triplicates to reduce PCR-stochasticity.

A	B	C	D	E
Reagent	Stock conc.	Mastermix conc.	Volume 1 reaction (µl)	Volume X reactions (µl)
5X Q5 Reaction Buffer	5X	1X	5	
Forward primer	10 µM	0,2 µM	0,5	
Reverse primer	10 µM	0,2 µM	0,5	
dNTPs	10 mM	200 µM	0,5	
Q5 High-Fidelity DNA Polymerase	0,25			
Nuclease-Free Water	16,25			
Template DNA			2*	
Total			25	

Table 1:

PCR cycling conditions has been optimized with an Applied Biosystems QuantStudio 3. Add up to 2 µl DNA template from environmental samples. For positive control, add 1 µl mock community template and 1 µl Nuclease-Free water. For negative control, add 2 µl Nuclease-Free water.

7.1 **Remember to wear gloves and lab coat. Change gloves whenever necessary.**

7.2 Clean the working area and equipment with 10% chlorine, wait 10 minutes before continuing.





- 7.3 Clean the working area and equipment with 70% ethanol, wait 10 minutes before continuing.
- 7.4 Briefly vortex and spin down all reagents in a microcentrifuge, except the polymerase. Mix the polymerase by pipette mixing 10 times and spin briefly down in a microcentrifuge.
- 7.5 Prepare a Mastermix according to Table 1. Remember to make 10% extra to account for pipetting loss and potential errors.
- 7.6 Vortex and spin down Mastermix in a microcentrifuge.
- 7.7 Add 23 µl Mastermix to each well according to plate layout.
- 7.8 Add 2 µl template DNA according to plate layout.
- 7.9 Add 2 µl Negative control according to plate layout.
- 7.10 Add 1 µl mock community + 1 µl PCR water according to plate layout.
- 7.11 Seal the plate, vortex and spin the plate down with a microcentrifuge.
- 7.12 Run PCR using the cycling conditions described in Table 2.
- 7.13 **PCR cycling conditions (first round of PCR)**

	A	B	C	D
	Step	CYCLES	TEMP	TIME
	Initial Denaturation		98°C	2 minutes
	Denaturation	10 Cycles*	98°C	10 seconds

	A	B	C	D
	Annealing	15 Cycles**	44°C	30 seconds
	Extension		72°C	20 seconds
	Denaturation		98°C	10 seconds
	Annealing		62°C	30 seconds
	Extension		72°C	20 seconds
	Final Extension		72°C	5 minutes
	Hold		4°C	∞

Table 2:

PCR cycling conditions have been optimized with an Applied Biosystems QuantStudio 3. **Aim to use as few cycles as possible**, we recommend to run a few test samples with the reaction conditions showed above.

*If you need to increase the number of cycles, do not increase the number of cycles at 44 degrees annealing temperature.

** If you need to increase the number of cycles, you can increase up to 30 cycles at 62 degrees annealing temperature.


PCR visualization

8 Conduct in post-PCR lab

After the first PCR reaction, pool the replicated samples and inspect the PCR products using Agarose gel electrophoresis (1%).

Successful amplification should give bands at ~530bp (See picture below).

Note that some primer dimers may occur, and are removed in the purification step.



Purification

9 Conduct in post-PCR

Perform purification with magnetic beads (Cytiva Sera-Mag Select) using 0.8x beads concentration.

DOI: [dx.doi.org/10.17504/protocols.io.n2bvj393wlk5/v1](https://doi.org/10.17504/protocols.io.n2bvj393wlk5/v1)

Second PCR reactions

- 10 A second PCR is conducted for attaching standard illumina handles (bold) and index primers

Multiplex_fwd:

AATGATACGGCGACCACCGAGA{TCTACAC}-[i5 index] ACACTCTTTCCTACACGACG

Multiplex_rev:

CAAGCAGAAGACGGCATACGAGAT-[i7 index]-
GTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT

We use 20 different forward index/barcode primers and 20 different reverse index/barcode primers.

Combining both primers (20X20) makes it possible to generate 400 unique tags to run 400 samples of each unique primer set in one final pool for sequencing.

(See Guidelines for index sequences)

Second PCR cycling conditions

- 11 Prepare Mastermix and mix with forward and reverse indices according to plate layout in the pre-PCR lab.

Bring the plate(s) with Mastermix and primers to the post-PCR lab, where samples from the 1st PCR reaction is added according to plate layout.

Calculate the volumes of each reagent needed to create the Mastermix for the number of samples you are running (+10% extra to account for pipetting lost and potential errors) (Table 3).

Note: Not necessarily to run the samples in duplicates

	A	B	C	D	E
	Reagent	Stock conc.	Mastermix conc.	Volume 1 reaction (μl)	Volume X reactions (μl)
	5xQ5 Reactio Buffer	5X	1X	4,00	
	Forward index (i5, illu-N5xx)	5μM	0,25 μM	1,00	
	Reverse index (i7, illu-N7xx)	5μM	0,25 μM	1,00	
	dNTPs	2mM	200 μM	2,00	
	Q5 HF DNA polymerase	2 U/μl	0.02 U/μl	0,20	
	Template from 1st PCR			2,00	
	Nuclease-Free water			9,80	
	Σ			20,00	

Table 3:

PCR cycling conditions has been optimized with an Applied Biosystems QuantStudio 3. Add 2 μl of purified PCR product of both environmental samples, positive and negative controls.

- 11.1 **Remember to wear gloves and lab coat. Change gloves whenever necessary.**
- 11.2 Clean the working area and equipment with 10% chlorine, wait 10 minutes before continuing.
- 11.3 Clean the working area and equipment with 70% ethanol, wait 10 minutes before continuing.

- 11.4 Briefly vortex and spin down all reagents in a microcentrifuge, except the polymerase. Mix the polymerase by pipette mixing 10 times and spin briefly down in a microcentrifuge.
- 11.5 Prepare a Mastermix according to Table 3. Remember to make 10% extra to account for pipetting loss and potential errors.
- 11.6 Vortex and spin down Mastermix in a microcentrifuge.
- 11.7 Add 16 µl Mastermix to each well according to plate layout.
- 11.8 Add 1 µl Forward and Reverse primer (unique combination per sample) to each well on the 96-well plate or PCR strip.
- 11.9 **Move to Post-PCR**
Seal the plate while transferring to Post-PCR room.
- 11.10 Add 2 µl template DNA according to plate layout.
- 11.11 Seal the plate, vortex and spin the plate down with a centrifuge.
- 11.12 Run the second PCR using the cycling conditions described in **Table 4**:
- 11.13 PCR cycling conditions (second round of PCR)

	A	B	C	D
	CYCLES	STEP	TEMP	TIME
		Initial Denaturation	98 C	30 sec
	15 cycles	Denaturation	98 C	10 sec
		Annealing	66 C	30 sec

	A	B	C	D
		Extension	72 C	30 sec
		Final Extension	72 C	2 min
		Hold	6 C	∞

Table 4:

PCR cycling conditions have been optimized with an Applied Biosystems QuantStudio 3. Number of cycles can be reduced to 10.


PCR visualization

12 Conduct in Post-PCR

After the second PCR reaction, pool the replicated samples and inspect the PCR products using Agarose gel electrophoresis (1%).

Successful amplification should give bands at ~600 bp (See picture below).

Note that some primer dimers may occur, and are removed in the purification step.



Purification

13 Conduct in post-PCR

Perform purification with magnetic beads (Cytiva Sera-Mag Select) using 0.75x beads concentration following the manufacturer recommendations.

DOI: dx.doi.org/10.17504/protocols.io.n2bvj393wlk5/v1



Quantify and pool in equimolar concentration

14 Conduct in post-PCR:

Quantify and pool the samples to create a final library using PicoGreen, if pooling samples from different marker lengths, see Sheet 3.

DOI: [dx.doi.org/10.17504/protocols.io.dm6gpdb1gzp/v1](https://doi.org/10.17504/protocols.io.dm6gpdb1gzp/v1)

Make sure to reach the requirements for illumina sequencing, with at least 20 μL of 10 nM library pool (Sheet 3).

Read more about the requirements here:

<https://ngisweden.scilifelab.se/methods/illumina-user-provided/>

14.1 Conduct in post-PCR

Calculation conducted in the Google Sheets.

Pool the PCR samples in equal DNA amount (ng) or for unequal length amplicons, in equal molecule amount (mol).

You will get one tube with a mix of all the samples in it.

To calculate the volume of each sample to be pooled (DNA amount mixing):

Use the lowest concentration sample to define the minimum amount of DNA (ng) that you have available from a single sample: the DNA concentration (ng/ μL) of the lowest concentration sample multiplied with its volume (μL). This will be your target DNA amount for each sample.

Calculate how many μL s of each sample you need to achieve the target DNA amount: divide the target DNA amount with the concentration of each sample.

Pipette into one tube the calculated volume of each sample. Aim to use the same pipette for all samples (dilute or pipette multiple times) to avoid pipette calibration errors.

Quality control

- 15 Ensure that primer dimers have been removed from the final pool by using Agarose gel electrophoresis (1%). If primer dimers are still present, conduct a new bead purification as described above, or conduct gel extraction.

If the DNA concentration is too low, <10nM, you can either conduct another bead purification and elute with a lower elution volume (recommended).

Or use a SpeedVac to further up-concentrate the DNA (Remember that SpeedVac also up-concentrates salts from the elution buffer), which may be unfavorable for



sequencing.