

Sep 15, 2025

Illumina Dual Index Amplicon Sequencing Sample Preparation COI Leray-XT primer set

DOI

<https://dx.doi.org/10.17504/protocols.io.e6nvw4z2dlmk/v1>

Eivind Stensrud¹, Alexander Eiler²

¹eDNA solutions AB; ²Alex private

eDNAsolutions

Tech. support phone: +46 (704) 139-481 email: eivind@ednasolutions.se



Eivind Stensrud

eDNA solutions AB

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.e6nvw4z2dlmk/v1>

Protocol Citation: Eivind Stensrud, Alexander Eiler 2025. Illumina Dual Index Amplicon Sequencing Sample Preparation COI Leray-XT primer set. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.e6nvw4z2dlmk/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: September 15, 2025

Last Modified: September 15, 2025

Protocol Integer ID: 227286

Keywords: leray xt primer, xt primer, preparation of pcr product, pcr product, amplicon

Abstract

Preparation of PCR products for amplicon sequencing using COI-Leray XT primers.

Guidelines

	Index 2 (i5)	Index 2 (i5) Sequence	Index 1 (i7)	Index 1 (i7) Sequence
	IIIu_N50_1F	TAG ATC GC	IIIu_N70_1R	TCG CCT TA
	IIIu_N50_2F	CTC TCT AT	IIIu_N70_2R	CTA GTA CG
	IIIu_N50_3F	TAT CCT CT	IIIu_N70_3R	TTC TGC CT
	IIIu_N50_4F	AGA GTA GA	IIIu_N70_4R	GCT CAG GA
	IIIu_N50_5F	GTA AGG AG	IIIu_N70_5R	AGG AGT CC
	IIIu_N50_6F	ACT GCA TA	IIIu_N70_6R	CAT GCC TA
	IIIu_N50_7F	AAG GAG TA	IIIu_N70_7R	GTA GAG AG
	IIIu_N50_8F	CTA AGC CT	IIIu_N70_8R	CCT CTC TG
	IIIu_N52_1F	CTT GCT TT	IIIu_N70_9R	AGC GTA GC
	IIIu_N52_2F	GG CTT CAA	IIIu_N71_0R	CAG CCT CG
	IIIu_N52_3F	AAT CG GCA	IIIu_N71_1R	TGC CTC TT
	IIIu_N52_4F	GGT TCA AA	IIIu_N71_2R	TCC TCT AC

	IIIu_N52 5F	ACT TCG AC	IIIu_N73 3R	GGT ATA AG
	IIIu_N52 6F	TGA CTT GC	IIIu_N73 5R	CAG CTA GA
	IIIu_N52 7F	TAG GAC CT	IIIu_N73 6R	CCA TAG CA
	IIIu_N52 8F	GGA GAC TT	IIIu_N73 8R	GGT ATA GC
	IIIu_N52 9F	AGG TTA CG	IIIu_N73 9R	GGT TAT GC
	IIIu_N53 0F	AAT TCG CT	IIIu_N74 0R	TAG GCA AG
	IIIu_N53 1F	TCA GCT AA	IIIu_N74 1R	TTG TCC AT
	IIIu_N53 2F	GC GAT ATG	IIIu_N74 3R	TCT AGG CA

Index primers from Sinclair et al. 2015

Materials

Thermo Scientific, dNTP Mix (2 mM each): R0242

New England Biolabs, Q5U Hot Start High-Fidelity DNA Polymerase: M0515LVIAL

New England Biolabs, Q5U Reaction Buffer: B9037SVIAL

Eurofins, NGSgrade Oligos (HPLC purified) COI-Leray XT primer pair (Wangesteen et al. 2018)

Thermo Fisher Scientific, Invitrogen, TE, pH 8.0, RNase-free: AM9858

Thermo Fisher Scientific, Invitrogen, UltraPure DNase/RNase-Free Distilled Water: 10977035

Eurofins, NGSgrade Oligos (HPLC purified), Index primers (Sinclair et al. 2015)

Abstract

- 1 Preparation of PCR products for amplicon sequencing using COI Leray-XT primer set.

Introduction

- 2 Standard procedure for metabarcoding using COI Leray-XT primers (Wangenstein et al., 2018).

The primer set is near eukaryotic universal as it successfully amplifies a wide range of eukaryotic organisms, and is favorable to use when metazoan eDNA studies should be conducted.

Preparation of PCR products for amplicon sequencing using COI Leray-XT primers (Wangenstein et al. 2018. <https://doi.org/10.7717/peerj.4705>), which incorporates extra inosine bases on variable nucleotides in the forward primer compared to the original COI Leray primer set (Leray et al. 2013. <https://doi.org/10.1186/1742-9994-10-34>). This decreases the primer bias, but prevents amplification using normal High-Fidelity polymerases.

To improve the quality of the analysis, eDNA solutions adapted the workflow to work with **New England Biolabs Q5U Hot Start High-Fidelity DNA Polymerase**, a two-step PCR workflow and lower cycling numbers.

Note: The primer pair amplifies a lot of prokaryotic DNA, and requires a greater sequencing depth than other primers.

Safety warnings

- 3 The lab work should be conducted in designated rooms for DNA-extractions, pre-PCR and post-PCR.

Perform plate set-up in the pre-PCR lab, and all PCR products should only be worked with in the post-PCR.

These rooms should have a strict one-way flow, meaning no equipment or samples should be moved from post-PCR to pre-PCR.

Always include (a) positive control(s) with a known mock-community, and add (a) negative plate control(s) in addition to the negative control(s) from the DNA extraction.

Before starting



4 Before you start in the lab, ensure that the metadata files are finalized.

Necessary equipment:

Pipettes P10, P100/P200 + P1000

Pipette tips (w + w/o filters)

Eppendorf tubes

96-well plate / PCR strips w/ lids

96-well plate seal

Seal applicator

Microplate rack / PCR strip rack

Racks

Microcentrifuge

Vortex

Clean bench and equipment (ex. pipettes and racks) with 10% chlorine, wait 10 minutes, then proceed to clean with 70% Ethanol and wait additional 10 minutes.

If working in an UV-hood, to sterilize the hood, turn on the UV for 30 minutes.

Avoid using chlorine on metal, and ethanol on hood plastics.

DNA extraction

5 Conduct in designated extraction room:

Conduct the DNA extraction, we recommend using the **Qiagen DNeasy PowerWater Sterivex Kit** (remember to add negative DNA extraction(s)). Follow the manufacturer recommendation with the exception of a dual elution of 100 µl volume.

Quantify the extracted DNA using PicoGreen:

DOI: dx.doi.org/10.17504/protocols.io.dm6gpdb1gzp/v1

Primers

6 COI Leray-XT primer set (mIColintF and jgHCO2198)

NGSgrade Oligos (HPLC purified) (Eurofins), the primer set has attached Illumina adapters from Sinclair et al. 2015. <https://doi.org/10.1371/journal.pone.0116955> and Juttonen et al. <https://doi.org/10.1111/1462-2920.15058> incorporate random nucleotides to reduce amplification bias in the forward primer.

Resuspend the primers with TE Buffer or Nuclease-Free water to obtain a 100 µM stock solution (store in -20 °C).

To reduce risk of cross-contamination and primer degradation, we recommend to make several aliquots of 10 µM primers by diluting with either TE-buffer or Nuclease-Free

Water.

Forward primer (65mer):

mIColintF (5'-3'): **ACA CTC TTT CCC TAC ACG ACG CTC TTC CGA TCT NNN NNN**

GGW ACW RGWTGR ACW [I] T [I] TAY CCY CC

Reverse primer (53mer):

igHCO2198 (5'-3'):

AGA CGT GTG CTC TTC CGA TCT NNN NNN TA [I] ACY TC [I] GGR TG [I] CCR AARAAY

CA

Illumina adapters and spacers are marked in bold, while the taxa specific primers are in roman.

First PCR reactions

7 Conduct in pre-PCR:

Calculate the volumes of each reagent needed to create the Mastermix for the number of samples you are running (+10% extra to account for pipetting lost and potential errors). Remember to run the samples in duplicates/triplicates to reduce PCR-stochasticity.

	A	B	C	D	E
	Reagents	Stock conc.	Mastermix conc.	Volume Per reaction (µL)	Total volume X reactions (µl)
	5X Q5U Reaction Buffer		1X	5	
	Forward Primer	10 µM	0,5 µM	1,25	
	Reverse Primer	10 µM	0,5 µM	1,25	
	dNTP Mix	2 mM each	40 µM each	0,5	
	Nuclease-Free Water			12,75	
	Q5U Hot Start High-Fidelity DNA Polymerase		0,02 U/µl	0,25	
	Template DNA			4	
	Total			25	

Table 1:

PCR cycling conditions has been optimized with an Applied Biosystems QuantStudio 3. Add up to 4 μ l DNA template from environmental samples. For positive control, add 1 μ l mock community template and 3 μ l Nuclease-Free water. For negative control, add 4 μ l Nuclease-Free water.

7.1 Remember to wear gloves and lab coat. Change gloves whenever necessary.

7.2 Clean the working area and equipment with 10% chlorine, wait 10 minutes before continuing.



7.3 Clean the working area and equipment with 70% ethanol, wait 10 minutes before continuing.



7.4 Briefly vortex and spin down all reagents in a microcentrifuge, except the polymerase. Mix the polymerase by pipette mixing 10 times and spin briefly down in a microcentrifuge.



7.5 Prepare a Mastermix according to Table 1. Remember to make 10% extra to account for pipetting loss and potential errors.

7.6 Vortex and spin down Mastermix in a microcentrifuge.

7.7 Add 21 μ l Mastermix to each well according to plate layout.

7.8 Add 4 μ l template DNA according to plate layout.

7.9 Add 4 μ l Negative control (Nuclease-Free Water) according to plate layout.

7.10 Add 1 μ l mock community + 3 μ l Nuclease-Free Water according to plate layout

7.11 Seal the plate, vortex and spin the plate down with a centrifuge.



7.12 Run PCR using the cycling conditions described in the Table 2.



7.13 PCR cycling conditions (first round of PCR)

	A	B	C	D
		STEP	TEMP	TIME
		Initial Denaturation	98°C	2 minutes
	25 Cycles*	Denaturation	98°C	10 seconds
		Annealing	45°C	30 seconds
		Extension	72°C	20 seconds
		Final extension	72°C	5 minutes
		Hold	4°C	


Table 2:
PCR cycling conditions have been optimized with an Applied Biosystems QuantStudio 3. Aim to use as few cycles as possible, we recommend to run a few test samples with 25 cycles.

PCR visualization

8 Conduct in post-PCR:

After the first PCR reaction, pool the replicated samples and inspect the PCR products using Agarose gel electrophoresis (1%).

Successful amplification should give bands at ~450 bp (See picture below).



Note: We rarely see primer dimers here, so expect a two bands a clear band of around 450 bp and a band at 50-60 bp for the residual primers

Purification

- 9 **Conduct in post-PCR:**
Perform purification with magnetic beads (Cytiva Sera-Mag Select) using 0.8x beads concentration following the manufacturer recommendations.
DOI: dx.doi.org/10.17504/protocols.io.n2bvj393wlk5/v1

Second PCR reactions

- 10 A second PCR is conducted for attaching standard illumina handles (bold) and index primers
- Multiplex_fwd:**
AATGATACGGCGACCACCGAGA{TCTACAC}-[i5 index] ACACTCTTCCCTACACGACG
- Multiplex_rev:**
CAAGCAGAAGACGGCATACGAGAT-[i7 index]-
 GTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT

We use 20 different forward index/barcode primers and 20 different reverse index/barcode primers.

Combining both primers (20X20) makes it possible to generate 400 unique tags to run 400 samples of each unique primer set in one final pool for sequencing.

(See Guidelines for index sequences)

Second PCR cycling conditions

11 Conduct in post-PCR:



Prepare Mastermix and mix with forward and reverse indices according to plate layout in the pre-PCR lab.

Bring the plate(s) with Mastermix and primers to the post-PCR lab, where samples from the 1st PCR reaction is added according to plate layout.

Calculate the volumes of each reagent needed to create the Mastermix for the number of samples you are running (+10% extra to account for pipetting lost and potential errors).

	A	B	C	D	E
	Reagents	Stock conc.	Mastermix conc.	Volume 1 reaction (µl)	Volume X reactions (µl)
	5xQ5U Reaction Buffer	5X	1X	4,00	
	Forward index (i5, illu-N501-N508)	5µM	0,25 µM	1,00	
	Reverse index (i7, illu-N701-N712)	5µM	0,25 µM	1,00	

	A	B	C	D	E
	dNTP mix	2mM each	200 μ M each	2,00	
	Q5U HF DNA polymerase	2 U/ μ l	0.02 U/ μ l	0,20	
	Template from 1st PCR			2,00	
	Nuclease-Free water			9,80	
	Σ			20,00	

Table 3:
 PCR cycling- and reaction conditions has been optimized with an Applied Biosystems QuantStudio 3.
 Add 2 μ l of purified PCR product of both environmental samples, positive and negative controls.

- 11.1 **Remember to wear gloves and lab coat. Change gloves whenever necessary.**
- 11.2 Clean the working area and equipment with 10% chlorine, wait 10 minutes before continuing.
- 11.3 Clean the working area and equipment with 70% ethanol, wait 10 minutes before continuing.
- 11.4 Briefly vortex and spin down all reagents in a microcentrifuge, except the polymerase. Mix the polymerase by pipette mixing 10 times and spin briefly down in a microcentrifuge.
- 11.5 Prepare a Mastermix according to Table 1. Remember to make 10% extra to account for pipetting loss and potential errors.
- 11.6 Vortex and spin down Mastermix in a microcentrifuge.



- 11.7 Add 16 µl Mastermix to each well according to plate layout.
- 11.8 Add 1 µl Forward and Reverse primer (unique combination per sample) to each well on the 96-well plate.
- 11.9 **Move to Post-PCR**
Seal the plate while transferring to Post-PCR room.
- 11.10 Add 2 µl template DNA according to plate layout.
- 11.11 Seal the plate, vortex and spin the plate down with a centrifuge.
- 11.12 Run the second PCR using the cycling conditions described in the Table 4:

11.13 **PCR cycling conditions (second round of PCR)**

STEP	TEMP.	TIME
Initial Denaturation	98 C	2 min
	98 C	10 sec
15 cycles	66 C	30 sec
	72 C	30 sec
Final Extension	72 C	2 min
Hold	6 C	∞

PCR cycling conditions have been optimized with an Applied Biosystems QuantStudio 3. Number of cycles can be reduced to 10.


PCR visualization

12 Conduct in Post-PCR:

After the second PCR reaction, inspect the PCR products using Agarose gel electrophoresis (1%).

Successful amplification should give bands at 550 bp (See picture below).

Note that some primer dimers ~240 bp are expected, and is removed in the purification step.



Purification

13 Conduct in post-PCR:



Perform purification with magnetic beads (Cytiva Sera-Mag Select) using 0.75x beads concentration following the manufacturer recommendations.

DOI: dx.doi.org/10.17504/protocols.io.n2bvj393wlk5/v1

Quantify and pool in equimolar concentration

14 Conduct in post-PCR:

Quantify and pool the samples to create a final library using PicoGreen, if pooling samples from different marker lengths, see Sheet 3.

DOI: dx.doi.org/10.17504/protocols.io.dm6gpdb1gzp/v1

Make sure to reach the requirements for illumina sequencing, with at least 20 μL of 10 nM library pool (Sheet 3).

Read more about the requirements here:

<https://ngisweden.scilifelab.se/methods/illumina-user-provided/>

14.1 Conduct in post-PCR:

Calculation conducted in the Google Sheets.

Pool the PCR samples in equal DNA amount (ng) or for unequal length amplicons, in equal molecule amount (mol).

You will get one tube with a mix of all the samples in it.

To calculate the volume of each sample to be pooled (DNA amount mixing):

Use the lowest concentration sample to define the minimum amount of DNA (ng) that you have available from a single sample: the DNA concentration (ng/ μL) of the lowest concentration sample multiplied with its volume (μL). This will be your target DNA amount for each sample.

Calculate how many μL s of each sample you need to achieve the target DNA amount: divide the target DNA amount with the concentration of each sample.

Pipette into one tube the calculated volume of each sample. Aim to use the same pipette for all samples (dilute or pipette multiple times) to avoid pipette calibration errors.

Quality control

15 Conduct in post-PCR:

Ensure that primer dimers have been removed from the final pool by using Agarose gel electrophoresis (1%). If primer dimers are still present, conduct a new bead purification as described above, or conduct gel extraction.



If the DNA concentration is too low, <10nM, you can either conduct another bead purification and elute with a lower elution volume (recommended).
Or use a SpeedVac to further up-concentrate the DNA (Remember that SpeedVac also up-concentrates salts, which may be unfavorable for sequencing).