

Dec 14, 2023

Image analysis of plasma membrane contacts

DOI

dx.doi.org/10.17504/protocols.io.n2bvj364plk5/v1

Chase Amos¹, Pietro De Camilli¹

¹Departments of Neuroscience and of Cell Biology, Howard Hughes Medical Institute, Program in Cellular Neuroscience, Neurodegeneration and Repair, Yale University School of Medicine, New Haven, Connecticut 06510, USA 2Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815



Chase Amos

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n2bvj364plk5/v1>

Protocol Citation: Chase Amos, Pietro De Camilli 2023. Image analysis of plasma membrane contacts. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.n2bvj364plk5/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: August 01, 2023

Last Modified: June 01, 2024

Protocol Integer ID: 85739

Keywords: ASAPCRN, plasma membrane contact, plasma membrane contact site, plasma membrane, k562 cell, membrane, cell, plasma, image analysis

Funders Acknowledgements:

ASAP

Grant ID: ASAP-000580

Abstract

This protocol details the image analysis of overexpressed VPS13A^{Halo} at plasma membrane contact sites in K562 cells.

Troubleshooting

Image analysis of plasma membrane contacts

- 1 Using the image analysis program FIJI, split the channels of the image containing the mitochondria (TMRE stain or overexpressed mito-BFP) and overexpressed VPS13A^{Halo}.
- 2 Threshold the mitochondria channel to include the mitochondria. Convert to a 32-bit image and convert the threshold to binary, setting values within the mitochondria to be NaN (not a number). This step generates a negative mitochondrial mask.

Note

Due to the close proximity of mitochondria to the cell edge in K562 cells, a negative mask is used to exclude signal originating from the mitochondria.

- 3 Multiply the VPS13A^{Halo} channel and the negative mitochondrial mask into a 32-bit result image. This generates a mitochondria-subtracted VPS13A^{Halo} image.
- 4 Using the mitochondria-subtracted VPS13A^{Halo} image, apply a threshold to include signal present at the plasma membrane. Convert to a 32-bit image and convert threshold to binary, setting values outside of the VPS13A signal to be NaN.

The resulting image includes plasma membrane and excludes signal originating from the mitochondria or background.

- 5 Trace the cell edge using the selection brush tool with a pixel size of 1. Apply the ROI to the image generated in step 4 and measure the "% area" (selected as "Area fraction" in the "Set measurement" window) of the cell edge length. The recorded value corresponds to the percent of plasma membrane length that contains VPS13A^{Halo} signal.