

Jan 22, 2025

Version 2

Imaging and quantification V.2

DOI

<https://dx.doi.org/10.17504/protocols.io.yxmvm31mbl3p/v2>

Tony Hsiao^{1,2}, Felicia Xaveria Suteja^{1,2}, Barkha Kishor Goswami^{1,2}, Anahid Ansari Mahabadian^{1,2}, YuHong Fu^{1,2}, Glenda Halliday^{1,2}

¹Brain and Mind Centre & Faculty of Medicine and Health School of Medical Sciences, The University of Sydney, Sydney, NSW 2050, Australia;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815



Tony Hsiao

Usyd

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.yxmvm31mbl3p/v2>

Protocol Citation: Tony Hsiao, Felicia Xaveria Suteja, Barkha Kishor Goswami, Anahid Ansari Mahabadian, YuHong Fu, Glenda Halliday 2025. Imaging and quantification. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.yxmvm31mbl3p/v2> Version created by Tony Hsiao

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 17, 2024

Last Modified: January 22, 2025

Protocol Integer ID: 115713

Keywords: ASAPCRN, using olympus vs200 scanner, olympus vs200 scanner, quantification imaging, quantification by qupath, qupath, imaging

Funders Acknowledgements:

igning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase,
Grant ID: 20815

Abstract

Imaging using Olympus VS200 scanner and quantification by Qupath ver 3.8

Cropping (OlyVIA ver 3.8) and Quantification (QuPath ver 0.4.3)

- 1 Cropping and section selection
 - 1.1 Images are cropped into groups following the Bregma sections (Bregma -2.69mm to -4.03mm):
 - 1.2 a. Rostral
 - 1.3 b. Intermediate
 - 1.4 c. Caudal
 - 2 Channel separation and Annotation
 - 2.1 Channels are renamed based on their antibody (DAPI, TH, NeuN, mCherry, α Syn)
 - 2.2 Annotations (categorized as SNC and VTA) is created surrounding TH+ cells clusters.
- 3 Cell detection
 - 3.1 Cell detection is applied with the parameters as followed:
 - a. Detection channel: DAPI/TH/NeuN/mCherry/ α Syn (depending on what channel we would like to detect)
 - b. Background radius: 8 μ m
 - c. Median filter radius: 0 μ m
 - d. Median filter radius: 0 μ m
 - e. Sigma: 1.5 to 2 μ m (depending on the size of the targeted stain)
 - f. Minimum area: X μ m² (depending on the size of the targeted stain)
 - g. Maximum area: 1000 μ m² (depending on the size of the targeted stain)

- h. Threshold: X (depending on the average targeted stain intensity on each detection)

3.2 Single Measurement Classifier

1. Add class for the targeted channels of the classifier:
 - a. TH+
 - b. TH-
 - c. NeuN+
 - d. NeuN-
 - e. mCherry+
 - f. mCherry-
 - g. α Syn+
 - h. α Syn-

3.3 After the necessary detections are made, a new 'Single measurement classifier' is created with the parameters as followed:

- a. Object filler: Detections (all)
- b. Channel filter: TH/NeuN/mCherry/ α Syn (depending on the targeted stain, but was not detected using 'cell detection')
- c. Measurement: Cell/Nucleus/Cytoplasm: DAPI/TH/NeuN/mCherry/ α Syn mean (depending on where the targeted channel is concentrated)
Threshold: X (depending on the average targeted stain intensity on each detection)
- d. Above Threshold: TH+/NeuN+/mCherry+/ α Syn+
- e. Below Threshold: TH-/NeuN-/mCherry-/ α Syn-

3.4 The new classifier was saved after given a name according to the detected channel and threshold.

3.5 Export measurement

- a. 'Annotation measurement' results of each image are then exported to excel and arranged accordingly.

3.6 DAPI count

- a. Apply 'cell detection' using the channel: DAPI

3.7 NeuN+ with TH+/TH- and α Syn+/ α Syn- count

- a. Apply 'cell detection' using the channel: NeuN
- b. Apply 'single measurement classifier' using the channel: TH
- c. Apply 'single measurement classifier' using the channel: α Syn