



Dec 11, 2025

## Imaging using 2-photon microscope through GRIN lens

DOI

<https://dx.doi.org/10.17504/protocols.io.4r3l21d2jg1y/v1>

Ines Rodrigues-Vaz<sup>1,2,3</sup>

<sup>1</sup>Zuckerman Mind Brain Behavior Institute, Departments of Neuroscience and Neurology, Columbia University;

<sup>2</sup>Champalimaud Neuroscience Programme, Champalimaud Centre for the Unknown;

<sup>3</sup>Aligning Science Across Parkinson's (ASAP) Collaborative Research Network



Ines Rodrigues-Vaz

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.4r3l21d2jg1y/v1>

**Protocol Citation:** Ines Rodrigues-Vaz 2025. Imaging using 2-photon microscope through GRIN lens. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.4r3l21d2jg1y/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working



**Created:** December 02, 2025

**Last Modified:** December 11, 2025

**Protocol Integer ID:** 234019

**Keywords:** ASAPCRN, photon microscope through grin lens, photon imaging with bruker microscope, photon microscope, photon imaging, bruker microscope, grin lens, imaging, photon

**Funders Acknowledgements:**

**Aligning Science Across Parkinson's**

Grant ID: ASAP-020551

## Abstract

This is the standard protocol for 2-photon imaging with Bruker microscope through GRIN lens used for Rodrigues-Vaz and Athalye et al, 2025.

## Materials

- Behavior setups
- 70% ethanol
- 10% ethanol
- 2-photon laser scanning microscope (Ultima In Vivo, Bruker)
- Femtosecond Ti:sapphire laser (Chameleon Vision-S, Coherent)
- Low repetition photostimulation laser (Monaco 1035-40, Coherent)
- Lens paper
- Distilled water

## Safety warnings

! -Wear appropriate PPE as required by your institution.

## Ethics statement

This protocol was approved by Columbia University IACUC. Please do not perform any of these procedures unless there is prior approval from the institution's animal ethics committee.

## Before start

- Mice were implanted with head plate and GRIN lens and habituated to being head fixed.
- Behavior setup was prepared and placed under microscope for *in vivo* experiments.
- 2-photon laser and Pockels cell controller were turned on.
- If doing holographic stimulation, calibration step must have been performed before. See "Holographic calibration" protocol for this.

## General preparation

- 1 Confirm laser is aligned with target at objective output.
- 1.1 If doing holographic stimulation, confirm that stimulation laser is also aligned
- 2 Head fix mouse in behavior setup under microscope.
- 3 Find top of GRIN lens using 4X Nikon objective. If necessary, clean top of the lens with wet cotton tip.
- 4 Swap objectives – Use 16X water-immersion Nikon objective – and find top of the GRIN lens again.
- 5 Turn off all lights in the room and switch to using 2-photon microscope for imaging.
- 6 Find desired field of view (FOV). If matching same FOV across multiple days follow steps below:
  - a. Load BOT file (brightness over time) that labeled ROIs from previous session – these ROIs were extracted using Suite2p with custom trained classifier to identify medium spiny neurons (MSNs).
  - b. Find XYZ coordinates that better match last session's FOV
  - c. To continuously monitor the imaging power, we used a beam pickoff (Thorlabs BSF10-B) and power meter (Thorlabs PDA100A2). Pockels cell bias was adjusted accordingly to maintain imaging power throughout session.
- 7 Record 1000 frames of functional and structural images – green and red PMTs, respectively. These are used to identify cell types.
- 8 Turn off red PMTs and red acquisition window.
- 9 Turn on imaging acquisition for 60000 frames (or for as long as behavior session lasts) followed by all behavior boards.



- 10 Record session, confirming that all TTLs are being recorded in Prairie's voltage recording as expected by design. 
- 11 At the end of the session, turn off all behavior boards and then the imaging recording.
- 12 Acquire another 1000 frames of functional and structural images – turn on red PMTs and acquisition window to do so.
- 13 Collect z-stack 50 um above and below imaged FOV, 2 um between each step. Each image is an average of at least 16 frames.
- 14 Turn off PMTs, rise objective and remove mouse from head-fixation.
- 15 Convert data using "Image Ripping" software. 
- 16 Transfer data from local computer to storage locations. 

## End of day procedure

- 17 Clean objectives with wet lens paper.
- 18 Turn off Pockels cell controller and laser.
- 19 Turn on lights in the room and take the mice to their holding room.