



Sep 11, 2025

Immunohistochemistry (Stadelmann Lab - single antibody, DAB-POX-Avidin-Biotin)

DOI

<https://dx.doi.org/10.17504/protocols.io.8epv5rnw6g1b/v1>

Jonas Franz¹, Katja Schulz¹, Christine Stadelmann¹

¹University Medical Center Göttingen, Neuropathology



Jonas Franz

University Medical Center Göttingen, Neuropathology

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.8epv5rnw6g1b/v1>

Protocol Citation: Jonas Franz, Katja Schulz, Christine Stadelmann 2025. Immunohistochemistry (Stadelmann Lab - single antibody, DAB-POX-Avidin-Biotin). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.8epv5rnw6g1b/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

This protocol is established over years and provides robust results.



Created: June 27, 2024

Last Modified: September 11, 2025

Protocol Integer ID: 102505

Keywords: ASAPCRN, primary antibody, single antibody, secondary antibody, coupled secondary antibody, using avidin, immunohistochemistry, avidin, biotin, coupled peroxidase, stadelmann lab

Abstract

This protocol describes staining of any primary antibody using a biotin coupled secondary antibody using Avidin coupled peroxidase.



Paraffin removal and rehydration

14h 35m

- 1 Immunohistochemistry was performed on 2-3 μm thick paraffin sections.
- 2 Remove Excess Paraffin at 56° for 1h-24h 1h
- 3 Rehydration
- 3.1 4x Xylol à 5 min 20m
- 3.2 1x Isopropanol:Xylol 1:1 à 3 min 3m
- 3.3 2× 100% Isopropanol à 3 min 6m
- 3.4 1× 90% Isopropanol à 3 min 3m
- 3.5 1× 70% Isopropanol à 3 min 3m
- 3.6 1× 50% Isopropanol à 3 min 3m
- 3.7 3x wash in aqua dest 1m 30s

Antigen Retrieval

14h 35m

- 4 Pretreatment
- 4.1 heat antigen retrieval (1x citrate buffer 10 mM , pH 6) - 30 min steamer



adopt according to tissue and antibody

4.2 3x wash PBS

4.3 formic acid (98%) (for Amyloid beta immunohistochemistry) - 1minute at RT

4.4 3x wash PBS

4.5 3% hydrogen peroxide in PBS 10 min @ RT

4.6 3x wash PBS

Blocking

14h 35m

5 Blocking in 10% normal goat serum in PBS for 10 min @ RT

10m

6 Decant 10% NGS in PBS and move directly to the next step without washing

Staining

14h 35m

7 Incubation of primary antibodies over night with respective dilution in 10% NGS @RT in a humid chamber

12h



	target antigen - human	host species	Clone	Company	Catalog Number	Dilution
	pTau	mouse	at8	Thermo Fisher Scientific	MN1020	1:100

Selected Antibodies

8 3x wash with PBS

1m 30s



9 secondary antibody

1h

Incubation with secondary antibody coupled to biotin at RT for 1h in a humid chamber

	target antigen	host species	Clone	Company	Catalog Number	Dilution
	mouse IgG	sheep		GE Healthcare	RPN1001	1:200

s

10 3x wash PBS

1m 30s

10.1 ExtrAvidin-peroxidase 1:1000 in 10% NGS for 1h in humid chamber at RT

1h

	target	coupled enzyme	Company	Catalog Number	Dilution
	biotin	Peroxidase			1:1000

Counterstaining + Mounting

14h 35m

11 3x wash PBS

1m 30s

12 Slides were developed using DAB (1gin40mL PBS → 1mL)
1 mL Aliquots + 50 mL PBS + 20 µL 30% H2O2)

10m

13 3x PBS

1m 30s

or optional

3x TBS → move to AP coupled secondary staining step 9
([dx.doi.org/10.17504/protocols.io.3byl4952rgo5/v1](https://doi.org/10.17504/protocols.io.3byl4952rgo5/v1))

14 10 min tap water

10m

15 Hematoxylin 30s



30s

16 Differentiation dip in HCL

5s

17 Dehydration reverse see **step 3**

40m

18 Entellan mounting with coverslip