




Nov 16, 2023

Immuno Fluorescence staining of human FFPE (formalin fixed paraffin embedded) gut mucosal biopsy

DOI

<https://dx.doi.org/10.17504/protocols.io.n2bvj3jrnlk5/v1>



Ran RZ Zhou¹

¹University of Chicago

Human Cell Atlas Metho...

Helmsley project_Basu lab



Ran RZ Zhou

University of Chicago

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n2bvj3jrnlk5/v1>

Protocol Citation: Ran RZ Zhou 2023. Immuno Fluorescence staining of human FFPE (formalin fixed paraffin embedded) gut mucosal biopsy. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.n2bvj3jrnlk5/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

I use this protocol and it's working

Created: November 16, 2023

Last Modified: November 16, 2023

Protocol Integer ID: 91039

Keywords: gut, FFPE, IF, fluorescence, staining, gut mucosal biopsy this protocol, human gut mucosal biopsy, gut mucosal biopsy, detail instruction on the immuno fluorescence, immuno fluorescence, human ffpe, fluorescence, fixation for the tissue, buffered formalin

Disclaimer

This protocol is for research only.

Abstract

This protocol is to provide a detail instruction on the immuno fluorescence staining on human gut mucosal biopsy that has been preserved in FFPE and sectioned at 5 micron thickness. The fixation for the tissues was 24 hour at room temperature in 10% neutral buffered formalin. Target retrieval needs optimization is fixation, sectioning conditions change.

Guidelines

The human gut mucosal tissues was retrieved by colonoscopy. The size of tissues range from $1 \times 1 \times 1$ mm, to $2 \times 4 \times 4$ mm.

Smaller tissues were embedded in 1% agarose gel before FFPE.



Materials

1x Citrate buffer 200 ml
20 ml 10x Citrate buffer
180 ml distilled water

1x PBST 1000 ml
100 ml 10x PBS
0.5ml 10% Tween20
900ml distilled water

Permeabilization buffer 5 ml
125ul 10% TritonX-100
3.3 1x PBST
0.5ml 10% BSA

Blocking buffer 5 ml
2.5ml 10% BSA
2.5ml 1x PBST

Histoclear II (National diagnostics HS-202)
Histology grade ethanol (VWR 89370-084)
10xCitrate buffer (MilliepporeSigma C9999)
10xPBS (FisherSci BP399500)
10%Tween20 (Teknova T0710)
10%TritonX (MilliepporeSigma 93443)
10%BSA (Miltenyi Biotec 130-091-376)

Safety warnings

 Boiling buffer, microwave and hotplate are used in the protocol. Please proceed with caution and use thermal gloves at these steps.

Ethics statement

The human intestinal tissue are obtained after patients' consents and approval from Institutional Review Board at the University of Chicago (IRB Number: 15573A). All the samples are processed for research use only.

Before start

Check the tissue block quality with H&E histology.



Deparaffinization

15m

- 1 Immerse tissue sections on slides in Histoclear II in a staining dish Room temperature 00:05:00 . Move the slide holder up and down x 5 times.
Repeat this step two more times with fresh Histoclear II. There should be no visible paraffine remaining on the slide by the end of this step.

Rehydration

16m

- 2 Move to another staining dish. Shake off Histoclear II from former step. Immerse tissue sections on slides in histology ethanol in a staining dish Room temperature 00:02:00 . Repeat this step one more times with fresh ethanol.
- 3 Move to another staining dish. Shake off ethanol from former step. Immerse tissue sections on slides in 70% histology ethanol in a staining dish Room temperature 00:02:00 . Repeat this step one more times with fresh 70% ethanol.
- 3.1 Prepare fresh 1 x Citra buffer 200 ml in a glass beaker, and place 200 ml distilled water in the second glass beaker.
- 4 Move to another staining dish. Shake off ethanol from former step. Immerse tissue sections on slides in 50% ethanol in a staining dish Room temperature 00:02:00 . Repeat this step one more times with fresh 50% ethanol.
- 4.1 Microwave the liquid in the beakers x 1 minute. Place the beakers with liquid on the hot plate and bring the liquid to the boiling point.
- 5 Move to another staining dish. Immerse tissue sections on slides in distilled water in a staining dish Room temperature 00:02:00 . Repeat this step one more times with fresh water.

Target retrieval

26m 30s



6 Place the slides to a slide basket. Immerse the tissue sections in the boiling water

🌡️ 100 °C .

30s

7 Move the slides from boiling water. Immerse the tissue sections in boiling 1 x citrate buffer 🌡️ 100 °C . There should be no agarose gel remaining on the slides.

15m

8 Move the slides from citrate buffer. Cool the tissues in distilled water dish

🌡️ Room temperature .

1m

8.1 Remove all the beakers from the hot plate.

9 Airdry the slides 🌡️ Room temperature .

5m

10 Draw a hydrophobic barrier around the tissue sections on the slides.

5m

10.1 Prepare the humidity chamber with the heated water from step 13.

Permeabilization

22m

11 Place the slides in a humidity chamber. Add 1 x PBS to the tissues

🌡️ Room temperature .

2m

12 Tap off water. Add permeabilization buffer to the tissues 🌡️ Room temperature

🕒 00:10:00 . Repeat one more time with fresh permeabilization buffer.

20m

Blocking and primary antibody incubation

17h 30m

13 Tap off the buffer. Add blocking buffer to the tissues 🌡️ Room temperature .

1h

13.1 Prepare the primary antibody dilutes in 1 x PBST.



- 14 Tap off the blocking buffer. Add antibody dilutes to the tissues 4 °C Overnight 16h
- 15 Tap off buffers. Immerse the slides in a staining dish with 1 x PBST Room temperature 00:10:00 on a horizontal shaker. Repeat the rinse for two more times. 30m
- 15.1 Prepare the secondary antibody dilutes in 1 x PBST.

Secondary antibody incubation and mounting

1h 15m

- 16 Tap off PBST. Place the slides in the humidity chamber. Add antibody dilutes to the tissues Room temperature 01:00:00 . 1h
- 17 Tap off buffers. Immerse the slides in a staining dish with 1 x PBST Room temperature 00:05:00 on a horizontal shaker. Repeat the rinse for two more times. 15m
- 17.1 If nuclei staining is need, incubate tissues with DAPI for 2 minutes at room temperature between the second and third rinses.
- 18 Tap off buffer. Add prolong Gold to the tissues and cover with cover glass.