

Jul 22, 2025

Immuno-isolation of endogenously tagged BLTP3A^{V5}

DOI

<https://dx.doi.org/10.17504/protocols.io.5jyl8qyp9l2w/v1>

Michael Hanna¹, Hely Rodriguez Cruz¹, Pietro De Camilli¹

¹Yale University



Michael G Hanna

Yale University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5jyl8qyp9l2w/v1>

Protocol Citation: Michael Hanna, Hely Rodriguez Cruz, Pietro De Camilli 2025. Immuno-isolation of endogenously tagged BLTP3A^{V5}. protocols.io <https://dx.doi.org/10.17504/protocols.io.5jyl8qyp9l2w/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 05, 2025



Last Modified: July 22, 2025

Protocol Integer ID: 219604

Keywords: ASAPCRN, isolate bltp3a, immuno, cell

Abstract

This protocol explains how to immuno-isolate BLTP3A (BLTP3A^{V5}) from endogenously edited A549 cells

Guidelines

This protocol is designed to use the following cells lines: A549 wild-type (RRID:CVCL_0023); A549 BLTP3A^{V5}, B6 clone (RRID:CVCL_E9W3); A549 BLTP3A^{V5}, C5 clone (RRID:CVCL_E9W4)

Materials

15 cm plates (Corning 353025)
10 cm plates (Corning 430293)
dPBS (Sigma-Aldrich D8537)
anti-V5 resin (ChromoTek v5tma)
Trypsin (Gibco 25300054)
Eppendorf tubes (Fischer 05-408-138)
Cell scraper (Fisher 353086)
Magnetic rack (Fischer 12321D)

Lysis buffer:

25 mM Tris pH 7.4 (Millipore Sigma T1503)
150 mM NaCl (Millipore Sigma S5886)
1% Triton (American Bioanalytical ab02025)
Protease inhibitors (Sigma-Aldrich 05056489001)

Complete media for A549 cells:

DMEM (Gibco 11965118)
10% FBS (Gibco/LifeTech 16140-071)
Pen-Strep (Gibco/LifeTech 15140-122)
Plasmocin (Invivogen ant-mpp)



Expanding cells

- 1 From one or more 10 cm plates at ~80-90% confluency, wash once with PBS before trypsinizing cells at 37C for 3-5 mins
- 2 Collect cells lifted off of the dish using complete media into a 15 mL Falcon tube 4m
- 3 Count the number of cells per mL of media using a hemocytometer or your choice of counting method 5m
- 4 Plate 1.5E6 cells onto to two 15 cm plates per condition per replicate: A549 wild-type and BLTP3A^V5 KI, in triplicate; for a total of 12 × 15 cm plates (six wild-type and six knock-in) 10m
- 5 Let the cells grow until ~90% confluent (typically 3-4 days) 4d

Lysate and anti-V5 preparation

3h 17m

- 6 Remove media from 15 cm plate and wash with at least 10 mL dPBS and then remove 5m
- 7 At 1 mL dPBS to each 15 cm plate (it works best to work with only two plates at a time to avoid cells drying out) 1m
- 8 Using a cell scraper, scrape up as many cells as you can 5m
- 9 Pipette up as many cells as you can using the 1 mL of dPBS added to the dish and collect 2 × 15 cm plates in a single 2 mL Eppendorf tube on ice 5m
- 10 Spin down cells at 500g for 5 min in a refrigerated (4C) table top centrifuge 5m
- 11 Remove supernatant and resuspend cells in ice cold 1 mL lysis buffer 1m
- 12 Let cells rotate at 4C for 30 mins 30m



12.1 While cells rotate, prepare anti-V5 resin: add 75 uL resin to 1 mL ice cold lysis buffer in a new 2 mL eppendorf tube, invert 3 times, place on magnetic rack; leave at 4C until ready

5m

13 Spin down lysate at 17000g for 10 mins in a refrigerated (4C) table top centrifuge

10m

14 Collect supernatant in a new, ice cold 2 mL eppendorf tube

5m

15 Remove lysis buffer from resin and add lysate to resin

5m

16 Let lysate rotate at 4C with resin for 2 hours

2h

Resin washing and storage

22m

17 Place tubes on magnetic rack at 4C for 1 min to capture resin

1m

18 Remove supernatant and replace with 1 mL ice cold lysis buffer

5m

19 Rotate resin for 5 min at 4C

5m

20 Repeat this wash one more time (total of 2 washes)

10m

21 After final wash, place tubes on magnetic rack at 4C for 1 min to capture resin and then remove supernatant

1m

22 Spin down tube at 500g for 2 mins to pellet all buffer

2m

23 Put tubes back on magnetic rack at 4C for 1 min to capture resin

1m



- 24 Using a 200 uL pipette, remove the last little bit of buffer from the bottom of the tubes while still on the magnetic rack
- 25 Store resin in absence of buffer at -80C

5m