



Nov 13, 2024

Immunoblotting protocol for dopamine and cortical iPSC-derived neurons with different Parkinsonian mutations

DOI

dx.doi.org/10.17504/protocols.io.yxmvm9xpol3p/v1

Nicolas Dzamko¹, Jessica Chedid¹

¹University of Sydney



courtney.wright Wright

University of Sydney

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.yxmvm9xpol3p/v1>

Protocol Citation: Nicolas Dzamko, Jessica Chedid 2024. Immunoblotting protocol for dopamine and cortical iPSC-derived neurons with different Parkinsonian mutations. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.yxmvm9xpol3p/v1>

Manuscript citation:

<https://doi.org/10.1101/2024.10.07.617117>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: November 13, 2024

Last Modified: November 13, 2024

Protocol Integer ID: 111984

Keywords: cortical cells, immunoblot, western blot, PD, alpha synuclein, P62, GBA, Rab-10, B-actin, neurons with different parkinsonian mutation, different parkinsonian mutation, different parkinson, parkinson, cortical ipsc, derived neuron, glucocerebrosidase protein, immunoblotting protocol, method for immunoblotting, differentiated vmda, mutations for quantification, mutation, protein

Funders Acknowledgements:

Michael J Fox Foundation

Grant ID: ASAP-000497

Abstract

This protocol details the method for immunoblotting of differentiated VMDA and cortical iPSC's from different Parkinson's mutations for quantification of Parkin, α -syn, Rab-10 and glucocerebrosidase proteins.



Materials

Solutions

- *1 x MES buffer*

50 mL 20x MES Buffer

950 mL dH₂O

*make 1 L at a time

- *10 x Transfer buffer (1L):*

105 g Glycine

22.5 g Tris

1L dH₂O

- *1 x Transfer Buffer (15% methanol)*

100 mL 10x Transfer buffer

750 mL dH₂O

150 mL Methanol

Note: Should be made up fresh (place in cold room or on ice to cool prior to transfer)

- *10x TBS buffer (1L):*

87 g Sodium Chloride

12 g Tris

* pH to 7.4

- *1x TBS + 0.1% Tween 20 (TBS-T)*

100 mL 10x TBS buffer

1 mL Tween-20

900 mL dH₂O

*use magnetic stirrer to dissolve Tween into solution

- *5% Skim Milk in TBST*

5 g skim milk powder

100 mL TBS-T

- *Stripping buffer (1L)*

15 g glycine

1 g SDS

10 ml Tween20

- Adjust pH to 2.2

- Bring volume up to 1 L with ultrapure water.

Antibodies



	Antibody	Company	Cat. No.	Dil. Factor	Molecular weight
	alpha synuclein	BD Biosciences	610787	1:1000	14kDa
	P62/SQSTM1	Abcam	ab56416	1:1000	62 kDa
	GBA	Abnova	H00002629-M01	1:1000	62 kDa
	Rab-10	Cell Signalling Technologies	8127S	1:1000	23 kDa
	B-actin	Abcam	ab6276	1:50000	37 kDa

Primary antibodies


	Antibody	Company	Cat. No.	Dil. Factor
	Anti-Mouse IgG (H+L) HRP Conjugate	Bio-Rad	1706516	1:5000
	Anti-Rabbit IgG (H+L) HRP Conjugate	Bio-Rad	1706515	1:5000

Secondary antibodies

Troubleshooting



Safety warnings

 For hazard information and safety warnings, please refer to the SDS (Safety Data Sheet) for each of the raw materials used.

Ethics statement

Experiments using human iPSCs were approved by the Human Research Ethics Committee at the University of Sydney (2017/094). Application of this protocol may need prior approval by the users' institutional review board (IRB) or equivalent ethics committee(s).

Before start

Ventral midbrain progenitors and Cortical neural progenitors were generated from iPSC according to the methods of Gantner and colleagues (1, 2).



BCA Assay and Sample preparation

- 1 Thaw samples and centrifuge at 12000 rpm, 4°C , 4°C, 00:20:00
- 2 Transfer supernatant to new microcentrifuge tubes and discard pellet.
- 3 Perform BCA assay as per manufacturer's protocol to determine protein concentration of lysates
- 4 Prepare 30 µL sample for loading by adding 20 ug of protein lysate to required volume of dH2O and 7.5 µL of 4X NuPage LDS sample Buffer
- 5 Store samples at -20 °C until ready to run western blot

Running Gels



1h 10m

- 6 Thaw samples from freezer on ice and place on heating block for 00:10:00 at 70 °C 10m
- 7 While waiting for samples to heat up remove the white sticker and gel comb carefully from gel(s) and rinse with dH2O.
- 8 Place 12-well gel(s) into casket and clip in place. Fill up the centre part with 1 X MES buffer and check for any leaks. Once in place, add 1 X MES to the outer space until half way full.
- 9 After samples have heated, quickly vortex and spin down.
- 10 Load 4 µL of the ladder and 25 µL of the samples prepared at 1 mg/mL (or as needed)
- 11 Run gel at 180 V (constant voltage) for 01:00:00 at Room temperature 1h








Transfer and blocking

2h 30m

- 12 Obtain large container and fill with 1 X transfer buffer
- 13 For each gel obtain a sandwich apparatus, 2 X sponges, 2 X filter paper, 1 X nitrocellulose membrane (0.45 um)
- 14 Open sandwich apparatus and place components in the following order: sponge, filter paper, gel, nitrocellulose membrane, filter paper, sponge
- 15 Perform transfer at 90 V (constant voltage) for  01:30:00 in cold room 1h 30m
- 16 Place membrane in container and cover with 5% skim milk in TBS-T solution, place on rocking platform for  01:00:00 1h

Primary and secondary antibody incubations

1h 20m

- 17 Prepare primary antibody in 5% skim milk in TBS-T solution (dilution 1:1000)
- 18 Cut blots using scissors into smaller strips as desired
- 19 Incubation at  4 °C overnight with gentle rocking
- 20 Wash blots with TBS-T every 3 times for  00:10:00 10m
- 21 Prepare secondary antibody (anti-mouse or anti-rabbit HRP conjugated antibody) in 2.5% skim milk in TBS-T solution
- 22 Incubation at  Room temperature for  01:00:00 with gentle rocking 1h
- 23 Wash blots with TBS-T every 3 times for  00:10:00 10m






Blot imaging

- 24 Visualise bands using BioRad ChemiDoc MP Imaging System chemiluminescence mode
- 25 Band analysis and quantitation is performed in ImageLab, all proteins are normalised to B-actin used as a loading control

Stripping and

1h 20m

- 26 If stripping membrane is required to reprobe with another antibody, add Abcam mild stripping buffer incubate for  00:10:00 at Room temperature with rocking. 10m
- 27 Wash membrane 3 times with TBS-T for  00:10:00 each, then block in 5% skim milk in TBS-T solution for  01:00:00 1h 10m
- 28 Repeat step 17 onwards with new primary antibody as required.

Protocol references

1. Gantner, C.W., et al., An optimized protocol for the generation of midbrain dopamine neurons under defined conditions. *Star Protocols*, 2020. 1(2): p.100065.
2. Gantner, C.W., et al., FGF-MAPK signaling regulates human deep-layer corticogenesis. *Stem cell reports*, 2021. 16(5): p. 1262-1275