

Apr 17, 2025

Immunocytochemistry

 In 2 collections

DOI

<https://dx.doi.org/10.17504/protocols.io.dm6gpdoz1gzp/v1>

Lucas Blasco-Agell^{1,2}, Veronica Testa^{1,2}, Valentina Baruffi^{1,2}, Antonella Consiglio^{1,2}, Miquel Vila³

¹Department of Pathology and Experimental Therapeutics, Bellvitge University Hospital- IDIBELL, 08908 Hospitalet de Llobregat, Spain;

²Institute of Biomedicine of the University of Barcelona (IBUB), Carrer Baldori Reixac 15-21, Barcelona 08028, Spain;

³VHIR-CIBERNED-ASAP

Vilalab Public



Miquel Vila

VHIR-CIBERNED-ASAP

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.dm6gpdoz1gzp/v1>



Protocol Citation: Lucas Blasco-Agell, Veronica Testa, Valentina Baruffi, Antonella Consiglio, Miquel Vila 2025.

Immunocytochemistry. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.dm6gpdoz1gzp/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 17, 2025

Last Modified: April 17, 2025

Protocol Integer ID: 130937

Keywords: detection of key microglia marker, ipscderived human microglia cell, key microglia marker, human microglia cell, immunocytochemistry immunocytochemistry

Funders Acknowledgements:

ASAP

Grant ID: ASAP-020505

Abstract

Immunocytochemistry (ICC) detection of key microglia markers in iPSCderived human microglia cells (hMG)

- 1 Cells plated in 24-well plates are fixed before any ICC staining using 4% of Paraformaldehyde (PFA; EMS) at RT for 15 minutes. Excess of PFA is removed from the culture wells with three washes of 15 minutes with DPBS.
- 2 For whole vmO IN TOTO ICC, samples are fixed in PFA at 4°C for 2 hours and rinsed 3 times with TBS for 20 minutes.
- 3 For ICC of cryosectioned vmO, samples are fixed in 4% PFA at 4°C Overnight (ON). Fixed vmO are transferred to 30% Sucrose (Sigma-Aldrich) solution at 4°C ON.
- 4 vmO are placed on plastic molds with O.C.T.TM compound (Tissue-Tek[®]) at -80°C until processing. Cryopreserved vmO were sectioned (20 μ m thick) on a cryostat (Leica) and mounted in Menzel-Glaser Superfrost[®] PLUS coverslips (Thermo Fisher Scientific), dried at RT ON and stored at -20°C for further ICC staining.
- 5 To stain hMG, we used a low Triton X-100 (Triton; Sigma-Aldrich) protocol. Coverslips are placed on Menzel-Glaser Superfrost[®] Microscope slides (Thermo Fisher Scientific) inside a humid chamber.
- 6 Samples are then rinsed in 1X Tris-Buffered Saline (TBS) before being blocked at RT for 2 hours with Blocking Solution (BS), which consisted of TBS with 0.01% of Triton and 3% of Donkey Serum (DS). After that, samples are incubated with primary antibodies at 4°C for 48 hours.
- 7 Samples are rinsed and blocked again with BS at RT for 1 hour and incubated with secondary antibodies for 2 hours at RT protected from light. All secondary antibodies employed are from Alexa Fluor Series (Jackson ImmunoResearch Europe) at a concentration of 1:250.
- 8 After that, samples are rinsed TBS and nuclei are counterstained with 0.5 μ g/mL of 4,6-Diamidino-2-phenylindole (DAPI; Abcam) at RT for 10 minutes.
- 9 Samples are mounted within a glass coverslip (Corning[®]) with Polyvinyl Alcohol-1,4-Diazabicyclo-Octane (PVA-DABCO; Sigma-Aldrich).
- 10 Mounted slides are dried at RT for 30 minutes protected from light and stored at 4°C before analysis. To stain vmDAn, 2D co cultures or vmO cryosections, we used a high Triton protocol except that rinses are done with TBS with 0.1% of Triton and BS consists of TBS with 0.3% of Triton and 3% of DS.
- 11 For the whole organoid structures, *IN TOTO* ICC staining is performed. vmO are rinsed with TBS and incubated with an antigen retrieval buffer composed by Sodium Citrate Tribasic (Sigma-Aldrich) with a pH of 9 at 60°C for 1 hour.



- 12 vmO are rinsed in TBS at RT in agitation (Labnet International) and permeabilized with TBS plus 0.5% Triton and 6% DS at 4°C ON in agitation. Samples are then incubated with primary antibodies, diluted in BS (TBS with 0.1% Triton and 6% DS) at RT at 4°C over-the-weekend in agitation.
- 13 After 3 days, samples are rinsed with BS at RT and incubated with secondary antibodies (diluted in BS) at 4°C ON in agitation.
- 14 The day after, vmO are rinsed TBS and nuclei are counterstained with 2.5 µg/mL DAPI (diluted in TBS). Samples are mounted with PVA-DABCO in a glass coverslip, dried and stored at 4°C for further imaging. Images are acquired using Zeiss AXIOIMAGER Z1 with an ApoTome (Zeiss Microscopy) microscope using a Cascade CCD camera (Photometrics). For confocal images, Zeiss AXIOIMAGER Z1 with an SPE confocal system or a Zeiss LSM 880 are employed (Zeiss Microscopy). Image acquisition is realized with ZEN pro software (Zeiss Microscopy) and subsequently analyses are performed using ImageJ (NIH), Photoshop® (Adobe) and IMARIS (Bitplane copyright).