



Aug 16, 2023

Immunocytochemistry (ICC) for Progenitor markers FOXA2 and OTX2

DOI

<https://dx.doi.org/10.17504/protocols.io.q26g7ppp3gwz/v1>

Ning Xia^{1,2}

¹Massachusetts General Hospital; ²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network



Daniel Choo

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.q26g7ppp3gwz/v1>

Protocol Citation: Ning Xia 2023. Immunocytochemistry (ICC) for Progenitor markers FOXA2 and OTX2. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.q26g7ppp3gwz/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 16, 2023

Last Modified: May 31, 2024

Protocol Integer ID: 86555

Keywords: ASAPCRN, progenitor markers foxa2, protocol for immunocytochemistry, immunocytochemistry, icc, otx2, protocol

Funders Acknowledgements:

Aligning Science Across Parkinson's (ASAP)

Grant ID: ASAP-000312

Abstract

This is the protocol for Immunocytochemistry (ICC) for Progenitor markers FOXA2 and OTX2.

Attachments



Immunofluorescence

P...

17KB

Materials

	A	B	C
	Blocking Buffer (BP)	1x PBS 5% BSA 0.3% Triton X-100	To make 10mL of Blocking Buffer: 0.5g of BSA 30µl of Triton X-100 Fill to 10mL with 1x PBS
	1° antibody (-20°C fridge)	FOXA2 Rb OTX2 Mus IgG	
	2° antibodies (4°C fridge)	A11008 Alexa Flour 488 goat anti-rabbit IgG (HHH+L) A10036 Alexa Flour546 donkey anti-mouse IgG (H+L)	Thermo Fisher Cat # A11008 Thermo Fisher Cat # A10036
	DAPI (-20°C fridge)	Stock 1mg/mL	Final (1µg/mL)
	- Each well of both 8µm Chamber slide and 48-well plate should use 200µL of solution		



- 1 Grow the cells on either 8µm chamber slide or 48-well plate. (1×10⁵cells/well)
 - Coating and plating densities can be found in the differentiation kit
- 2 Aspirate medium and rinse with 1x PBS
- 3 Fixation: Cover cells with warm 4% PFA and fix for 20 min at RT.
 - Cover the vessel and place in the hood
- 4 Rinse 3X in 1x PBS for 5 min each at the lab bench
- 5 Blocking and permeabilize for 30 min at RT with Blocking Buffer
- 6 1° antibody incubation dilute in 1/3 blocking buffer + 1x PBS and incubate overnight in the cold room
 - 20° (box labelled antibodies for iPSCs)
 - FOXA2 Rb à (1:250)
 - OTX2 Mus IgG à (1:250)
 - o Add 4µl of antibody per 1mL of prepared antibody solution
- 7 Rinse three times in 1x PBS for 5 min each
- 8 2° antibody incubation for 1hr at RT in the dark
 - A11008 Alexa Flour 488 goat anti-rabbit IgG (HHH+L) 1:400 in PBS
 - A10036 Alexa Flour546 donkey anti-mouse IgG (H+L) 1:400 in PBS
 - o Add 2.5µL of antibody per 1mL of prepared antibody solution
 - *Combine these two antibodies in one tube. Ex. If you were making 1.6mL (for 8 wells) you would add 4µL of each antibody
- 9 Rinse 3X in 1X PBS for 5 min each
- 10 Incubate with DAPI at (1µg/mL) for 10min at RT (Stock is 1mg/mL which is 1000x)
 - a. Add 1µl of Stock DAPI to 1mL of prepared solution.



11 11. Rinse 2X in 1X PBS and image.