

May 27, 2025

Immunocytochemistry of autophagosome-associated proteins

DOI

<https://dx.doi.org/10.17504/protocols.io.rm7vzqpb8vx1/v1>

Sierra Palumbos¹, Erika Holzbaur¹

¹University of Pennsylvania



Sierra Palumbos

University of Pennsylvania

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.rm7vzqpb8vx1/v1>

Protocol Citation: Sierra Palumbos, Erika Holzbaur 2025. Immunocytochemistry of autophagosome-associated proteins. protocols.io <https://dx.doi.org/10.17504/protocols.io.rm7vzqpb8vx1/v1>

Manuscript citation:

<https://doi.org/10.1101/2024.11.07.621551>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: May 27, 2025

Last Modified: May 27, 2025

Protocol Integer ID: 219002

Keywords: immunocytochemistry of autophagosome, autophagosome associated protein, autophagosome, proteins protocol for immunocytochemistry, associated proteins protocol, associated protein, primary cortical neuron, immunocytochemistry, protein, antibody, neuron, other cultured cell type, cortical neuron

Funders Acknowledgements:

National Institute of Neurological Disorders and Stroke

Grant ID: R01-NS060698

National Institute of Neurological Disorders and Stroke

Grant ID: F32-NS129586

Michael J. Fox Foundation for Parkinson's Research (MJFF)

Grant ID: MJFF-021130

Michael J. Fox Foundation for Parkinson's Research (MJFF)

Grant ID: MJFF-15100

Michael J. Fox Foundation for Parkinson's Research (MJFF)

Grant ID: MJFF-019411

Abstract

Protocol for immunocytochemistry used on primary cortical neurons in Palumbos et al., 2025. Protocol can be adapted to other cultured cell types and antibodies. Described protocol works for autophagosome associated proteins



Fix Cells

- 1 Heat 4% paraformaldehyde (PFA) to 37 °C
- 2 Add 4% sucrose by volume
- 3 Wash cultured cells 2X with warmed 1X PBS
- 4 Fix with 4% PFA + sucrose for 00:08:00 at 37 °C
- 5 Wash with 1X PBS for 5 minutes
- 6 Repeat wash 2X [go to step #5](#)
- 7 Stopping point - store at 4°C in 1X PBS or continue



Permeabilize Cells

- 8 Add permeabilizing agent and incubate cells.

Note

Permeabilization method should be determined based on protein(s) to be detected.
Follow Step 8.2 (MeOH permeabilization) for autophagosome or lysosome associated proteins
Follow Step 8.3 (Triton-X permeabilization) for other visualization
*You may need to try both methods to find best fit for a given antibody

- 8.1 MeOH permeabilization - Add 1 mL cold (-20°C) pure MeOH to well and place in -20°C for 8 minutes.



- 8.2 Triton X permeabilization - Add 0.2% Triton X-100 in PBS to well for 15 minutes at room temperature (not ideal for autophagosomes or lysosomes).
- 9 Wash 3X for 5 minutes with 1X PBS.

Block Cells

- 10 Prepare block (5% goat serum, 1% BSA, + sodium azide) in advance. Filter sterilize and keep at 4°C.
- 11 Add block to plates and incubate for  01:00:00 at  Room temperature .

Primary Antibody Incubation

- 12 Prepare primary antibody dilution in block (1:250 good starting dilution if unknown).
- 13 Incubate overnight with primary antibody at 4°C.

Secondary Antibody Incubation

- 14 Wash 3X for 5 minutes with 1X PBS.
- 15 Prepare secondary antibody dilution in block (1:1000 good starting dilution).
- 16 Incubate at  Room temperature for  02:00:00 .
- 17 Wash 3X for 5 minutes with 1X PBS.
- 18 Add Hoechst 1:1000 for 10 minutes.



19 Wash 3X for 5 minutes with 1X PBS.

20 Mount coverslip for imaging.