



Immunocytochemistry Staining Protocol V.3

DOI

dx.doi.org/10.17504/protocols.io.tdmei46



BioLegend, Inc., Kelsey Knight¹

¹BioLegend

BioLegend

Tech. support email: tech@biolegend.com



Sam Li

BioLegend

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account



DOI: <https://dx.doi.org/10.17504/protocols.io.tdmei46>

External link:

http://www.biolegend.com/media_assets/support_protocol/Immunofluorescence_Microscopy_Protocol_050514.pdf

Protocol Citation: BioLegend, Inc., Kelsey Knight . Immunocytochemistry Staining Protocol. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.tdmei46>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Created: September 10, 2018

Last Modified: August 29, 2019

Protocol Integer ID: 15501

Keywords: immunocytochemistry staining protocol, staining protocol, immunocytochemistry, staining

Troubleshooting

Before start


Reagent List:

<ul style="list-style-type: none">▪ Chamber slides, cover slips, or 12-well plates▪ Phosphate-buffered saline (PBS)	<ul style="list-style-type: none">▪ Fixation solution: 1% Paraformaldehyde, in PBS▪ Permeabilization solution: 0.5% Triton X-100 in PBS▪ Blocking buffer: 5% FBS in PBS
--	---



Sterilization (for 12-well plates with coverslips)


- 1 Transfer a single cover slip into a 12-well plate. Then add 1mL of 70% Ethanol into a well for 20 minutes at room temperature.

 00:20:00

- 2 Wash quickly three times with PBS.

Poly-Lysine Coating for 12-Well Plates (optional; for loosely attached cells)

- 3 Add 1 mL of 0.1 mg/mL Poly-D-lysine solution into a well for 15 minutes at room temperature.

 00:15:00


- 4 Wash quickly three times with PBS and let dry before plating cells.

Sample Preparation

- 5 Grow cultured cells on cover slips or in wells overnight at 37°C. At the time of fixation, cells should be ~70-80% confluent in single layer.


- 6 Rinse cells briefly in PBS.

- 7 Fix cells by incubation with freshly made 1% Paraformaldehyde in PBS for 10 minutes at room temperature.

 00:10:00


- 8 Rinse three times quickly in PBS.

- 9 For intracellular staining, add permeabilization solution and incubate at room temperature for 10 minutes. Then wash quickly three times in PBS.

 00:10:00

Sample Blocking

- 10 Block samples in 1 mL of blocking buffer at room temperature for 30 minutes.


 00:30:00




Sample Staining


- 11 Dilute the primary antibody to the recommended concentration/dilution in blocking buffer.
- 12 For 8-well chamber slides, add 200 μ L per well. For 12-well plates, add 500 μ L per well. Incubate two to three hours at room temperature or overnight at 4°C. If using conjugated antibodies, perform this step in the dark.
- 13 For surface staining, rinse 3 times quickly in PBS. For intracellular staining, quickly wash once followed by incubation with wash buffer for 5-10 minutes. Then quickly wash additional two times.


Note: If using primary antibodies directly conjugated to fluorochromes, then skip to step 17.

 00:10:00

- 14 Prepare fluorochrome-conjugated secondary antibody in blocking buffer according to the manufacturer's specification data sheet, and add 200 μ L per well to the 8-well chamber slides. For 12-well plates, add 500 μ L per well.
- 15 Incubate the samples for one hour, at room temperature, in the dark.
- 16 For surface staining, rinse three times quickly in PBS. For intracellular staining, quickly wash once followed by incubation with wash buffer for 5-10 minutes, then quickly wash additional two times.
- 17 Optional: To stain F-actin, prepare a working solution of Flash Phalloidin™ by diluting it 1:20-1:100 in PBS. Add 200 μ L per well for an 8-well plate or 500 μ L per well for a 12-well plate. Stain for 20 minutes at room temperature in the dark.

 01:00:00

 00:10:00

 00:20:00

- 18 Apply anti-fade mounting medium to the cover slip.
- 19 Seal slides with nail polish.