

Dec 17, 2024

Immunocytochemistry with Primary Cultured Hippocampal Neurons

 In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.eq2ly6j2egx9/v1>

Mukesh Kumar¹, Timothy A. Ryan^{1,2}

¹Department of Biochemistry, Weill Cornell Medicine, New York, NY 10065;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, Maryland 20815, USA



Mukesh Kumar

Weill Cornell Medical College

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.eq2ly6j2egx9/v1>

Protocol Citation: Mukesh Kumar, Timothy A. Ryan 2024. Immunocytochemistry with Primary Cultured Hippocampal Neurons. protocols.io <https://dx.doi.org/10.17504/protocols.io.eq2ly6j2egx9/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: November 27, 2024

Last Modified: December 17, 2024

Protocol Integer ID: 114395

Keywords: ASAPCRN, immunocytochemistry with primary cultured hippocampal neuron, primary cultured hippocampal neuron, hippocampal neurons this protocol, hippocampal neuron, immunocytochemistry, performing immunocytochemistry, neuron

Funders Acknowledgements:

NIH

Grant ID: NS036942

NIH

Grant ID: NS11739

ASAP

Grant ID: ASAP-024404

Abstract

This protocol describes a detailed method for performing immunocytochemistry (ICC) on hippocampal neurons cultured from post-natal day 0-1 pups.

Guidelines

The protocol needs prior approval by the users' Institutional Animal Care and Use Committee (IACUC) or equivalent ethics committee.

Materials

- Hippocampal neurons: DIV-14
-  Poly-D-lysine hydrobromide (PDL) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P7886**
- 4% paraformaldehyde (ThermoScientific, J19943.K2) + 10% sucrose (SigmaAldrich, S0389)
 -  Sucrose **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S0389**
-  Triton X-100 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T8787**
-  Goat serum **Merck MilliporeSigma (Sigma-Aldrich) Catalog #G9023**
- Primary antibodies:
 -  DDHD2 Polyclonal antibody **Proteintech Catalog #25203-1-AP** ,
 -  Synapsin1/2 antibody **Synaptic Systems Catalog #106 004** ,
 -  Monoclonal Anti-TOMM20 antibody produced in mouse **Merck MilliporeSigma (Sigma-Aldrich) Catalog #WH0009804M1**
 - ,  Neuron-specific beta -III Tubulin Antibody **R&D Systems Catalog #MAB1195** ,
- NeuN (Abcam, ab104224, ab177487)
 -  Anti-NeuN antibody [1B7] - Neuronal Marker **Abcam Catalog #ab104224** ,
 -  Anti-NeuN antibody [EPR12763] - Neuronal Marker **Abcam Catalog #ab177487**
- GFAP (Abcam, ab7260)  Anti-GFAP antibody **Abcam Catalog #ab7260**
 -  GFAP Antibody **Novus Biologicals Catalog #NB300-141**),
 -  GM130 (E9Z6S) Rabbit mAb #70767 **Cell Signaling Technology Catalog #70767T**
- Alexa Fluor-conjugated secondary antibodies (All from Invitrogen)
-  BODIPY™ 493/503 **Life Technologies Catalog #D3922**
- MDH (Abgent, SM1000a)
-  ProLong Diamond Antifade Mountant **Thermo Fisher Scientific Catalog #P36970**
-  PBS, pH 7.4 **Thermo Fisher Catalog #10010031**

Equipment

- Dissection tools
- Humidified chamber (Custom made)
- Frosted microscope slides (Avantor, 16004-422)
- Confocal Microscope (Zeiss 880)



Neuronal Culture Preparation

- 1 Isolate hippocampal neurons from post-natal day 0-1 pups using standard dissection techniques.
- 2 Plate the neurons sparsely onto glass coverslips coated with poly-D-lysine.
- 3 Maintain neurons in culture media under standard conditions ( 37 °C , 5% CO₂) until they reach 14 days in vitro (DIV-14). 

Fixation

25m

- 4 Prepare a fixative solution containing 4% PFA and 10% sucrose in PBS.
- 5 Fix the neurons by gently adding the fixative solution and incubating at  Room temperature for  00:10:00 .  
- 6 Rinse the fixed neurons 3 times with PBS (5 minutes each wash).
- 6.1 Rinse the fixed neurons with PBS for  00:05:00 each wash (1/3). 
- 6.2 Rinse the fixed neurons with PBS for  00:05:00 each wash (2/3). 
- 6.3 Rinse the fixed neurons with PBS for  00:05:00 each wash (3/3). 

Permeabilization

20m



7 Permeabilize the fixed neurons by incubating in 0.25% Triton X-100 in PBS for

00:05:00 at Room temperature .

5m



8 Wash the permeabilized neurons 3 times with PBS (5 minutes each wash).



8.1 Wash the permeabilized neurons with PBS for 00:05:00 each wash (1/3).

5m



8.2 Wash the permeabilized neurons with PBS for 00:05:00 each wash (2/3).

5m



8.3 Wash the permeabilized neurons with PBS for 00:05:00 each wash (3/3).

5m



Blocking

1h 10m

9 Prepare a blocking solution containing 5% goat serum in PBS.

10 Incubate the neurons in the blocking solution for 01:00:00 at

Room temperature .

1h



11 Gently wash the neurons once with PBS for 00:10:00 .

10m



Primary Antibody Incubation

10m

12 Dilute the primary antibodies in a solution containing 5% goat serum and 0.05% Triton X-100 in PBS, following the manufacturer's instructions.

13 Apply the primary antibody solution to the neurons and incubate Overnight at

4 °C in a humidified chamber.

10m





Washing

15m

14 Wash the neurons 3 times with PBS (5 minutes each wash).



14.1 Wash the neurons with PBS for 00:05:00 each wash (1/3).

5m



14.2 Wash the neurons with PBS for 00:05:00 each wash (2/3).

5m



14.3 Wash the neurons with PBS for 00:05:00 each wash (3/3).

5m



Secondary Antibody Incubation

1h

15 Dilute the Alexa Fluor-conjugated secondary antibodies (Alexa 488, Alexa 546, Alexa 647) to a 1:500 dilution in 0.05% Triton X-100 in PBS.

16 Apply the secondary antibody solution to the neurons and incubate for 01:00:00 at Room temperature in a humidified chamber.

1h



Washing and Lipid Droplet Staining

30m

17 Wash the neurons twice with PBS (5 minutes each wash).



17.1 Wash the neurons with PBS for 00:05:00 each wash (1/2).

5m



17.2 Wash the neurons with PBS for 00:05:00 each wash (2/2).

5m



18 For lipid droplet staining: Prepare a staining solution containing 1 μ L BODIPY or 100 micromolar (μ M) MDH in PBS. Perform the third wash using this solution,

10m





incubating for  00:10:00 .

- 19 Gently wash the neurons with PBS for an additional  00:10:00 .

10m



Mounting

- 20 Carefully mount the glass coverslips onto glass slides using Prolong antifade mounting media containing BODIPY or MDH to prevent photobleaching.
- 21 Seal the coverslip with nail-polish to prevent drying.

Imaging

- 22 Use a confocal or fluorescence microscope to acquire images of the stained neurons.



Note

- Bovine Serum Albumin (BSA) can be used in place of Goat serum.
- Maintain gentle handling of neurons throughout the protocol to prevent damage of synapses.
- Protect fluorescently labeled samples from light to preserve signal intensity.
- Imaging conditions should be optimized to capture fluorescence signals without significant photobleaching.