

Jun 06, 2025

Version 1

Immunofluorescence and confocal imaging in brain sections and image analysis V.1

DOI

<https://dx.doi.org/10.17504/protocols.io.yxmvm6q6v3p/v1>

Chuyu Chen¹, Loukia Parisiadou¹

¹Northwestern University, Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD 20815



Chuyu Chen

Northwestern University, Aligning Science Across Parkinson's...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.yxmvm6q6v3p/v1>

Protocol Citation: Chuyu Chen, Loukia Parisiadou 2025. Immunofluorescence and confocal imaging in brain sections and image analysis. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.yxmvm6q6v3p/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: May 30, 2025

Last Modified: June 06, 2025

Protocol Integer ID: 219236

Keywords: ASAPCRN, Irrk2 expression in snc dopamine neuron, vulnerable ventral dopamine neuron, distinct dopamine neuronal subtype, dopamine neuron, snc dopamine neuron, distinct dopamine, resilient dorsal neuron, fibers in the dorsolateral striatum, mouse brain section, confocal imaging in brain section, neuronal subtype, aldh1a1 marker, neuron, image analysis single cell profiling study, important insights into cell, dorsolateral striatum, preclinical mouse model, brain section, single cell profiling study, confocal imaging, Irrk2, postmortem pd tissue, immunofluorescence, calbindin1, imaging, dopamine, mouse brain

Funders Acknowledgements:

Aligning Science Across Parkinson's [ASAP-020600] through the Michael J. Fox Foundation for Parkinson's Research (MJFF)
Grant ID: Grant ID: ASAP-020600

Abstract

Single cell profiling studies and research on both postmortem PD tissues and preclinical mouse models demonstrated that the Sox6/Aldh1a1 markers molecularly define the vulnerable ventral dopamine neurons, while Calbindin1 (Calb1) marks the resilient dorsal neurons. These distinct dopamine neuronal subtypes not only differ in their molecular profiling but also exhibit unique projection patterns and functional responses. Therefore, studying these neurons in a subtype-specific manner can yield important insights into cell-intrinsic mechanisms of vulnerability in PD. We first addressed whether LRRK2 is expressed in dopamine neurons, as this is a prerequisite for any cell-autonomous effects. We performed immunolabeling on mouse brain sections and confirmed LRRK2 expression in SNc dopamine neurons. We also employed confocal imaging and 3D reconstructions to assess Th+ fibers in the dorsolateral striatum.

Safety warnings

- ! These protocols need prior approval by the users' Institutional Animal Care and Use Committee (IACUC) or equivalent ethics committee

Immunofluorescence

- 1 Mice were perfused with 50 ml of PBS, followed by 50 ml of 4% paraformaldehyde in PBS. The brains were dehydrated using 30% sucrose in PBS for 48 hours and then sectioned coronally at 30 μ m thickness with a cryostat (Leica Biosystems, CM305). The slices were collected in PBS containing 0.1% sodium azide and stored at 4°C for immunohistochemistry.
- 2 Sections were incubated in 5% goat serum with 0.2% Triton X-100 for 2 hours.
- 3 Antigen retrieval was performed by incubating the sections in antigen unmasking solution (Vector, H-3300) at 100°C for 30 minutes
- 4 cooling to room temperature for another 30 minutes
- 5 wash with PBS for 5 mins
- 6 sections were incubated overnight at 4°C in 5% goat serum with 0.2% Triton X-100
- 7 incubated for an additional 48 to 72 hours at 4°C with primary antibodies in 5% goat serum with 0.2% Triton X-100
- 7.1 ***Lrrk2 signal in SNc sections:*** anti-Th (1:1000, SYSY), anti-Aldh1a1 (1:300, Abcam), anti-LRRK2/Dardarin clone N241A/34 (1:200, NeuroMab), anti-GFP (1:1000, invitrogen).
- 7.2 ***Striatal Th expressing fibers ::*** anti-Th (1:1000, SYSY), anti-Aldh1a1 (1:300, Abcam), and anti-Vmat2 (1:200, Immunostar).
- 8 washing with PBS for 5 mins
- 9 sections were incubated with secondary antibodies—Alexa FluorTM 488, Alexa FluorTM 568, and Alexa FluorTM 647 (1:300, invitrogen) in 5% goat serum with 0.2% Triton X-100 for 3hrs RT
- 10 all slices were washed with PBS and mounted using ProLongTM Diamond Antifade Mountant (invitrogen)

Confocal imaging

- 11 Confocal images of the fixed 30 μm thick striatal and midbrain sections were obtained using the Nikon A1R microscope
- 12 Fluorescence images were captured with a 20x objective for Lrrk2 expression, Th area, intensity (Supplementary Figure 6) and APEX2 expression, and 100x objective for Th fluorescence at a resolution of $1,024 \times 1,024$ pixels, maintaining constant laser power for each channel across genotypes.
Images for volume measurements of Th fibers were captured with Nikon CSU-W1 SoRa with a 60x objective at 0.1 mm intervals at $1,024 \times 1,024$ pixel resolution

Image analysis: The protein signals for Th, Aldh1a1, and Lrrk2 were measured using Imaris 10.1 software (Bitplane, Concord, USA).

13 *Lrrk2 signal in SNc sections*

- 13.1 The surface rendering function was used to segment Th/Aldh1a1 cells
- 13.2 Background subtraction was enabled; the diameter for the largest sphere was set to 10 μm and automatically thresholded, with a smoothing surface set to 3 μm
- 13.3 Mean intensity of Aldh1a1 protein within the Th surface above 1.5 times the average Th protein channel mean intensity was considered positive
- 13.4 The Lrrk2 intensity within Th-positive cells was measured automatically by the software

14 *Striatal Th expressing fibers analysis:*

- 14.1 Volume measurements of Th fibers
Surface rendering function was used to segment the Th fiber. Background subtraction was enabled, the diameter of the largest sphere was set at 0.805 μm and threshold set at 2000.
- 14.2 Segments were filtered with "Number of voxels $\text{Img}=1$ " set above 50. Th volume was measured automatically by software.

14.3 Intensity and area of Th expressing axons.

Rolling ball background subtraction was performed using Nikon Elements software.

14.4 For the high magnification analysis

The surface rendering function was used to segment TH, Aldh1a1, and Vmat2 proteins. Background subtraction was disabled and automated thresholds were applied with manual adjustments made as needed.

14.5 :

Segments were filtered with "Number of voxels lmg=1" set above 50.

14.6 For the Aldh1a1 positive Th axon, filtering was applied based on the "overlapped area ratio to surface surface=TH-aldh1a1" above 0.4.

14.7 Vmat2 within Th axonal terminal was filtered with "overlapped area ratio to surface surface= Th-VMAT2" above 0.4.

14.8 Area and intensity were measured automatically by software. For the low magnification analysis

14.9 For the low magnification analysis

The surface rendering function was used to segment TH. Background subtraction was disabled and the threshold was set to 1500 with manual adjustments made if needed

14.10 Segments were filtered with "Number of voxels lmg=1" set above 50. Th area and intensity were measured automatically by software.