

Aug 25, 2023

## Immunofluorescence and live-cell Imaging

DOI

[dx.doi.org/10.17504/protocols.io.bp2l6x9b1lqe/v1](https://dx.doi.org/10.17504/protocols.io.bp2l6x9b1lqe/v1)

Shenjie Wu<sup>1</sup>, Nancy C. Hernandez Villegas<sup>2</sup>, Iona Thomas-Wright<sup>3</sup>, Richard Wade-Martins<sup>3</sup>, schekman<sup>1</sup>

<sup>1</sup>Department of Molecular and Cell Biology, Howard Hughes Medical Institute, University of California, Berkeley, Berkeley, United States;

<sup>2</sup>Helen Wills Neuroscience Institute, University of California, Berkeley, Berkeley, United States;

<sup>3</sup>Oxford Parkinson's Disease Centre, Department of Physiology, Anatomy and Genetics and Kavli Institute for Nanoscience Discovery, University of Oxford, Oxford, United Kingdom



Nancy C. Hernandez Villegas

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.bp2l6x9b1lqe/v1>

**Protocol Citation:** Shenjie Wu, Nancy C. Hernandez Villegas, Iona Thomas-Wright, Richard Wade-Martins, schekman 2023. Immunofluorescence and live-cell Imaging. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bp2l6x9b1lqe/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



**Protocol status:** Working

**We use this protocol and it's working**

**Created:** August 02, 2023

**Last Modified:** August 25, 2023

**Protocol Integer ID:** 85900

**Keywords:** cell imaging this protocol, ipscs cell, immunofluorescence, cell imaging, cell imaging procedure, confocal microscope, microscope, zeiss lsm900, imaging, cell, different cell type, u2o

**Funders Acknowledgements:**

**Nancy C Hernandez Villegas**

Grant ID: NIH training program T32GM139780

**Randy Schekman**

Grant ID: Howard Hughes Medical Institute

## Abstract

This protocol contains a detail description of how to perform immunostaining on two different cell types, U2OS and iPSCs cells.

It also describes how to perform live-cell imaging procedure using a Zeiss LSM900 confocal microscope in a temperature-controlled environment.

## Troubleshooting



## Immunofluorescence of U2 cells

50m

- 1 U2OS cells were washed once with PBS and immediately fixed by 4% EM-grade paraformaldehyde for 00:10:00 at Room temperature 10m
- 2 Cells were washed three times with PBS for 00:10:00 each time. 10m
- 3 Blocked and permeabilized for 00:30:00 in permeabilization buffer (5% FBS and 0.1% saponin in PBS) 30m
- 4 Cells were then incubated with 1:100 dilution of primary antibodies at 4 °C Overnight 1h
- 5 Cells were washed three times with PBS for 00:10:00 each time. 10m
- 6 Cells were incubated with 1:500 dilution of fluorophore-conjugated secondary for 00:30:00 at Room temperature 30m
- 7 Prolong Gold with DAPI was used as mounting solution
- 8 Images were acquired with a Zeiss LSM900 confocal microscope and analyzed with Fiji/ImageJ software

## Immunofluorescence of hiPSC dopamine neurons

2h 40m

- 9 Cells were fixed with 4% paraformaldehyde in PBS and 0.1% Triton-X was used for permablization 00:10:00 10m
- 10 Blocked in 10% normal donkey serum for 01:00:00 Room temperature 1h
- 11 Cells were then incubated with 1:100 dilution of primary antibodies at 4 °C Overnight 1h



12 Cells were washed three times with PBS for ⌚ 00:10:00 each time. 10m

13 Cells were incubated with 1:500 dilution of fluorophore-conjugated secondary for ⌚ 00:30:00 at 🌡 Room temperature 30m

14 Prolong Gold with DAPI was used as mounting solution

15 Images were acquired with a Zeiss LSM900 confocal microscope and analyzed with Fiji/ImageJ software

## Live-cell imaging

15m

16 Cells were cultured in 35 mm glass bottom dishes (MatTek).

17 HaloTag fluorescent ligands were added according to the manufacturer's protocol (Promega).

18 After incubation for ⌚ 00:15:00 in the incubator ( 🌡 37 °C and 5% CO<sub>2</sub> ), the cells were quickly washed twice with PBS. The medium was replaced with Opti-MEM supplemented with 10% FBS. 15m

19 Imaging was performed using a Zeiss LSM900 confocal microscope in a temperature-controlled ( 🌡 37 °C and 5% CO<sub>2</sub> ) environment.