

Jun 03, 2024

# 🌐 Immunofluorescence labelling and imaging of cholinergic interneurons in post-mortem human brain tissues

📁 In 2 collections

DOI

[dx.doi.org/10.17504/protocols.io.eq2lywn5qvx9/v1](https://dx.doi.org/10.17504/protocols.io.eq2lywn5qvx9/v1)

Alan K. L. Liu<sup>1,2</sup>, Laura Parkkinen<sup>1,2</sup>

<sup>1</sup>Oxford Parkinson's Disease Centre, University of Oxford, Oxford OX1 3PT, UK;

<sup>2</sup>Nuffield Department of Clinical Neurosciences, University of Oxford, Oxford OX3 9DU, UK



Cláudia C. Mendes

University of Oxford



## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.eq2lywn5qvx9/v1>

**Protocol Citation:** Alan K. L. Liu, Laura Parkkinen 2024. Immunofluorescence labelling and imaging of cholinergic interneurons in post-mortem human brain tissues. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.eq2lywn5qvx9/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



**Protocol status:** Working

**We use this protocol and it's working**

**Created:** May 09, 2024

**Last Modified:** June 03, 2024

**Protocol Integer ID:** 99495

**Keywords:** imaging of cholinergic interneuron, cholinergic interneuron, striatal astrocytes in formalin, striatal astrocyte, astrocyte, immunofluorescence labelling, brain, formalin fixed paraffin embedded, imaging, immunofluorescence, protocol details immunofluorescence labelling

## Abstract

This protocol details immunofluorescence labelling and imaging of cholinergic interneurons along with striatal astrocytes in formalin-fixed paraffin-embedded (FFPE) post-mortem human brain tissues.

## Troubleshooting



## Collection and fixation of post-mortem human brain tissues

- 1 Collect 6  $\mu\text{m}$ -thick sections of formalin-fixed paraffin-embedded (FFPE) tissues containing the anterior basal ganglia at the level of nucleus accumbens.
- 2 Place tissue sections in an oven at 70°C for 30 mins.
- 3 Dewax tissue sections in xylene (3  $\times$  5 mins).
- 4 Rehydrate tissue sections through decreasing concentration of industrial denatured alcohol (IDA) (100%, 100%, 90%, 70%; 5 mins each) and subsequently in distilled water (5 mins).

## Primary Antibody

- 5 Perform heat-induced epitope retrieval by placing tissue sections in autoclave (121°C, 20 mins) in 0.01 M sodium citrate buffer (pH 6.0).
- 6 Rinse in PBS (3  $\times$  5 mins).
- 7 Incubate samples at 4°C overnight with a mixture of the following primary antibodies:
  - anti-ChAT antibodies (1:50, AB144P; Millipore, UK; RRID:AB\_2313845)
  - anti-GFAP antibodies (1:2000; Z0334; Agilent Dako, Santa Clara, United States; RRID:AB\_10013382).

### Note

Antibodies were diluted in 0.3% Triton X-100, 2% fetal bovine serum (FBS) and PBS.

## Secondary Antibody

- 8 Rinse samples in PBS (3  $\times$  5 mins).

- 9 Prepare Secondary Antibody mixture of Alexa Fluor 488-conjugated donkey anti-rabbit secondary antibody (1:200; A-11055; ThermoFisher Scientific, UK) and Alexa Fluor 594-conjugated donkey anti-goat secondary antibody (1:200; ThermoFisher Scientific, UK) diluted in 0.1% Triton-X in PBS.
- 10 Incubate sections in Secondary Antibody mixture in the dark (covered with foil) for 1 hour at room temperature.

## Mounting

- 11 Incubate sections with TrueBlack® Lipofuscin Autofluorescence Quencher (1:20 with 70% ethanol; Biotium, Fremont, CA, United States) for 30 seconds to block endogenous fluorescence signal.
- 12 Rinse in PBS (3 × 2 mins).
- 13 Mount and coverslip with Vectashield antifade mounting medium with DAPI (H1900, Vector Laboratories, Peterborough, UK).

## Confocal Imaging and Analysis

- 14 Visualise immunofluorescent-stained slides using Zeiss 880 Airyscan inverted confocal microscope (Carl Zeiss).
- 15 Acquire images using 10x objectives (Plan-Apochromat/0.45 NA) and 63x objectives (Plan-Apochromat, oil immersion; DIC M27/1.4 NA) with laser excitation at 405 nm (solid state), 488 nm (Argon), and 561 nm (solid state).
- 16 Perform image capture and processing using the Zen Black and Zen Blue (Carl Zeiss, Germany) software.