

Apr 01, 2023

Immunofluorescence microscopy of R1441C or VPS35 D620N MEF cells

DOI

dx.doi.org/10.17504/protocols.io.j8nlkw1n5l5r/v1

Chloe A Hecht^{1,2}, Herschel Dhekne^{1,2}, Wondwossen Yeshaw^{1,2}, Suzanne Pfeffer^{1,2}

¹Department of Biochemistry, Stanford University School of Medicine;

²Aligning Science Across Parkinson's Disease



Chloe A Hecht

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlkw1n5l5r/v1>

Protocol Citation: Chloe A Hecht, Herschel Dhekne, Wondwossen Yeshaw, Suzanne Pfeffer 2023. Immunofluorescence microscopy of R1441C or VPS35 D620N MEF cells. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.j8nlkw1n5l5r/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: March 21, 2023

Last Modified: April 01, 2023

Protocol Integer ID: 79201

Keywords: vps35 d620n mef cell, immunofluorescence microscopy, transfected myc, stain r1441c, cell, rilpl1, r1441c

Funders Acknowledgements:

Aligning Science Across Parkinson's Disease

Grant ID: 000463

Abstract

We present here a method to culture, fix, permeabilize, and stain R1441C or VPS35 D620N MEF cells to visualize transfected Myc-RILPL1 and endogenous TMEM55b or pRab10.

Materials

Paraformaldehyde 16% (Electron Microscopy Sciences 50-980-487)

PBS (Cold Spring Harbor protocol)

TMEM55B polyclonal antibody (Proteintech 23992-1-AP)

Recombinant Anti-Rab10 (Phospho T73) antibody (Abcam ab241060)

Myc antibody (Biolegend 626802)

Donkey anti-Rabbit 568 Alexa Fluor (Thermo Fisher A10042)

Donkey anti-Mouse 488 Alexa Fluor (Thermo Fisher A-21202)

DMEM high glucose (Cytiva SH30243.01)

Fetal bovine serum (Sigma #F0926)

Microscope cover glass (Fisher 12-545-81)

FuGENE 6 transfection reagent (Promega E2691)

Troubleshooting

Culturing, plating, and treating cells

- 1 Seed R1441C MEF cells or VPS35 D620N MEF cells at 50-60% confluency in a six well plate in 2 mL of complete DMEM (DMEM containing 10% FBS and 1% penicillin-streptomycin) 24 hours before transfection.
- 2 Transfect with Myc-RILPL1 plasmid with FuGENE 6 transfection reagent (E2691) at a 3:1 FuGENE:plasmid ratio using 2µg per well according to the manufacturer's guidelines for 24 hours.
- 3 Trypsinize cells and plate onto 12 mm glass coverslips resting in a six well plate at 60% confluency and leave cells to attach to coverslips in a 37°C incubator with 5% CO₂ for 24 hours.
- 4 If needed for MLi-2 conditions, add 200nM MLi-2 and leave in the incubator for 1.5 -2 hours.
- 5 For Nocodazole, add 20µM Nocodazole and leave in the incubator for 1.5-2 hours.

Fixing, permeabilizing, and mounting coverslips

- 6 Solutions needed:
 - 4% PFA in 1x PBS
 - 0.2% Saponin in 1x PBS
 - 2% BSA in 1x PBS (blocking solution)
- 7 Transfer each coverslip from the 6-well plate to a 24 well plate filled with ~400 µl 4% PFA per well and fix for 10 minutes at room temperature (all subsequent steps and solutions are at room temperature).
- 8 Wash 3x with 1x PBS.
- 9 Permeabilize cells with 0.2% Saponin in 1x PBS for 10 minutes.
- 10 Wash 3x in 1x PBS.
- 11 Block with 2% BSA in 1x PBS for 30 minutes.

- 12 Incubate primary antibody in blocking buffer for 2 hours.
- 12.1 Dilute antibodies (see Materials) as follows:
Rabbit anti-TMEM55B PA5-61760 (1:4000)
Rabbit anti-pRab10 (1:1000)
Mouse anti-Myc (1:1000)
- 13 Wash 3x with 1x PBS.
- 14 Incubate with 2° Antibody for 1 hour in the dark in blocking buffer.
- 14.1 Dilute antibodies (see Materials) as follows:
Donkey anti-Rabbit 568 for TMEM55b and pRab10 (1:2000) and
Donkey anti mouse 488 for Myc (1:2000)
DAPI (1:1000)
- 15 Wash 3x with 1x PBS.
- 16 Use a Kimwipe held on the edge of the coverslip to remove excess liquid; Invert the coverslip onto 3µl Mowiol on a Gold Seal glass slide (size 1×3", 1 mm thick).
- 17 Allow coverslips to dry at RT overnight.
- 18 Image