

Dec 06, 2023

Immunofluorescence Multi-label Protocol for Free-floating Fixed Tissue

DOI

<https://dx.doi.org/10.17504/protocols.io.j8nlkorn5v5r/v1>

Yaping Chu¹, Jeremy Molina¹, Scott Muller¹, Jeffrey H Kordower¹

¹Arizona State University

Team Kordower



Scott Muller

Arizona State University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlkorn5v5r/v1>

Protocol Citation: Yaping Chu, Jeremy Molina, Scott Muller, Jeffrey H Kordower 2023. Immunofluorescence Multi-label Protocol for Free-floating Fixed Tissue. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.j8nlkorn5v5r/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 05, 2023

Last Modified: May 31, 2024

Protocol Integer ID: 91873

Keywords: Immunohistochemistry, Immunofluorescence, ASAPCRN, floating fixed tissue immunofluorescence, immunofluorescence multi, immunofluorescence, fixed tissue, tissue in the kordower laboratory, tissue

Abstract

Immunofluorescence multi-label protocol for staining free-floating fixed tissue in the Kordower Laboratory.

Attachments



Immunofluor Multi-La...

19KB

Guidelines

HISTO- NOTES:

- Primate tissue staining dishes use  100 mL solution per dish
- Rodent tissue staining dishes  50 mL solution per dish
- If staining a large number of primate cases, incubate 1' & 2' Ab in individual cups to conserve volume of Ab used.
- Be conscious of tissue saturation while washing and incubating. i.e. Check that tissue is fully submerged in solution & not clumping. This will ensure proper penetration of antibodies & other reagents.
- Always include Positive & Negative Controls.
- Positive: Use relevant control tissue to confirm specific antibody detection. (i.e. pS129; control tissue should consist of nigral sections previously successfully stained for pS129).
- Negative: Ideally, use tissue that you know does not contain the targeted antigen. If not available, use a section of tissue not incubated in the 1' Ab (primary delete).
- When incubating 1' Ab overnight, leave on shaker in refrigerator.
- Can incubate in fridge on a shaker, covered in parafilm, over the weekend or up to 3 days.
- Select a secondary antibody directed against the species in which the primary antibody was raised (i.e. if a primary antibody raised in rabbit is used, an anti-rabbit secondary antibody raised in a species other than rabbit must be used).



Materials

- Dilution Media (DM) (1M 0.2 Mass Percent TBS plus 1M 0.05 % volume Triton X-100)
- 1M 0.2 Molarity (M) Tris-buffered saline (TBS)
- Sodium meta-periodate
- Normal Serum (species matching the host of the secondary antibody, e.g. horse, goat)
- Bovine Serum Albumin (BSA)
- Triton X-100
- Vectastain Elite ABC-HRP Kit (PK-6100)
- Imidazole
- Sodium Acetate
- 3,3-Diaminobenzidine Tetrahydrochloride (DAB)
- 1M 30 % (v/v) hydrogen peroxide
- 1M 0.2 Molarity (M) Phosphate-buffered saline (PBS)
- Household Bleach
- Primary antibody against the target antigen
- Secondary antibody directed against the species in which the primary antibody was raised (i.e. if a primary antibody raised in rabbit is used, an anti-rabbit secondary antibody raised in a species other than rabbit must be used).



DAY 1 (4 hrs)

1 Wash sections (6 x  00:10:00) in Dilution Media (DM) ( 0.2 Mass Percent TBS plus  0.05 % volume Triton X-100).

10m



1.1 Wash sections for  00:10:00 in DM (1/6).

10m

1.2 Wash sections for  00:10:00 in DM (2/6).

10m

1.3 Wash sections for  00:10:00 in DM (3/6).

10m

1.4 Wash sections for  00:10:00 in DM (4/6).

10m

1.5 Wash sections for  00:10:00 in DM (5/6).

10m

1.6 Wash sections for  00:10:00 in DM (6/6).

10m

2 Endogenous peroxidase inhibition ( 00:20:00).  0.1 Mass Percent Sodium meta-periodate in TBS.

20m

Note

Only necessary if using ABC, because it has HRP and could cause background.

-  100 mL  0.2 Molarity (M) Tris-buffered saline (TBS)
-  2.13 g Sodium meta-periodate

3 Wash (2 x  00:10:00) in DM.

10m

Note

Only necessary following endogenous peroxidase inhibition.





3.1 Wash for  00:10:00 in DM (1/2).

10m

3.2 Wash for  00:10:00 in DM (2/2).

10m

4 Serum blocking step ( 01:00:00 incubation):

1h

-  100 mL DM
-  3 mL Normal Serum (species matching the host of the secondary antibody, e.g. horse, goat)
-  2 g Bovine Serum Albumin (BSA)



5 Incubation in primary antibody ( 18:00:00 -  72:00:00). See antibody catalog for concentration of primary antibody.

3d 18h

-  100 mL DM
-  1 mL Normal Serum (species matching the host of the secondary antibody, e.g. horse, goat)
-  1 g BSA
-  0.5 mL Triton X-100



Note

****Optionally, refrigerate  4 °C to keep antibody stable****

DAY 2 (8 hrs)

2h 30m

6 Wash (6 x  00:10:00) in DM.

10m



6.1 Wash in DM for  00:10:00 (1/6).

10m

6.2 Wash in DM for  00:10:00 (2/6).

10m



- 6.3 Wash in DM for 00:10:00 (3/6). 10m
- 6.4 Wash in DM for 00:10:00 (4/6). 10m
- 6.5 Wash in DM for 00:10:00 (5/6). 10m
- 6.6 Wash in DM for 00:10:00 (6/6). 10m
- 7 Fluorophore-conjugated secondary antibody incubation. (01:00:00) Concentration of secondary antibody is always 1:200 in solvent. 1h
- 100 mL DM
 - 1 mL Normal Serum (species matching the host of the secondary antibody, e.g. horse, goat)
 - 1 g BSA
- 8 Wash (2 x 00:10:00) in TBS. 10m
- 8.1 Wash for 00:10:00 in TBS (1/2). 10m
- 8.2 Wash for 00:10:00 in TBS (2/2). 10m
- 9 Serum blocking step (01:00:00 incubation):
- 100 mL DM
 - [M] 5 % (v/v) Normal Serum (species **MATCHING THE HOST OF THE PRIMARY** antibody, saturates open binding sites on secondary antibody)
 - 1 g Bovine Serum Albumin (BSA)

Note

****DO NOT USE any detergent (i.e. Triton X-100, Tween-20, DM) from this step onward!
Detergent will wash away the fragment antibodies!****



10 Oversaturate with Fab antibody **against host of primary antibody** and from **same host species as secondary antibody**. (ex. If primary was mouse and secondary was goat anti-mouse, you would use a Fab-goat **anti-mouse** antibody). Working concentration: [M] 40 Mass Percent . (⌚ 01:00:00 incubation)

1h



- ⌚ 100 mL TBS
- [M] 40 Mass Percent Fab antibody (base concentration is 1.3 mg/mL for fab-goat anti mouse, so use $M_1V_1 = M_2V_2$ to find V_1/X).

11 Wash (2 x ⌚ 00:10:00) in TBS.



11.1 Wash for ⌚ 00:10:00 in TBS (1/2).

11.2 Wash for ⌚ 00:10:00 in TBS (2/2).

12 Incubation in **SECOND** primary antibody (⌚ 18:00:00 - ⌚ 72:00:00). See antibody catalog for concentration of primary antibody.



- ⌚ 100 mL DM
- ⌚ 1 mL Normal Serum (species matching the host of the secondary antibody, e.g. horse, goat)
- ⌚ 1 g BSA

Note

****Optionally, refrigerate ⌚ 4 °C to keep antibody stable** CAN ALSO ADD THIRD PRIMARY ANTIBODY FROM DIFFERENT SPECIES, IF NEEDED.**

DAY 3 (2 hrs)

13 Wash (6 x ⌚ 00:10:00) in TBS.



13.1 Wash in TBS for ⌚ 00:10:00 (1/6).



13.2 Wash in TBS for  00:10:00 (2/6).

13.3 Wash in TBS for  00:10:00 (3/6).

13.4 Wash in TBS for  00:10:00 (4/6).

13.5 Wash in TBS for  00:10:00 (5/6).

13.6 Wash in TBS for  00:10:00 (6/6).

14 Fluorophore-conjugated secondary antibody incubation against second (and third, if used) primary antibody. ( 01:00:00) Use different fluorophores for each secondary. Concentration of secondary antibody is always 1:200 in solvent.

-  100 mL DM
-  1 mL Normal Serum (species matching the host of the secondary antibody, e.g. horse, goat)
-  1 g BSA

15 Wash (3 x  00:10:00) in TBS.

15.1 Wash for  00:10:00 in TBS (1/3).

15.2 Wash for  00:10:00 in TBS (2/3).

Note

****Can add DAPI (1:15,000) during the second TBS washing step, if desired.****



- 15.3 Wash for  00:10:00 in TBS (3/3).
- 16 Store tissue in TBS in the refrigerator at  4 °C until mounted.
- 17 Control for Fragment antibody (Fab): Control tissue should be processed alongside experimental tissue through Day 2, Step 11. Skip second primary incubation all together (Step 12), and complete Day 3. Check under microscope to ensure there is no co-labeling between the two chosen fluorophores.
Use appropriate +/- controls.