

Jul 07, 2023

Immunofluorescence of macrophages and microglia

DOI

<https://dx.doi.org/10.17504/protocols.io.81wgby7ovpk/v1>

Narayana Yadavalli^{1,2,3,4,5}, Shawn M. Ferguson^{1,2,3,4,6,5}

¹Department of Cell Biology, Yale University School of Medicine, New Haven, Connecticut 06510, USA;

²Neuroscience, Yale University School of Medicine, New Haven, Connecticut 06510, USA;

³Program in Cellular Neuroscience, Neurodegeneration and Repair;

⁴Wu Tsai Institute Yale University School of Medicine, New Haven, Connecticut 06510, USA;

⁵Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815, USA;

⁶Kavli Institute for Neuroscience, Yale University School of Medicine, New Haven, Connecticut 06510, USA



Berrak Ugur

Yale University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.81wgby7ovpk/v1>

Protocol Citation: Narayana Yadavalli, Shawn M. Ferguson 2023. Immunofluorescence of macrophages and microglia. protocols.io <https://dx.doi.org/10.17504/protocols.io.81wgby7ovpk/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 25, 2023

Last Modified: May 31, 2024

Protocol Integer ID: 82457

Keywords: Immunofluorescence, antibody, ASAPCRN, immunofluorescence of macrophage, microglia, microglia this protocol, paraformaldehyde solution, staining of cultured cell, cultured cell, immunofluorescence, macrophage, cell

Funders Acknowledgements:

ASAP

Grant ID: 000580

Abstract

This protocol describes the fixation and staining of cultured cells in paraformaldehyde solution.

Attachments



[mqjybj3up.docx](#)

16KB



Materials

Solutions to prepare

Fixative solution: 4% paraformaldehyde in 0.2 M phosphate buffer pH 7.4.

Permeabilization and locking solution: 0.1% (v/v) saponin + 5% (w/v) BSA in PBS

Primary antibody solution: 0.1% (v/v) saponin + 5% (w/v) BSA in PBS solution

Secondary antibody solution: 0.1% (v/v) saponin + 5% (w/v) BSA in PBS solution

⊗ LAMP-1 (human) **Developmental Studies Hybridoma Bank Catalog #ID4B**

⊗ Anti-TFE3 antibody produced in rabbit **Merck MilliporeSigma (Sigma-Aldrich) Catalog #HPA023881**

⊗ ProLong[®] Gold Antifade Mountant with DAPI **Thermo Fisher Catalog #P36935**



Immunofluorescence of macrophages and microglia

3h 55m

- 1 Plate 50,000 iPSC derived macrophages or microglia or Neurons or BMDM's cells on 22 x 22 mm glass coverslips in 24-well dishes in respective culture media.
- 2 Incubate  Overnight at  37 °C in 5% CO₂.  
- 3 Pre-warm a solution of 4% paraformaldehyde to  Room temperature .
- 4 Aspirate media from cells and immediately proceed to step 5.
- 5 Add pre-warmed 4% paraformaldehyde to cells.
- 6 Incubate for  00:15:00 at  Room temperature . 

- 7 Remove 4% paraformaldehyde solution and dispose of it in an appropriate chemical waste container.
- 8 Rinse each well with 3X  1 mL PBS. Remove rinse and dispose of it in an appropriate chemical waste container. 
- 9 Add 0.1% (v/v) saponin + 5% (w/v) BSA in PBS and incubate for  01:00:00 at  Room temperature to permeabilize cells. 
 
- 10 Aspirate 0.1% (v/v) saponin + 5% (w/v) BSA in PBS solution.
- 11 Prepare respective antibody dilutions in 0.1% (v/v) saponin + 5% (w/v) BSA in PBS solution. mouse LAMP1 (DSHB #ID4B) 1:200 dilution, TFE3 (Sigma #HPA023881), 1:200 dilution.
- 12 Incubate at  4 °C  Overnight . 



13 Aspirate antibody solution and wash cells with PBS for 00:05:00 .

5m



14 Repeat wash.



14.1 Wash with PBS for 00:05:00 . (1/3)

5m

14.2 Wash with PBS for 00:05:00 . (2/3)

5m

14.3 Wash with PBS for 00:05:00 . (3/3)

5m

15 Add filtered PBS containing 5% BSA and secondary antibodies (1:600, Alexa fluorophores 488 and 555, Invitrogen).

16 Incubate at Room temperature for 01:00:00 .

1h



17 Aspirate antibody solution and wash cells.



17.1 Wash cells with PBS for 00:05:00 . (1/3)

5m

17.2 Wash cells with PBS for 00:05:00 . (2/3)

5m

17.3 Wash cells with PBS for 00:05:00 . (3/3)

5m

18 Rinse coverslips in milliQ water.



- 19 Mount coverslips onto slides using ProLong Gold Antifade Mountant with DAPI (ThermoFisher P36935).
- 20 Allow slides to cure  Overnight in darkness.
- 21 Image.

5m

