

Nov 22, 2024

Immunofluorescence of paraffin-embedded tissue sections

DOI

<https://dx.doi.org/10.17504/protocols.io.n92ld8mk8v5b/v1>

Carolina Lopez¹

¹Washington University



Cecilia Escudero

wustl

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n92ld8mk8v5b/v1>

Protocol Citation: Carolina Lopez 2024. Immunofluorescence of paraffin-embedded tissue sections. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.n92ld8mk8v5b/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 28, 2024

Last Modified: November 22, 2024

Protocol Integer ID: 100735

Keywords: immunofluorescence of paraffin, immunofluorescence, embedded tissue section, tissue sections protocol, embedded tissue sections protocol, tissue section, paraffin, tissue

Abstract

Protocol for immunostaining of paraffin-embedded tissue sections

Materials

	A	B	C
	Reagent ID	Vendor	Catalog #
	Antigen Unmasking Solution, Tris-Based (100X)	Vector	H3301
	Antigen Unmasking Solution, Citric Acid Based (100X)	Vector	H-3300
	Xylene	Fisher	X3P-1GAL



Protocol

2h 25m

1 De-Paraffinization and re-hydration of tissue sections

30m

For the following incubations, prepare coplin jars or any container that can completely enclosure your slide rack, one per solution. The solutions should completely cover the slides. Between incubations, just move the slide to the next container.

a) Xylene (3 Changes) ⌚ 00:05:00 each

b) 100% EtOH ⌚ 00:05:00

c) 95% EtOH ⌚ 00:05:00

d) 90% EtOH ⌚ 00:05:00

e) 70% EtOH ⌚ 00:05:00

f) dH₂O ⌚ 00:05:00

2 Antigen Retrieval (if using a steamer)

1h 20m

1. Preheat steamer in the "cooking" or "boiling" function. Careful to not put a lot of dH₂O in the pot so it will overflow to the slide container.

2. Prepare the correct amount of 1x antigen unmasking solution in dH₂O and fill a coplin jar or plastic container where the slides will be incubated. Example: 🧴 200 mL for a plastic container jar (12 slide rack). Fill up the container, submerge the slide rack and transfer the container + unmasking solution + slide rack to the preheated steamer. Incubate for ⌚ 00:20:00 at max temperature.

3. Turn of steamer and take off the coplin jar containing the slides. Let it cool down at 🌡 Room temperature for ⌚ 01:00:00 .

4. Wash the slides in 1x PBS

OBS.: Troubleshoot the best unmasking solution for your staining. Examples: The Citric Acid Based solution works good for staining Foxp3 in mouse lungs, while Tris-based solution works better for detecting both virus (SARS-CoV-2 N protein) and host cytoplasmic proteins in mouse lungs.

3 Antigen Retrieval (If using a hot plate)

1) Place a 2L beaker containing 1X antigen unmasking solution on a hot plate, cover with foil, and turn plate on high for ⌚ 00:15:00 , then turn low to maintain a mild boil at

🌡 102 °C . Check temperature with a thermometer.

2) With forceps, slowly submerge slide rack into mildly-boiling beaker. Boil the slides for ⌚ 00:20:00 .



3) Turn off the hot plate and let the slides cool down inside the beaker containing antigen unmasking solution for  01:00:00 , or until the temperature drops to around

 40 °C .

4) Wash slides 3-5 times by moving rack up and down in PBS.

4 Immunofluorescence Staining

- Once sections are de-paraffinized, re-hydrated and after antigen retrieval the staining can be performed using the **PROTOCOL FOR TISSUE IFA** (<https://www.protocols.io/edit/protocol-for-tissue-ifa-cryopreserved-sections-demmm3c46>) described previously.