

Jan 23, 2025

## Immunofluorescence on acute hippocampal slices following electrophysiology

DOI

<https://dx.doi.org/10.17504/protocols.io.4r3l2923qv1y/v1>

Quinn Pauli<sup>1</sup>, Robert Bonin<sup>1</sup>

<sup>1</sup>Leslie Dan Faculty of Pharmacy, University of Toronto, Toronto, Ontario, Canada



Quinn Pauli

University of Toronto

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

[Create free account](#)

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.4r3l2923qv1y/v1>

**Protocol Citation:** Quinn Pauli, Robert Bonin 2025. Immunofluorescence on acute hippocampal slices following electrophysiology . **protocols.io** <https://dx.doi.org/10.17504/protocols.io.4r3l2923qv1y/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** December 05, 2024

**Last Modified:** January 23, 2025

**Protocol Integer ID:** 114319

**Keywords:** Immunofluorescence , Staining , Hippocampal slices, Electrophysiology, IF, immunofluorescence on acute hippocampal slice, immunofluorescence on hippocampal slice, acute hippocampal slice, labeled protein, hippocampal slice, spatial quantification of protein, immunofluorescence, using immunofluorescence, following electrophysiology recording, protein detection, electrophysiology recording, electrophysiology, glutamate receptor, following electrophysiology, protein, cryopreservation

**Funders Acknowledgements:**

**Natural Sciences and Engineering Research Council of Canada**

Grant ID: RGPIN-2024-06215

## Abstract

This protocol describes protein detection using immunofluorescence on hippocampal slices obtained from adult mice following electrophysiology recordings. Here, we detail the fixation, cryopreservation, slicing and staining methods to produce reliable detection of fluorescently labeled proteins. It is intended to enable spatial quantification of proteins (e.g. glutamate receptors) following electrophysiology recordings.

## Materials

### Reagents:

Phosphate-buffered saline (PBS)  
Paraformaldehyde 4% (in PBS)  
Triton-X  
Horse serum (or other blocking agent)  
Sucrose  
Optimal cutting temperature compound (OCT)  
Dry ice  
Mountant (e.g. ProLong<sup>TM</sup> Glass Antifade, ThermoFischer Scientific P36980)

### Equipment:

Cryostat, chucks and blades  
Microscope slides  
Culture dishes  
Pasteur pipettes  
Cryomolds  
Parafilm  
Tissue spatula  
Tweezers  
PAP pen  
Coverslips  
Clear nail polish

## Troubleshooting

## Safety warnings

 If you choose to use this protocol, you do so at your own risk and must ensure that any local guidance is adhered to.

## Ethics statement

Experiments involving animals must be conducted according to internationally-accepted standards and should always have prior approval from an Institutional Animal Care and Use Committee (IACUC) or equivalent ethics committee(s). Animal care and experimental procedures were reviewed and approved by the Animal Ethics and Compliance Program at the University of Toronto and conducted in accordance with the Canadian Council on Animal Care (CACC) guidelines.



## Fixation

1w 0d 2h 15m

- 1 Immediately after electrophysiology recording is complete, transfer hippocampal slice from recording solution to  2-3 mL fixative solution (4% PFA in PBS) in a small culture dish using a pasteur pipette. Leave on rocker at  Room temperature for  02:00:00  
- 2 Wash 3 x  00:05:00 with 1 X PBS 
- 3 Replace PBS with 30% sucrose solution (w/v in PBS) and seal culture dish with parafilm
- 4 Store at  4 °C for 1-2 weeks 

## Cryo-embedding

5m

- 5 Transfer slice to a cryomold, placing it flat at the bottom using a spatula. Include orientation markers on outside of mold
- 6 Using a rolled kimwipe, soak up any excess sucrose solution without directly touching the slice
- 7 Pour Optimal Cutting Temperature compound (OCT) otop of slice and fill the mold, being careful not to displace the slice. Remove any bubbles in the OCT carefully using a pair of tweezers or spatula
- 8 Embed cryomold in crushed dry ice to flash freeze 
- 9 When frozen, wrap cryomold in parafilm and/or a small plastic bag, label and store at  -80 °C until processing 

## Cryo-sectioning

2h

- 10 Place sample and tools in cryostat for  00:30:00 to acclimatize to  -22 °C 30m
- 11 Pour OCT on specimen chuck and allow it to freeze. When nearly frozen, place a weight on top to flatten surface of OCT
- 12 Remove OCT-embedded sample from cryomold. Place a small amount of OCT on flat surface of chuck and immediately place sample flat on top, with the tissue facing upwards.
- 13 Tighten the chuck in the chuck holder, orienting the tissue with the longest side facing the blade
- 14 Adjust the chuck holder knobs such that the OCT is flush with the blade before slicing
- 15 Slice  20 μm sections, ensuring that the blade remains flush with the sample (not on an angle). Adjust as needed before proceeding.
- 16 Collect each slice from under the anti roll plate and orient it appropriately. Pick slice up swiftly with a  Room temperature microscope slide
- 17 Collect 16-20 slices across 4-5 microscope slides and allow slides to dry for  01:00:00 at  Room temperature 1h
- 18 Store slides at  -20 °C until staining \*

## Immunostaining - Primary Ab Incubation

2h 30m

- 19 Allow slides to acclimatize to  Room temperature for  01:00:00 1h
- 20 Use a hydrophobic barrier PAP pen to outline sections on microscope slide and allow to dry 5m
- 21 Gently wash slides 5 x  00:05:00 with 0.1% Triton-X in PBS (PBS-T) using a pasteur pipette (~  250 μL /slide) 25m  


22 Incubate in blocking solution (10% horse serum in PBS-T) for  01:00:00 at  Room temperature 

1h

23 Incubate in primary antibody cocktail (typically 1:1000 dilution) in blocking solution overnight at  4 °C 

## Immunostaining - Secondary Ab Incubation

4h 20m

24 Wash slides 5 x  00:05:00 with PBS-T 

25m

25 Incubate in secondary antibody cocktail (typically 1:500 dilution) in blocking solution for  02:00:00 at  Room temperature protected from light 

2h

26 Wash slides 5 x  00:05:00 with PBS-T, keeping slides protected from light in between washes 

25m

27 Incubate in DAPI (1:500 in PBS) for  00:10:00 at  Room temperature protected from light 

10m

28 Wash slides 2 x  00:05:00 in PBS, keeping slides protected from light in between washes 

10m

29 Let slides dry  00:10:00 on an angle protected from light 

10m

30 Mount slides with an appropriate mounting medium and coverslip

31 Wait  00:30:00 for mountant to solidify 

30m

32 Seal edges of coverslip with nail polish to prevent from drying out and allow to dry

30m

33 Store in fridge until imaging

\*

## Imaging

34 Clean slides with ethanol

35 Image slices using a confocal microscope. When possible, image and quantify slices that have visible electrode marks to identify the recording plane within the slice to assess differences between electrophysiological conditions

fc