

Feb 16, 2026

Immunofluorescence on Free-Floating Brain Sections

 In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.yxmvm1q5ov3p/v1>

Cristian González-Cabrera¹, Pablo Henny²

¹LIN; ²Universida de Chile



Cristian González-Cabrera

Leibniz Institute für Neurobiologie

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.yxmvm1q5ov3p/v1>

Protocol Citation: Cristian González-Cabrera, Pablo Henny 2026. Immunofluorescence on Free-Floating Brain Sections. protocols.io <https://dx.doi.org/10.17504/protocols.io.yxmvm1q5ov3p/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: February 16, 2026

Last Modified: February 16, 2026

Protocol Integer ID: 243360

Keywords: coronal brain sections for neurochemical characterization, floating coronal brain section, floating brain section, labeled neuron, immunofluorescence on free, neurochemical characterization, immunofluorescence, brain sections this protocol, antigen retrieval

Abstract

This protocol describes immunofluorescence staining of free-floating coronal brain sections for neurochemical characterization of juxtacellularly labeled neurons. Antigen retrieval was performed selectively for VGAT and VGlut2 staining. Confocal images were acquired using identical acquisition settings across experimental conditions.

Materials

****Solutions****

- Phosphate-buffered saline (PBS), 0.01 M.
- Sodium citrate buffer, 0.1 M, pH 6.0.
- Blocking solution: PBS with 3 percent normal horse serum and 0.3 percent Triton X-100.
- Fluorescence mounting medium.

****Primary Antibodies****

- Guinea pig anti-Tyrosine Hydroxylase (TH), 1:5000, Synaptic Systems.
- Rabbit anti-VGlut2, 1:2000, Synaptic Systems.
- Rabbit anti-VGAT, 1:5000, Synaptic Systems.

****Secondary Antibodies****

- Species-appropriate Alexa Fluor 488 or Alexa Fluor 635 conjugates.
- Dilution 1:1000, Jackson ImmunoResearch.

****Tissue Preparation****

- Free-floating coronal sections.
- Section thickness: 40 μ m.



Procedure

- 1 Rinse sections 3 times in PBS.
5 min per wash.
- 2 **Antigen retrieval**
Incubate sections in 0.1 M sodium citrate buffer, pH 6.0.
80 °C for 30 min.
Rinse 3 times in PBS.
Note: Antigen retrieval was performed only for VGAT and VGlut2 staining.
- 3 Incubate sections for 1 h at room temperature.
Blocking buffer: PBS + 3 percent normal horse serum + 0.3 percent Triton X-100.
- 4 Incubate overnight at 4 °C.
Primary antibodies diluted in blocking buffer.
- 5 Rinse sections 3 times in PBS.
- 6 Incubate 2 h at room temperature.
Protect from light.
Secondary antibodies diluted 1:1000.
- 7 Rinse sections 3 times in PBS.
Mount onto glass slides.
Coverslip using fluorescence mounting medium.

Imaging

- 8 Confocal microscopy.
Identical laser power and detector settings across conditions.
Settings kept constant within antibody sets.

Critical Steps

- 9 Perform antigen retrieval only for VGAT and VGlut2.
Maintain consistent antibody dilutions across experiments.
Protect fluorescent samples from light during and after secondary incubation.
Use identical confocal settings for comparisons.



Expected Outcome

- 10 Successful staining yields clear neurochemical identification of labeled neurons, allowing classification based on TH, VGluT2, or VGAT expression.