



Dec 03, 2024

Immunofluorescence - RAW 264.7

 [The Journal of Cell Biology](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.36wgqd7wkvk5/v1>

Amanda Bentley-DeSousa^{1,2,3,4,5}, Shawn Ferguson^{1,2,3,4,5,6}

¹Department of Cell Biology, Yale University School of Medicine, New Haven, Connecticut 06510, USA;

²Neuroscience, Yale University School of Medicine, New Haven, Connecticut 06510, USA;

³Program in Cellular Neuroscience, Neurodegeneration and Repair;

⁴Wu Tsai Institute Yale University School of Medicine, New Haven, Connecticut 06510, USA;

⁵Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815, USA;

⁶Kavli Institute for Neuroscience, Yale University School of Medicine, New Haven, Connecticut 06510, USA



Amanda Bentley-DeSousa

Yale University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.36wgqd7wkvk5/v1>

Protocol Citation: Amanda Bentley-DeSousa, Shawn Ferguson 2024. Immunofluorescence - RAW 264.7. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.36wgqd7wkvk5/v1>

Manuscript citation:

<https://doi.org/10.1101/2023.10.31.564602>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 02, 2024

Last Modified: December 03, 2024

Protocol Integer ID: 113418

Keywords: RAW 264.7, immunofluorescence, immunofluorescence in raw, protocol details step, cell, raw, protocol detail

Funders Acknowledgements:

Aligning Science Across Parkinson's Disease

Grant ID: ASAP-000580

Abstract

This protocol details steps taken to perform immunofluorescence in RAW 264.7 cells.

Materials

DMEM (Thermo Fisher Scientific, 11965-092)

FBS (Thermo Fisher Scientific, 16140-071)

Penicillin/Streptomycin (Thermo Fisher Scientific, 15140122)

CellStripper (Corning, 356230)

PBS [1.1 mM KH₂PO₄ (J.T. Baker, 3246-01), 155.2 mM NaCl (Sigma-Aldrich, 3624-05), and 3 mM Na₂HPO₄ (J.T. Baker, 3828-05)]

Paraformaldehyde (Electron Microscopy Sciences, 19202)

Sodium phosphate buffer (pH 7.3 / Buffer: 153.56 mM Sodium phosphate, dibasic, anhydrous, J.T. Baker 3828, 53.63 mM Sodium dihydrogen phosphate monohydrate, J.T. Baker 3818)

Sucrose (Sigma-Aldrich, S0389)

BSA (Sigma-Aldrich, A9647)

Saponin (Sigma-Aldrich, S4521)

Triton (Sigma-Aldrich, X100)

Coverslips (12 mm, Carolina Biological Supplies, 633029)

Microscope Slides (Thermo Fisher Scientific, 12-550-143)

Prolong Gold Mounting Media (Thermo Fisher Scientific, P36935)



Day 0

- 1 Plate 2×10^5 RAW 264.7 cells on poly-D-lysine coated coverslips.

Day 1

35m

- 2 Add the indicated compound at the provided concentration and time point in fresh DMEM media containing FBS and P/S.
- 2.1 Proceed to next step if not treating cells.
- 3 Fix cells in 4% paraformaldehyde supplemented with 4% sucrose for 00:30:00 at Room temperature .
- 4 Wash coverslips 3X for 00:05:00 with PBS at Room temperature .
- 5 Add primary antibodies Overnight at 4 °C in 3% BSA in PBS.

Day 2

1h 10m

- 6 Wash coverslips 3X for 00:05:00 with PBS at room temperature.
- 7 Add secondary antibody for 01:00:00 at Room temperature in PBS in the dark.
- 8 Wash coverslips 3X for 00:05:00 with PBS at Room temperature .
- 9 Mount coverslips onto microscope slides with Prolong Gold mounting media and store at 4 °C .



Day 3+

10 Image coverslips.