



Dec 26, 2025

Immunofluorescence staining and confocal microscopy in cultured cells

DOI

<https://dx.doi.org/10.17504/protocols.io.yxmvm177nv3p/v1>

Xin Xu¹, Magnus Lam¹

¹Princess Margaret Cancer Centre



Xin Xu

University Health Network

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.yxmvm177nv3p/v1>

Protocol Citation: Xin Xu, Magnus Lam 2025. Immunofluorescence staining and confocal microscopy in cultured cells. protocols.io <https://dx.doi.org/10.17504/protocols.io.yxmvm177nv3p/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 22, 2025



Last Modified: December 26, 2025

Protocol Integer ID: 235612

Keywords: immunofluorescence, confocal microscopy, protein localization, fluorescence imaging, cell biology, proximity ligation assay, PLA, cancer cells, confocal microscopy in cultured cell, transduced cultured cell, cultured cell, cell, cells this protocol, antibody

Disclaimer

This protocol is provided for research use only and may require optimization depending on experimental conditions.

Abstract

This protocol describes immunofluorescence staining and confocal microscopy in lentivirally transduced cultured cells, including fixation, permeabilization, antibody staining, and mounting for fluorescence imaging.

Guidelines

- Antibody dilutions and incubation times may require optimization depending on antibody performance.
- Avoid drying samples during staining and wash steps.
- PLA experiments should follow the manufacturer's protocol exactly.

Materials

- Lab-Tek II Chamber 8-well slides (Thermo Fisher Scientific, 154534)
- Alexa Fluor 647-conjugated anti-Rabbit IgG (Thermo Fisher Scientific, A-31573)
- Alexa Fluor 488-conjugated anti-Mouse IgG (Thermo Fisher Scientific, A-11001)
- Triton X-100 (Millipore Sigma, X100-500ML)
- 4% paraformaldehyde (Thermo Fisher Scientific, J61899-AK)
- ProLongTM Diamond Antifade Mountant containing DAPI (Thermo Fisher Scientific, P36971)
- glass coverslips (Fisher Scientific, 12541032CA)

Before start

- Perform all staining steps at room temperature unless otherwise indicated.
- Prepare all buffers fresh or thaw on the day of use.
- Protect samples from light after secondary antibody incubation.
- Use appropriate biosafety precautions when handling paraformaldehyde.



Cell seeding

- 1 Seed lentivirally transduced cells onto Lab-Tek II Chamber 8-well slides (Thermo Fisher Scientific, 154534).
- 2 Incubate cells for 24 h to reach 40–60% confluency at the time of fixation.

Fixation

20m

- 3 Rinse cells three times with PBS.
- 4 Fix cells with 4% paraformaldehyde for 15 min at room temperature.
- 5 Wash cells three times with PBS.

15m

5m

Permeabilization

10m

- 6 Permeabilize cells with 0.5% Triton X-100 in PBS for 5 min at room temperature.
- 7 Wash cells three times with PBS using gentle agitation.

5m

5m

Blocking

1h

- 8 Block cells with 3% BSA and 0.1% Tween-20 in PBS for 1 h at room temperature.

1h

Primary antibody incubation

2h

- 9 Dilute primary antibodies in 1% BSA and 0.1% Tween-20 in PBS.



10 Incubate cells with primary antibodies for 2 h at room temperature.

2h

Secondary antibody incubation

5h

11 Dilute secondary antibodies in 1% BSA and 0.1% Tween-20 in PBS:
- Alexa Fluor 647–conjugated anti-Rabbit IgG (1:1000; Thermo Fisher Scientific, A-31573)
- Alexa Fluor 488–conjugated anti-Mouse IgG (1:1000; Thermo Fisher Scientific, A-11001)

12 Incubate cells with secondary antibodies for 1.5 h at room temperature, protected from light.

5h

13 Wash cells three times with PBS.

Mounting

14 Mount cells with ProLong Diamond Antifade Mountant containing DAPI (Thermo Fisher Scientific, P36971).

15 Cover samples with glass coverslips (Fisher Scientific, 12541032CA).

16 Cure mounted slides at room temperature, protected from light, for at least 24 h before sealing with nail polish.



Imaging

17 Acquire images using confocal microscopy under appropriate laser and detector settings.

Proximity ligation assay (PLA) (optional)

18 Perform the proximity ligation assay (PLA) according to the manufacturer's instructions (Sigma Aldrich, DUO92101).