

Feb 03, 2025

Immunofluorescence staining and image analysis

 In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.n92ldr3k7g5b/v1>

Maria Jose Perez J.¹

¹Institut Imagine

Team Deleidi



Federico Bertoli

ASAP

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n92ldr3k7g5b/v1>

Protocol Citation: Maria Jose Perez J. 2025. Immunofluorescence staining and image analysis. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.n92ldr3k7g5b/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: November 06, 2024

Last Modified: February 03, 2025

Protocol Integer ID: 111683

Keywords: image analysis immunofluorescence, staining

Funders Acknowledgements:

ASAP

Abstract

Immunofluorescence staining and image analysis

2D cultures IF

- 1 Seed cells on coverslips and fix with 4% paraformaldehyde (PFA) for 10 minutes
- 2 Rinse cells with phosphate-buffered saline (PBS)
- 3 Block cells with 10% normal goat or donkey serum (NGS/NDS) in PBS containing 0.1% Triton X-100 (PBST) at room temperature for 1 hour
- 4 Incubate cells overnight at 4°C with primary antibodies diluted in PBST containing 5% NGS/NDS
- 5 The next day, after washing with PBS 3 times for 5', incubate cells for 2 hours at room temperature with Alexa 488/568/647-conjugated secondary antibodies at a 1:1000 dilution, diluted in PBST containing 5% NGS/NDS
- 6 Remove secondary antibodies and stain nuclei with 1 μ M DAPI for 5 minutes
- 7 Wash coverslips 3 times for 10'
- 8 Mount coverslips on glass slides with Dako mounting medium

Brain organoids and Assemblods IF

- 9 For brain organoids and assemblods, fix samples with 4% (w/v) PFA for 1 hour
- 10 Equilibrate organoids in 30% sucrose in PBS overnight at 4°C
- 11 Embed organoids in OCT compound the following day
- 12 Cryosection organoids to 20 μ m thickness and mount on SuperFrost Plus slides



- 13 Block sections with 10% NGS in PBS containing 0.5% Triton X-100
- 14 Incubate sections with primary antibodies diluted in PBST containing 5% NGS/NDS overnight
- 15 After washing sections with PBS 3 times for 5' Incubate with secondary antibodies diluted in PBST containing 5% NGS/NDS for 3 hours
- 16 (If BODIPY staining is needed, incubate cells or tissues with 1 μ M BODIPY 493/503 dye in PBS for 30 minutes at room temperature)
- 17 (For protein aggregate detection, stain sections with 10 μ M thioflavin-T for 30 minutes at room temperature)
- 18 Stain sections with 1 μ M DAPI for 5 minutes
- 19 Wash sections three times with PBS for 10'
- 20 Mount sections with Dako mounting medium

Mouse tissue IF

- 21 For mouse brain tissue, perform antigen retrieval on free-floating brain slices for 20 minutes in citrate buffer at 80°C
- 22 Block tissues with mouse-on-mouse (MOM) blocking reagent following manufacturer's instructions
- 23 Incubate with primary antibodies overnight using MOM protein concentrate
- 24 After washing in PBS 3 times 5', incubate with secondary antibodies for 3 hours
- 25 Stain sections with 1 μ M DAPI for 5 minutes



26 Wash sections three times with PBS for 10'

27 Mount sections with Dako mounting medium

Image analysis

28 Acquire images using a Leica TCS SP8 confocal microscope with a 60×/1.4 numerical aperture objective or an EVOS M7000 Imaging System with a 20×/0.4 numerical aperture

29 For quantification, analyze p21-, p16-, and γH2AX-positive nuclei in CellProfiler with DAPI as a counterstain

30 Calculate fluorescence intensity per area using Fiji