



Apr 03, 2023

Immunofluorescence staining

 In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.36wgg7oq5vk5/v1>

michela.deleidi¹, Pascale Baden¹, Maria Jose Perez J.¹, Hariam Raji¹, Federico Bertoli¹

¹German Center for Neurodegenerative Diseases (DZNE), Tübingen, 72076 Germany



Hariam Raji

Institut Imagine Paris - Inserm U1163, Aligning Science Acro...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.36wgg7oq5vk5/v1>

Protocol Citation: michela.deleidi , Pascale Baden, Maria Jose Perez J., Hariam Raji, Federico Bertoli 2023.

Immunofluorescence staining. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.36wgg7oq5vk5/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: April 26, 2022

Last Modified: May 31, 2024

Protocol Integer ID: 61418

Keywords: Fixation of the cells, HEK cells, neurons, Blocking and permeabilization, immunofluorescence Staining, ASAPCRN, immunofluorescence, cell

Abstract

This protocol describes the immunofluorescence staining of cells.

Attachments



[404-872.docx](#)

24KB



Materials

Recipes and products:

	A	B
	8% Paraformaldehyde (PFA)	
	PFA	20 g
	1M NaOH	0.5 ml
	1x PBS	100 ml

Note

- Heat to ~ 60°C to dilute. Then filter through folded filters into new cylinder.
- Adjust $\text{pH } 7.4$ (normally its ~ $\text{pH } 7.38$ without adding something)
- Fill up to 250 mL with 1x PBS.

Dako Mounting Medium  Fluorescence Mounting Medium **Agilent Technologies Catalog #S302380-2** :

NGS (normal-goat serum)  Normal Goat Serum Blocking Solution **BIOZOL Catalog #VEC-S-1000**

Note

- Prepare aliquots out of stock and store in -20°C .
- Filter before use to avoid contamination.

DAPI (DAPI (4',6-Diamidino-2-Phenylindole, Dilactate)

 DAPI (46-Diamidino-2-Phenylindole Dilactate) **BioLegend Catalog #422801** :

- Dissolve the content in 2 mL deionized water (DAPI concentration $10.9\text{ millimolar (mM)}$).



Fixation of the cells: Strategy 1, i.e., HEK cells

25m

1 Remove medium.

2 Wash with 1× PBS.



3 Remove PBS and add 4% PFA to the cells.



Note

Note: For a coverslip in a 24 multi-well use at least  300 µL /well.

4 Incubate  00:10:00 at  Room temperature .

10m



5 Collect PFA.

6 Wash 1× PBS.



6.1 Wash with 1× PBS for  00:05:00 at  Room temperature . (1/3)

5m

6.2 Wash with 1× PBS for  00:05:00 at  Room temperature . (2/3)

5m

6.3 Wash with 1× PBS for  00:05:00 at  Room temperature . (3/3)

5m

Fixation of the cells: Strategy 2, i.e., neurons

25m

7 Add the same volume of 8% PFA as medium is in the well to the well.





8 Incubate  00:10:00 at  Room temperature .

10m



9 Collect PFA in a falcon.

10 Wash in 1× PBS.



10.1 Wash with 1× PBS for  00:05:00 at  Room temperature . (1/2)

5m

10.2 Wash with 1× PBS for  00:05:00 at  Room temperature . (2/2)

5m

Note

Storage until ICC: Keep coverslips in 1× PBS, seal the plate with parafilm and store at  4 °C .

Blocking and permeabilization

25m

11 Remove 1X PBS.

12 Add  300 µL /well of blocking solution.



Note

Blocking solution: 10% NGS [normal goat serum] in PBS + Triton X-100 0,1%, filter the solution before using it.

13 Incubate at least  01:00:00 at  Room temperature .

1h





Staining: Day 1

25m

- 14 Prepare antibody in blocking solution containing 5% NGS.
- 15 Put a drop ( 50 μ L) of Primary antibody solution on the parafilm surface.
- 16 Remove the coverslips from the plate and gently put it upside-down on the antibody drop.
- 17 Incubate  Overnight at  4 °C .

1h



Staining: Day 2

25m

- 18 Wash: 
- 18.1 Wash for  00:05:00 in PBS + Triton X-100 0.1%. (1/3) 5m
- 18.2 Wash for  00:05:00 in PBS + Triton X-100 0.1%. (2/3) 5m
- 18.3 Wash for  00:05:00 in PBS + Triton X-100 0.1%. (3/3) 5m
- 19 Prepare secondary antibody in blocking solution containing 5% NGS.

Note

Note: Keep in the dark.

- 20 Put a drop ( 50 μ L) of secondary antibody solution on the parafilm.
- 21 Take the coverslips and put it upside-down on the antibody drop.



22 Incubate  01:00:00 at  Room temperature in the dark.

1h



23 Transfer the coverslip to a 24-well containing  500 μ L 1X PBS + Triton X-100 0,1%.

24 Incubate for  00:05:00 at  Room temperature (in the dark).

5m



25 Dilute DAPI 1:10000 in 1X PBS.

26 Remove PBS and incubate with DAPI for  00:05:00 at  Room temperature in the dark.

5m



27 Wash.



27.1 Wash with 1X PBS. (1/2)

27.2 Wash with 1X PBS. (2/2)

28 Mount the slides:

28.1 Put a drop ( 10 μ L) of DAKO mounting reagent on the slide.

28.2 Take out the coverslip from the plate.

28.3 Dry it by gently tapping the coverslip's edge on a lens-cleaner tissue.

28.4 Gently put the coverslips upside-down on the DAKO drop.



28.5 Leave it dry (🕒 24:00:00 , DARK, 🌡 Room temperature).

1d

29 Store in a slide-box at 🌡 4 °C .