

Apr 29, 2024

Immunofluorescence Staining in Mouse Brain Tissue Sections

DOI

<https://dx.doi.org/10.17504/protocols.io.14egn6zopl5d/v1>

madalynn.erb Erb¹

¹Van Andel Research Institute

Team Lee



Jane Balster

ASAP - Team Lee

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.14egn6zopl5d/v1>

Protocol Citation: madalynn.erb Erb 2024. Immunofluorescence Staining in Mouse Brain Tissue Sections. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.14egn6zopl5d/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 08, 2024

Last Modified: May 31, 2024

Protocol Integer ID: 98951

Keywords: ASAPCRN, immunofluorescence staining in mouse brain tissue section, immunofluorescence staining, mouse brain tissue section, immunofluorescence, staining, tissue

Funders Acknowledgements:

aligning science across parkinsons

Grant ID: 000592

Abstract

This protocol details the immunofluorescence staining in mouse brain tissue sections.

Materials

Pyrex® Staining Dish **Ted Pella Inc. Catalog #36754-60**

8-Section Staining Nets **Ted Pella Inc. Catalog #36154-64**

Fisherbrand™ Superfrost™ Plus Microscope Slides **Fisher Scientific Catalog #12-550-15**

ProLong™ Diamond Antifade Mountant **Invitrogen - Thermo Fisher Catalog #P36961**

Primary Antibodies:

	A	B	C	D
	Target	Species	Conc	Manufacturer
	TH	Rabbit	1:2000	Novus Biological N300109
	TH	Chicken	1:750	Abcam ab76442
	GFP	Chicken	1:2000	Aves Labs- GFP-1010
	GFP	Rabbit	1:500	Thermofisher A-11122
	GFP	Mouse	1:1000	Roche 11814460001
	P62	Guinea pig	1:2000	Progen GP62-C
	LAMP2	Rat	1:1000	Abcam ab13524
	GAD67	Mouse	1:500	Millipore Sigma MAB5406
	parvalbumin	Rabbit	1:500	Abcam ab11427
	TFE3	Rabbit	1:500	Abcam ab93808

Secondary Antibodies (From Thermofisher):



	A	B	C	D
	Antibody	Fluorophore	Concentration	Cat Number
	goat-anti rabbit	488	1:500	A-11008
	goat-anti rabbit	546	1:500	A-11010
	goat-anti rabbit	647	1:500	A-21245
	goat-anti mouse	488	1:500	A-11029
	goat-anti mouse	546	1:500	A-11003
	goat-anti rat	647	1:500	A-21247
	goat-anti chicken	488	1:500	A-11039
	goat-anti chicken	647	1:500	A-21449



Day 1

30m

- 1 Staining of 35µm free-floating mouse brain sections is performed in glass staining dishes (Pyrex 36754-60) using 8-section staining nets (Ted Pella 36154-64).
[dx.doi.org/10.17504/protocols.io.5jyl8pzk9g2w/v1](https://doi.org/10.17504/protocols.io.5jyl8pzk9g2w/v1)
- 1.1 The sections can be transferred between wells using a paint brush.
- 1.2 Volume of solution:
 -  20 mL -  30 mL for antibodies.
 -  50 mL for washing.
- 2 Wash sections 3 times (5 min each wash) in Phosphate Buffer Saline (PBS) at  Room temperature to remove cryo-protectant solution.
- 2.1 Wash sections for  00:05:00 in Phosphate Buffer Saline (PBS) at  Room temperature to remove cryo-protectant solution (1/3).
- 2.2 Wash sections for  00:05:00 in Phosphate Buffer Saline (PBS) at  Room temperature to remove cryo-protectant solution (2/3).
- 2.3 Wash sections for  00:05:00 in Phosphate Buffer Saline (PBS) at  Room temperature to remove cryo-protectant solution (3/3).
- 3 Wash sections 3 times (5 min each wash) in PBS + 0.1% triton X-100.
- 3.1 Wash sections for  00:05:00 in PBS + 0.1% triton X-100 (1/3).
- 3.2 Wash sections for  00:05:00 in PBS + 0.1% triton X-100 (2/3).
- 3.3 Wash sections for  00:05:00 in PBS + 0.1% triton X-100 (3/3).



5m



5m



5m



5m



5m



5m



4 Block sections for ⌚ 01:00:00 in PBS + 10% Normal Goat Serum (NGS) + 0.4% Bovine Serum Albumin (BSA) and 0.2% triton X-100 for ⌚ 01:00:00 at 🌡️ Room temperature .

2h

5 Incubate sections with primary antibodies in blocking solution (PBS + 0.1% triton X-100, 0.4% BSA) ⌚ Overnight at 🌡️ 4 °C .

1h



Day 2

2d 2h 5m

6 Wash sections 3 times (10 min each wash) in PBS + 0.1% triton X-100 at 🌡️ Room temperature .



6.1 Wash sections for ⌚ 00:10:00 in PBS + 0.1% triton X-100 at 🌡️ Room temperature (1/3).

10m



6.2 Wash sections for ⌚ 00:10:00 in PBS + 0.1% triton X-100 at 🌡️ Room temperature (2/3).

10m



6.3 Wash sections for ⌚ 00:10:00 in PBS + 0.1% triton X-100 at 🌡️ Room temperature (3/3).

10m



7 Incubate with fluorescent secondary antibodies (1:500) diluted in PBS + 0.1% triton X-100 + 0.4% BSA for 2 hours at 🌡️ Room temperature or ⌚ 24:00:00 at 🌡️ 4 °C .

1d



Note

Sections should be protected from light at this step and all proceeding steps.

8 Wash sections 3 times (10 min each wash) in PBS + 0.1% triton X-100 at 🌡️ Room temperature .





- 8.1 Wash sections for  00:10:00 in PBS + 0.1% triton X-100 at  Room temperature  
- 8.2 Wash sections for  00:10:00 in PBS + 0.1% triton X-100 at  Room temperature 
- 8.3 Wash sections for  00:10:00 in PBS + 0.1% triton X-100 at  Room temperature  
- 9 Incubate with DAPI in PBS (1:5000) for  00:30:00 at  Room temperature .  
- 10 Wash sections 3 times (10 min each wash) in PBS at  Room temperature . 
- 10.1 Wash sections for  00:10:00 in PBS at  Room temperature (1/3).  
- 10.2 Wash sections for  00:10:00 in PBS at  Room temperature (2/3).  
- 10.3 Wash sections for  00:10:00 in PBS at  Room temperature (3/3).  
- 11 Mount sections:
- 11.1 Use a petri dish filled with PBS to mount the sections on superfrost plus slides (Fisher Scientific 12-550-15).
- 11.2 This is easiest to do using a medium sized paint brush.
- 11.3 Let sections dry  00:02:00 -  00:05:00 . 

**Note**

Coverslip using ProLongTM Diamond Antifade Mountant (Thermofisher P36961).

12 Dry slides for  24:00:00 at  Room temperature .

1d

1. Protect slides from light at this step.

2. Store slides at  4 °C .

