



Jun 26, 2024

Version 1

# Immunofluorescence staining on larval and adult Drosophila gonads V.1

DOI

[dx.doi.org/10.17504/protocols.io.n2bvjnd7ngk5/v1](https://dx.doi.org/10.17504/protocols.io.n2bvjnd7ngk5/v1)

Samantha Goetting<sup>1</sup>

<sup>1</sup>Johns Hopkins University



Samantha Goetting

Johns Hopkins University

## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n2bvjnd7ngk5/v1>

**Protocol Citation:** Samantha Goetting 2024. Immunofluorescence staining on larval and adult Drosophila gonads. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.n2bvjnd7ngk5/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working

**Created:** June 26, 2024



**Last Modified:** June 26, 2024

**Protocol Integer ID:** 102463

**Keywords:** drosophila, immunofluorescence, staining protocol, larval, adult drosophila

## Abstract

Immunofluorescence staining protocol for *Drosophila* gonads

## Materials

1x PBS

0.3% PBTx (0.3% Triton-X in PBS)

1% PBTx (1% Triton-X in PBS)

Paraformaldehyde or formaldehyde

Normal serum (usually NGS)

Primary antibodies

Secondary antibodies

DAPI

## Troubleshooting



## Before start

All steps are done with gentle rotation.



## Day 1

2h 40m

- 1 Dissect tissue in 1x PBS. Transfer to a 1.5 mL tube containing **1x PBS**. If it is a quick dissection (<20 mins), no ice needed. If more time is needed, keep samples on ice.
- 2 Remove PBS and add fixative. Fix in 4% paraformaldehyde in **0.3% PBTx** for 20 min RT with gentle rotation.  
500 uL fixative = 125 uL of 16% paraformaldehyde + 375 uL 0.3% PBTx 
- 3 Aspirate the fixative and wash twice for 10 min in **1% PBTx**. Not getting rid of fix will affect your immunostaining. 
- 4 Aspirate the supernatant and block/permeabilize for at least 2 hours in **1% PBTx** + 5% normal serum, or overnight at 4°C.  
1 mL block solution = 50 uL NGS + 950 uL 1% PBTx
- 5 Primary antibodies are diluted in **0.3% PBTx** + 5% normal serum and incubate for 1 hour at RT or overnight at 4°C. Overnight will give better staining.  
Primary antibody mix in 0.3% PBTx + 15 uL NGS (300 uL total)

## Day 2

5h 20m

- 6 Remove the primary antibody mix and wash in **0.3% PBTx** three times for 20 min at RT.
- 7 Wash in **0.3% PBTx** + 5% normal serum twice for 30 min at RT.  
1 mL wash solution = 50 uL NGS + 950 uL 0.3% PBTx
- 8 Secondary antibodies are diluted in **0.3% PBTx** + 5% normal serum and incubated for ~2 hours at RT or overnight at 4°C. Keep tubes covered from light.
- 9 Aspirate the supernatant and add 500 uL DAPI to each tube. Incubate for 10 min at RT. Keep tubes covered from light.
- 10 Aspirate the supernatant and wash in **0.3% PBTx** three times for 20 min at RT.
- 11 Aspirate the supernatant and wash in **PBS** for 10 min at RT.
- 12 Store in **PBS** at 4°C or proceed with mounting. Keep tubes covered from light.



## Protocol references

Slaidina et al. (2020)