



Jan 16, 2026

Version 1

Immunofluorescence staining –Paraffin sections V.1

DOI

<https://dx.doi.org/10.17504/protocols.io.bp2l6ex9rgqe/v1>

Domenic Filingeri¹

¹University of Pittsburgh



Domenic DF Filingeri

University of Pittsburgh, Children's Hospital of Pittsburgh

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bp2l6ex9rgqe/v1>

Protocol Citation: Domenic Filingeri 2026. Immunofluorescence staining –Paraffin sections. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.bp2l6ex9rgqe/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 13, 2026



Last Modified: January 16, 2026

Protocol Integer ID: 238596

Keywords: paraffin section, immunofluorescence, dapi step

Abstract

IF staining for paraffin sections (DAPI step included)

Materials

- HistoClear
- Ethanol (100%, 95%, 90%, 80%, 70%, 50%)
- Tap water
- Triton X-100
- PBS
- Sodium citrate buffer (10 mM Sodium citrate, 0.05% Tween 20, pH 6.0)
- EDTA buffer pH 8.0
- Tris-EDTA buffer pH 9.0
- Blocking buffer
- Primary antibody
- Secondary antibody
- 1xDAPI
- Anti-fade mounting medium
- Sodium citrate (dihydrate) 2.94 g
- Distilled water
- 1N HCl
- Tween 20
- EDTA 0.37 g
- NaOH
- Tris 1.21 g
- BSA
- NP40
- Normal Serum



- 1 60°C for 25 min in dry incubator 25m 
- 2 Re-hydrate with HistoClear and Ethanol, create each as a 50mL conical tube, pour into glass slide box. Pour slowly and not directly onto slides. 45m 
HistoClear 10 min X 2 times
100% EtOH 5 min X 2 times
90% EtOH 3 min, (45mL in 5mL DI water)
80% EtOH 3 min, (40mL in 10mL DI water)
70% EtOH 3 min, (35mL in 15mL DI water)
50% EtOH 3 min, (25mL in 25mL DI water)
- 3 Wash in tap water, 5 min X 3 times 15m 
- 4 **Antigen retrieval.** Select one method from below. Use plastic slide box
 - 4.1 Sodium citrate (pH 2.0, **pH6.0**, or pH 9.0)
 - 4.2 Tris-EDTA buffer pH 9.0
 - 4.3 EDTA buffer pH 8.0
- 5 Heat them using a presser cooker for 20 minutes in plastic box. Add DI-water to pressure cooker up to metal plate. Select porridge option, remove at "LO 25min", then RT for 1 hour 45m 
- 5.1 Allow to cool to RT slowly. This can take up to 2 hours. (Minimum 30 min) 2h
- 6 Wash in tap water, 5 min X 2 times 10m 
- 7 Wipe away water, circle section with pap-pen, then add 200uL/slide of tap water. (Use 200uL for all future steps unless noted) 



- 7.1 cut P-200 pipette tip to prevent high pressure for subsequent steps
- 8 Permeabilize nuclear membrane with 0.25% Triton X-100 in PBS 5 min 5m 
- 9 Wash in tap water 10 min X 2 times 20m  
- 10 Incubate slides with **blocking buffer** for 30min - 1 hour at RT in hydration chamber (add tap water to wells between slides). 1h  
- 10.1 Blocking buffer is 5%BSA in TBST
- 11 Rinse in PBST briefly  
- 12 Incubate for **primary antibody** diluted in blocking buffer. Leave overnight at 4°C in hydration chamber. (add tap water to wells between slides).  
- 12.1 Antibody diluted in 3%BSA in TBST
- 13 Wash in PBST 3 min X 3 times. Double wash 10m  
- 14 Incubate with **secondary antibody** diluted in 1x PBS at RT for 30 min to one hour in hydration chamber (dilute in PBS). 1h  
- 14.1 Antibody diluted in 3%BSA in TBST
- 15 Wash in PBST 3 min X 3 times. Double wash 10m  



16 Incubate with 1xDAPI at RT for 15 min in a hydration chamber. (Dilute 1:100 in ddH₂O)

15m



17 Wash in PBST 3 min X 3 times. Double wash

10m



18 Add one small drop anti-fade mounting medium (10-15 μ L) to each section, cover, press with plastic pipette tip to remove bubbles, blot clean with paper towel. Store at least 15', then move to 4°C



1 mM EDTA, pH 8.0

19 EDTA 0.37 g
Distilled water 1 L
Adjust to pH 8.0 with NaOH
Store at room temperature for 3 months

Tris-EDTA buffer (10 mM Tris base, 1 mM EDTA solution, 0.05% Tween 20, pH 9.0)

20 Tris 1.21 g
EDTA 0.37 g
Distilled water 1 L
Mix to dissolve. Adjust pH to 9.0.
Add 0.5 mL of Tween 20 and mix well. Store at room temperature for 3 months or at 4°C for longer storage.

Blocking buffer

21 10% Normal Serum
1% BSA
0.25% NP40
In 1x PBS

PBST

22 1x PBS
0.1% Tween 20

4x Blocking Buffer stock (50 ml)

23 0.5g BSA in 40 ml of 1xPBS. Mix to dissolve. It can be warmed in a water bath.



Add 1.25 ml of 10% NP40 and mix by inverting. Do not vortexing. It causes a lot of bubbles.

Add 1x PBS up to 50 ml Vol.

Aliquot (250 ul in 1.5 ml tube) and store at -20 C.

Working blocking buffer

- 24 250ul of 4x blocking buffer
100 ul of normal serum (if all secondary antibodies are from one species)
650 ul of 1x PBS

Sodium Citrate

- 25 2.94g Tri-Sodium Citrate
800 mL distilled (pure) Water
Mix, dissolve
Adjust pH to 6.0 with 1M HCl
Top off to 200 mL using Pure diH₂O
Add 0.5 mL Tween 20, or 5mL 10% Tween

PBST

- 26 990 mL PBS
10 mL (10% Tween 20) or 1mL (100% tween)

.25% Triton X 100 in PBS

- 27 10x PBS (900 mL Pure, 100 mL PBS 10X)
2.5 mL Triton X 100
25 mL 10x Triton

Acknowledgements

Edited by Jeong Lee on 6/28/2022